



# FavorPrep™ Circulating Nucleic Acid Extraction HE Mini Kit

## ■ Kit Contents

| Cat. No.  | FACN0030<br>(4 Preps) | FACN0033<br>(50 Preps) | FACN0034<br>(100 Preps) |
|---|-----------------------|------------------------|-------------------------|
| CL Lysis Buffer   | 15 ml                 | 200 ml                 | 400 ml                  |
| CB Binding Buffer (Concentrate) ▲   | 25 ml                 | 135 ml × 2             | 270 ml × 2              |
| Wash Buffer 1 (Concentrate) ■   | 0.48 ml × 2           | 12 ml                  | 24 ml                   |
| Wash Buffer 2 (Concentrate) ●   | 3 ml                  | 25 ml                  | 50 ml                   |
| RNase-Free Water  | 0.5 ml                | 5 ml                   | 8 ml                    |
| Proteinase K (Liquid)   | 1050 µl × 2           | 10.5 ml × 2            | 10.5 ml × 4             |
| HE Column R   | 4 pcs                 | 50 pcs                 | 50 pcs × 2              |
| HE Collection Tubes   | 4 pcs                 | 50 pcs                 | 100 pcs                 |
| Elution Tubes   | 4 pcs                 | 50 pcs                 | 100 pcs                 |
| Column Extenders  | 4 pcs                 | 50 pcs                 | 100 pcs                 |
| User Manual   | 1                     | 1                      | 1                       |
| <b>Preparation of CB Binding Buffer, Wash Buffer 1 and Wash Buffer 2.</b> |                       |                        |                         |
| Volume of 100% Isopropanol for CB Binding Buffer ▲                        | 17 ml                 | 90 ml                  | 180 ml                  |
| Volume of 96~100% Ethanol for Wash Buffer 1 ■                             | 0.72 ml               | 18 ml                  | 36 ml                   |
| Volume of 96~100% Ethanol for Wash Buffer 2 ●                             | 12 ml                 | 100 ml                 | 200 ml                  |


All kit components are shipped at room temperature and should be stored at room temperature between 15~25°C upon receipt.

## ■ Specification


|                  |   |
|------------------|---|
| Format/Principle | Spin column (Silica matrix)   |
| Binding Capacity | ≤150 µg Nucleic acid/Column   |
| Operation Time   | <60 mins  |
| Sample Size      | 1~4 ml cell-free fluid sample such as serums, plasma, or body fluids. |
| Elution Volume   | 30 µl   |

## ■ Procedure Overview

Sample




- Add 1~4 ml cell-free fluid sample into a tube.
- Add 0.1x sample volume of **Proteinase K**.
- Add 0.8x sample volume of **CL Lysis Buffer**.
- Mix thoroughly, incubate at 56°C for 30 mins (vortex 2~3 times).



- Add 1.8x volume of **CB Binding Buffer** (isopropanol contained).
- Incubate the sample mixture on ice for 5 mins.




Vacuum  
-650 mmHg



- Connect an **HE Column R** with a **Column Extender**.
- Transfer all the sample mixture into the HE Column R assembly.
- Pass the sample mixture through the **HE Column R** by vacuum pressure.



Vacuum  
-650 mmHg



- Detach the **Column Extender**.
- Add 500 µl **Wash Buffer 1** (ethanol contained).
- Add 900 µl **Wash Buffer 2** (ethanol contained), repeat.



Centrifuge  
18,000 xg, 2 mins



- Place the **HE Column R** to an **HE Collection Tube**.
- Dry the column membrane.



Centrifuge  
18,000 xg, 2 mins



- Add 30 µl **RNase-Free Water**.
- Stand the column for 5 mins.
- Obtain purified nucleic acid.

## ■ Preparation Before Starting

1. Buffers provided in this system contain irritants. Wear gloves and lab coat when handling these buffers. Make sure all materials are nuclease-free when executing extraction.
2. Additional materials: 96~100% ethanol, 100% isopropanol, and 15 ml or 50 ml centrifuge tube.
3. Additional equipment: vacuum source capable of -650 mmHg and vacuum manifold with an adaptor for column tip.
4. Add the indicated volume of isopropanol (100%) into the **CB Binding Buffer**, mix well; add the indicated volume of ethanol (96~100%) into the **Wash Buffer 1** and **Wash Buffer 2**, mix well. Store the buffer at room temperature.
5. Specimen storage and handling:
  - a. To prevent the degradation of cell-free nucleic acids and the release of cellular nucleic acids, whole blood should be stored at 2~8°C for no more than 6 hours. (e.g., in EDTA tube or specialized cfDNA collection tubes) Separate plasma or serum promptly.
  - b. For short-term storage, the specimen can be stored at 2~8°C for up to 24 hours after collection.
  - c. For long-term storage, keep the specimen at -20°C or -80°C for up to 12 months for DNA, or at -80°C for up to 4 weeks for RNA.
  - d. Make sure the cell-free fluid samples are clear. **Centrifuge the samples for 2 mins at 400 xg if the debris are still visible.**
6. All centrifugation steps should be performed at **18,000 xg**.

## ■ General Protocol: For 1~4 ml cell-free fluid

1. Add the indicated volume of **Proteinase K** and cell-free fluid sample into a 15 ml or 50 ml centrifuge tube (Refer to **Table 1/Step 1**), then mix thoroughly.
2. Add the indicated volume of **CL Lysis Buffer** to mixture then mix thoroughly (Refer to **Table 1/Step 2**).
3. Incubate the sample mixture at 56°C for 30 mins. Vortex occasionally during incubation (2~3 times).
4. Add the indicated volume of **CB Binding Buffer** (isopropanol contained) to mixture then mix thoroughly (Refer to **Table 1/Step 4**).

**Table 1. Quick Setup Guide**

| Step 1 | Sample Volume (ml)                               | 1   | 2   | 3   | 4   |
|--------|--|-----|-----|-----|-----|
|        | Proteinase K (μl)                                | 100 | 200 | 300 | 400 |
| Step 2 | CL Lysis Buffer (ml)                             | 0.8 | 1.6 | 2.4 | 3.2 |
| Step 3 | Incubate the sample mixture at 56°C for 30 mins. |     |     |     |     |
| Step 4 | CB Binding Buffer (isopropanol contained, ml)    | 1.8 | 3.6 | 5.4 | 7.2 |

5. Incubate the sample mixture on ice for 5 mins.
6. Connect an **HE Column R** with a **Column Extender** tightly then place the assembly on a vacuum manifold. Transfer the sample mixture carefully to the HE Column R assembly.
7. Pass all the sample mixture through the HE Column R by applying vacuum pressure until the column is empty. Detach the Column Extender.
8. Add 500 µl **Wash Buffer 1** (ethanol contained) to the HE Column R. Apply vacuum pressure until the column is empty.
9. Add 900 µl **Wash Buffer 2** (ethanol contained) to the HE Column R. Apply vacuum pressure until the column is empty.
10. Repeat **step 9**.
11. Place the HE Column R to an **HE Collection Tube**. Centrifuge the HE Column R for an additional 2 mins to dry the HE Column R.
12. Place the HE Column R to an **Elution Tube**, then add 30 µl of **RNase-Free Water** to the membrane center of the HE Column R. Stand the HE Column R for 5 mins.
  - **IMPORTANT!** For effective elution, make sure that RNase-Free Water is dispensed onto the membrane center and absorbed completely.
13. Centrifuge the HE Column R for 2 mins to elute nucleic acid.

For more product information, please visit <https://www.favorgen.com/>  
For technical assistance, please email us at [Technical@favorgen.com](mailto:Technical@favorgen.com)

This product is for research use only. It is not intended to be used for therapeutic or diagnostic purposes.