



FavorPrep™ Gel/PCR Purification HE Mini Kit

■ Kit Contents

Cat. No.	FAGC1030 (4 Preps)	FAGC1033 (50 Preps)	FAGC1034 (100 Preps)
GPP Buffer	7 ml	90 ml	180 ml
Wash Buffer (Concentrate) ▲	1.5 ml	15 ml	30 ml
Elution Buffer	0.5 ml	5 ml	7 ml
HE Columns	4 pcs	50 pcs	50 pcs x 2
HE Collection Tubes	8 pcs	100 pcs	100 pcs x 2
Elution Tubes	4 pcs	50 pcs	100 pcs
User Manual	1	1	1
Preparation of Wash Buffer by adding 96~100% ethanol.			
Volume of ethanol for Wash Buffer ▲	6 ml	60 ml	120 ml


All kit components are shipped at room temperature and should be stored at room temperature between 15~25°C upon receipt.

■ Specification

Format/Principle	Spin column (Silica matrix)
Binding Capacity	≤125 µg DNA/Column
Operation Time	<20 mins for gel extraction <10 mins for PCR clean up
Sample Size	≤500 mg agarose gel ≤300 µl PCR product or enzymatic reaction
DNA Size	50 bp~10 kbp
Recovery	70~95% for gel extraction 90~95% for PCR clean up
Elution Volume	30 µl


■ Procedure Overview

For gel extraction



- Transfer up to 500 mg agarose gel into microcentrifuge tubes and add 2.5x gel volumes of **GPP Buffer**. (If the input gel exceeds 300 mg, transfer the gel and GPP Buffer to a 15 ml centrifuge tube.)
- Incubation at 55°C for 10 mins (invert per 2~3 mins).
- Incubate the sample mixture at room temperature for 2 mins.


For PCR clean up



- Transfer up to 300 µl PCR product into a microcentrifuge tube and add 5x gel volumes of **GPP Buffer**. (If the input PCR product exceeds 200 µl, transfer the PCR product and GPP Buffer to a 15 ml centrifuge tube.)

HE Column

↻ Centrifuge
12,000 xg, 1 min



- Transfer up to 1000 µl of the sample mixture into the **HE Column** for DNA binding.
- Repeat the step if the sample mixture exceeds 1000 µl.

↻ Centrifuge
12,000 xg, 1 min

↻ Centrifuge
12,000 xg, 2 mins



- Add 900 µl **Wash Buffer** (ethanol contained).
- Add 500 µl **ethanol** (96~100%) and dry the column membrane.

↻ Centrifuge
12,000 xg, 2 mins



- Add 30 µl **Elution buffer**.
- Stand the column for 5 mins.
- Obtain purified DNA.

■ Preparation Before Starting

1. Additional materials: 100% Isopropanol and 96~100% ethanol.
2. Add the indicated volume of ethanol (96~100%) into the **Wash Buffer**, mix well, and store at room temperature.
3. Set up a water bath or dry bath at 55°C and preheat the **Elution Buffer** to 55°C for the elution step.
4. All centrifugation steps should be performed at **12,000 xg** at room temperature.

■ Protocol: For extraction of DNA fragments from agarose gel

1. Excise the agarose gel with a clean scalpel. Weigh the gel slice in a 1.5 ml microcentrifuge tube (not provided).
 - Remove the extra agarose gel to minimize the size of the gel slice.
 - The maximum volume of the gel slice is 500 mg. If the input gel exceeds 300 mg, transfer the gel to a 15 ml centrifuge tube (not provided).
2. Add 2.5x gel volumes of **GPP Buffer** to 1x volume of gel (100 mg gel ~100 µl).
3. Incubate at 55°C for 5~10 mins and invert the tube every 2~3 mins until the gel slice is dissolved completely.
 - Make sure that the gel slice has been dissolved completely before proceeding to the next step.
 - If the color of the mixture is orange or violet, add 10 µl 3 M sodium acetate, pH 5.0, and mix. The color of the mixture will turn yellow.
4. Cool down the sample mixture to room temperature for 2 mins.
5. **(Optional)** If the desired DNA fragment is ≤ 200 bp, add 1x gel volume of 100% isopropanol to the sample mixture. Mix thoroughly.
6. Place an **HE Column** in an **HE Collection Tube**, then carefully transfer up to 1000 µl of sample mixture into the HE Column.
 - If the sample mixture is more than 1000 µl, repeat this step for the rest of the sample mixture.
7. Centrifuge for 1 min. Discard the flow-through and place the HE Column in a new HE Collection Tube.
8. Add 900 µl **Wash Buffer** (ethanol contained) to the HE Column. Centrifuge for 1 min, then discard the flow-through.
9. Add 500 µl ethanol (96~100%) to the HE Column. Centrifuge for 2 mins to dry the membrane. Discard flow-through and HE Collection Tube.
 - **Important step!** Don't let the tip of the column touch the flow-through. If in doubt, spin again to remove all residual ethanol.
10. Place the HE Column in an **Elution Tube**, then add 30 µl prewarmed Elution Buffer or ddH₂O (pH 7.5~9.0) directly onto the membrane. Stand the HE Column for 5 mins.
 - **Important step!** For effective elution, ensure that the elution solution is dispensed onto the membrane center and absorbed completely.
11. Centrifuge for 1 min to elute the DNA.

■ Protocol: For purification of PCR products or reaction mixtures

1. Transfer up to 300 µl of PCR product (excluding oil) to a 1.5 ml microcentrifuge tube (not provided) and add 5x volumes of **GPP Buffer**, mix well by vortexing.
 - If PCR product is more than 200 µl, transfer the PCR product to a 15 ml centrifuge tube (not provided).
 - If the color of the mixture is orange or violet, add 10 µl 3 M sodium acetate, pH 5.0, and mix. The color of the mixture will turn yellow.
2. **(Optional)** If the desired DNA fragment is ≤ 200 bp, add 1x PCR product volume of 100% isopropanol to the sample mixture. Mix thoroughly.
3. Place an **HE Column** in an **HE Collection Tube**, then carefully transfer all mixture into the HE Column.
 - If the sample mixture is more than 1000 µl, repeat this step for the rest of the sample mixture.
4. Centrifuge for 1 min. Discard the flow-through and place the HE Column in a new HE Collection Tube.
5. Add 900 µl **Wash Buffer** (ethanol contained) to the HE Column. Centrifuge for 1 min then discard the flow-through.
8. Add 500 µl ethanol (96~100%) to the HE Column. Centrifuge for 2 mins to dry the membrane. Discard flow-through and HE Collection Tube.
 - **Important step!** Don't let the tip of the column touch the flow-through. If in doubt, spin again to remove all residual ethanol.
7. Place the HE Column in an **Elution Tube**, then add 30 µl prewarmed Elution Buffer or ddH₂O (pH 7.5~9.0) directly onto the membrane. Stand the HE Column for 5 mins.
 - **Important step!** For effective elution, ensure that the elution solution is dispensed onto the membrane center and absorbed completely.
8. Centrifuge for 1 min to elute the DNA.

For more product information, please visit <https://www.favorgen.com/>
For technical assistance, please email us at Technical@favorgen.com

This product is for research use only. It is not intended to be used for therapeutic or diagnostic purposes.