



FavorPrep™ Urine DNA Extraction HE Mini Kit

■ Kit Contents

Cat. No.	FAUD1030 (4 Preps)	FAUD1033 (50 Preps)	FAUD1034 (100 Preps)
FAUD Buffer	3 ml	30 ml	60 ml
W1 Buffer (Concentrate) ▲	0.8 ml × 2	12 ml	24 ml
Elution Buffer	0.5 ml	5 ml	7 ml
Lysozyme Solution	55 µl	450 µl	450 µl × 2
Lyticase Solution	50 µl	450 µl	450 µl × 2
Proteinase K (Liquid)	100 µl	1050 µl	1050 µl × 2
HE Columns	4 pcs	50 pcs	50 pcs × 2
HE Collection Tubes	8 pcs	100 pcs	100 pcs × 2
Elution Tubes	4 pcs	50 pcs	100 pcs
User Manual	1	1	1
Preparation of W1 Buffer by adding 96~100% ethanol.			
Volume of Ethanol for W1 Buffer ▲	1 ml	18 ml	36 ml

1. All kit components are shipped at room temperature and should be stored at 15~25°C upon receipt, **except Lysozyme Solution** and **Lyticase Solution**.
2. Store Lysozyme Solution and Lyticase Solution at -20°C upon receipt.

■ Specification

Format/Principle	Spin column (Silica matrix)
Binding Capacity	≤125 µg DNA/Column
Operation Time	<45 mins
Sample Size	300 µl Urine
Elution Volume	30 µl

■ Procedure Overview

Urine Sample



- Add 300 μ l Urine, 8 μ l **Lysozyme Solution**, 8 μ l **Lyticase Solution** into a tube.
- Incubate at 37°C for 10 mins.
- Add 20 μ l **Proteinase K** and 500 μ l **FAUD Buffer** at 60°C for 10 mins.



- Stand at room temperature for 3 mins.
- (Optional) Add RNase A at room temperature for 2 mins.
- Add 150 μ l ethanol (96~100%).

HE Column



↻ Centrifuge
12,000 xg, 1 min

- Transfer the mixture into the **HE Column** for DNA binding.




↻ Centrifuge
12,000 xg, 1 min

- Add 500 μ l **W1 Buffer** (ethanol contained).



↻ Centrifuge
12,000 xg, 2 mins

- Add 500 μ l **ethanol** (96~100%) and dry the column membrane.



↻ Centrifuge
12,000 xg, 1 min

- Add 30 μ l **Elution Buffer**.
- Stand the column for 5 mins.
- Obtain purified genomic DNA.

■ Preparation Before Starting

1. Add the indicated volume of ethanol (96~100%) into the **W1 Buffer**, mix well, and store at room temperature.
2. Additional materials: 96~100% ethanol, RNase A (Cat. No. FARA2093; optional).
3. Set up two water baths or dry baths: one at 37°C for the **Lysozyme** & **Lyticase** incubation step; the other at 60°C for **Proteinase K** incubation and **Elution Buffer** preheating.

■ General Protocol

1. Add 300 µl urine sample, 8 µl **Lysozyme Solution** and 8 µl **Lyticase Solution** into an eppendorf (not provided). Mix thoroughly.
 - For DNA isolation from Gram-positive bacteria, urine sample requires an additional incubation at 95°C for 5 mins before lysozyme and lyticase reaction.
2. Incubate the mixture at 37°C for 10 mins to disrupt cell wall.
3. Add 20 µl **Proteinase K** and 500 µl **FAUD Buffer** into the sample mixture. Mix thoroughly.
 - **DO NOT** add Proteinase K directly into FAUD Buffer.
4. Incubate the mixture at 60°C for 10 mins to lyse the sample. During incubation, vortex the tube every 5 mins.
5. Incubate the sample mixture at room temperature for 3 mins.

(Optional) If RNA-free genomic DNA is required, add 12 µl of RNase A (50 mg/ml; not provided). Mix thoroughly and incubate at room temperature for 2 mins.
6. Add 150 µl ethanol (96~100%) to the sample mixture. Mix gently by pipetting or inverting.
7. Place an **HE Column** in an **HE Collection Tube**, then carefully transfer all mixture into the HE Column.
8. Centrifuge at 12,000 xg for 1 min. Discard the flow-through and place the HE Column in a new HE Collection Tube.
9. Add 500 µl **W1 Buffer** (ethanol contained) to the HE Column. Centrifuge at 12,000 xg for 1 min then discard the flow-through.
10. Add 500 µl ethanol (96~100%) to the HE Column. Centrifuge for 2 mins to dry the membrane directly. Discard flow-through and HE Collection Tube.
11. Place the HE Column in an **Elution Tube**, then add 30 µl prewarmed **Elution Buffer** or ddH₂O (pH 7.5~9.0) directly onto the membrane. Stand the HE Column for 5 mins.
 - **Important step!** For effective elution, ensure that the elution solution is dispensed onto the membrane center and absorbed completely.
12. Centrifuge at 12,000 xg for 1 min to elute the DNA.

For more product information, please visit <https://www.favorgen.com/>
 For technical assistance, please email us at Technical@favorgen.com