



TSL-MINI Incubator

Instruction manual

Catalog No. MINI-INCBASIC



www.thistlescientific.co.uk
enquiries@thistlescientific.co.uk

Version 05A
Revised on: 2023/05/31

Packing List

- CLV-MINI-INCBASIC Incubator x1
- Power Cord x1
- Stainless Steel Shelf x1
- Instruction Manual x1

Signed by:

Date:

Thistle Scientific is liable for all missing or damaged parts / accessories within 7 days after customer received this instrument package. Please contact Thistle Scientific immediately regarding this issue. If no response within such time period from consignee party, that will be consignee party's whole responsibility.

Table of Contents

Packing List.....	2
Warning.....	5
Section 1 Introduction	8
1.1 Overview	8
1.2 Feature	8
1.3 Components Guide	9
Section 2 Technical Specification.....	10
Section 3 Installation Instructions	11
Section 4 Operation Instructions	11
4.1 Control Interface	11
4.1.1 Touch Keypad	12
4.1.2 Color Display	13
4.2 Start the operation	14
4.3 System Setup	15
4.4 Temperature Calibration	16
Section 5 Function Control Software	20
5.1 Installation Instruction	20
5.2 Operation Instruction	20
Section 6 Trouble shooting and Maintenance.....	39
6.1 Trouble Shooting.....	39
6.2 Warning Illustration	40
6.3 Cleaning and Maintenance	42
Section 7 Ordering Information	43

Section 8 Warranty	44
---------------------------------	-----------

Warning

This equipment has been tested and verified to comply with safety limits. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment may generate, use, and radiate radio frequency energy, and if not installed and used in accordance with the instruction manual, may cause harmful interference with radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at their expense. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. It is strongly recommended that the user carefully read the following prints before operating any equipment.

1. Read and follow the manual instructions carefully.
2. Do not alter the equipment. Failure to follow these directions could result in personal and/or laboratory hazards, as well as invalidate equipment warranty.
3. Use a properly grounded electrical outlet with correct voltage and current handling capacity.
4. Disconnect from power supply before maintenance and servicing. Refer servicing to qualified personnel.
5. Never use this instrument series without having the safety cover correctly in position.
6. Do not use the unit if there is any sign of damage to the external tank or cover. Replace damaged parts.
7. Do not use in the presence of flammable or combustible material; fire or explosion may result. This device contains components which may ignite such materials.
8. Refer maintenance and servicing to qualified personnel.
9. Ensure that the system is connected to electrical service according to local and national electrical codes. Failure to make a proper connection may create fire or shock hazard.
10. Use appropriate materials and operate correctly to avoid possible hazards of explosion, implosion or release of toxic or flammable gases arising from overheated materials.
11. The unit shall be operated only by qualified personnel.

Safety Information

Use high level of precaution against any electrical device. Before connecting the electrical supply, check to see if the supply voltage is within the range stated at the rating label, and see to it that the device be seated firmly. Place the unit in a safe and dry location; it must NOT touch the surrounding. Follow the safety precautions for chemicals / dangerous materials. If needed, please contact qualified service representative or enquiries@thistlescientific.co.uk

Environmental Conditions

Ensure the instrument is installed and operated strictly under the following conditions:

1. Indoor use only
2. $\leq 95\%$ RH
3. 75 kPa – 106 kPa
4. Altitude must not exceed 2000 meters
5. Ambient to 40°C operating temperature
6. Pollution degree: 2
7. Mains supply voltage fluctuations up to $\pm 10\%$ of the normal voltage

Avoiding Electrical Shock

Follow the guidelines below to ensure safe operation of the unit.

The MS MINI Incubator has been designed to utilize shielded wires thus minimizing any potential shock hazard to the user. Thistle Scientific recommends against the use of unshielded wires.

To avoid electrical shock:

1. In the event of solution spilling on the instrument, it must be dried out for at least 2 hours and restored to NORMAL CONDITION before each operation.
2. Never connect or disconnect wires loading from the power jacks when the red indicator light of power switch is on.
3. WAIT at least 5 seconds after stopping a run before handling output leads or any connected apparatus.
4. ALWAYS make sure that your hands, work area, and instruments are **clean** and **dry** before making any connections or operating the power supply.
5. ONLY connect the power cord to a properly grounded AC outlet.

Avoiding Damage to the Instrument

1. Do not attempt to operate the device if damage is suspected.
2. Protect this unit from physical damage, corrosive agents and extreme temperatures (direct sunlight, etc.).
3. For proper ventilation and safety concerns, keep at least 10 cm of space behind the instrument, and at least 5 cm of space on each side.
4. Use high level of precaution against the damages on the unit.
5. Do not operate the unit out of environmental conditions addressed above.
6. Prior to applying any cleaning or decontamination methods other than manufacturer's recommendation, users should check with the manufacturer's instruction to see if the proposed method will damage the equipment.

Equipment Operation

Follow the guidelines below to ensure safe operation of the unit:

1. NEVER access dangerous chemicals or other materials to prevent possible hazard of explosion and damage.
2. Do not operate the unit without lids or covers to prevent possible hazards.
3. A temporary conductivity caused by condensation might occur even though this series is rated Pollution Degree 2 in accordance with IEC 664.

Symbol

The symbol used on MS Mini Incubator is explained below.



Indicates an area where a potential shock hazard may exist. Consult the manual to avoid possible personal injury or instrument damage.



Indicates disposal instruction.

DO NOT throw this unit into a municipal trash bin when this unit has reached the end of its lifetime. To ensure utmost protection of the global environment and minimize pollution, please recycle this unit.

Section 1 Introduction

1.1 Overview

Thistle Scientific's Mini Incubator is designed for personal use and small laboratories, saving much space. The unit features a broad temperature range to meet a variety of microbiology or hematology applications.

It is compact and economically priced, yet offers unique features not typically found in other incubators. The housing is all corrosion resistant metal, so is the door frame. The glass window of the door offers full visibility to the interior. One stainless steel shelf is included, and can be adjusted to two different levels – additional shelves can be purchased to increase storage capacity. The color touch panel presents vivid and intuitive user interface and the backlit LED ensures enough brightness even in dark environment.


The heating structure incorporated in MS Mini Incubator performs evenly heating and precise temperature control. The microprocessor ensures temperature stability within the chamber. A thermometer port on the top of the unit allows the user to monitor temperatures. In addition, a designed spare port is on the rear of the chamber which allows for an electrical cord to be routed out if needed. Optional Function Control software allows Mini Incubator to perform real-time data recording through a PC.

1.2 Feature


- Ideal for microbiology or hematology applications
- Corrosion resistant metal chamber
- Door with large viewing area
- Backlit color touched panel
- One stainless steel shelf is included
- Functional control software available

1.3 Components Guide

Front view



Rear view



Section 2 Technical Specification

Display	LCD
Temperature Range	Ambient +5°C to 70°C
Capacity	17L
Exterior Dimension (WxLxH)	310x306x380 mm
Interior Dimension (WxLxH)	261x255x255 mm
Weight	Approx. 12.65 kg
Power	100-240V~, 50/60 Hz, 2A
Material	Metal
PC Connection	USB


Section 3 Installation Instructions

The MS Mini Incubator is a pre-installed instrument. Place the unit on a sturdy, level, safe and dry place; it is ready for operation.

Note: DO NOT hold the doorknob while lifting up the machine. Please use both hands from the bottom to move the machine.

Section 4 Operation Instructions

4.1 Control Interface



(1) Touch Keypad: Parameters such as temp. or timer can be set through touch screen interface. Detailed descriptions are mentioned below.


(2) Color Display: It shows present value, setting value, and timer.

4.1.1 Touch Keypad

Icon	Description	Function
	Selection Key (Upward)	- Press  to move upward.
	Value-changing Key (Increase)	- Press  to increase the value. Long press the key to speed up changing.
	Selection Key (Leftward)	- Press  to move leftward.
	Selection Key (Downward)	- Press  to move downward.
	Value-changing Key (Decrease)	- Press  to decrease the value. Long press the key to speed up changing.
	Selection Key (Rightward)	- Press  to move rightward.
	Enter Key	- Press  to enter selected value.
	Confirmation Key	- Press  to confirm the adjusted value.
	System Setup Key	- Press  to enter system setup page to set temperature offset value and touch keypad sensitivity. Password is needed to enter the page.
	Start / Stop Key	- Press  to start or stop the operation.

4.1.2 Color Display

The system is equipped with a color display. Parameters such as temperature or timer can be viewed here.



(1) PV: Shows the present temperature status in the chamber.

(2) SV: Shows the set temperature value user designed.

(3) Timer: Use the touch keypad to set up time control. **The timer starts count down when present value reaches the set value.** If time is set as “00:00:00” it means **continuous**; icon  at the right will show in blue.

4.2 Start the operation

Step 1 Prepare the samples according to lab protocol. Place the sample into the oven chamber. Make sure the door closes all the way.

Step 2 To set the temperature or time value, press   to move between the 2 parameter. The selected one will display in blue. Press , the figure will start flashing, use   to set the value. Press  again to confirm.

Note 1: If users do not press  to confirm for more than 3 seconds, system will automatically confirm.

For Timer:

1. Timer is for temperature control only.
2. Timer activates when the target temperature has been reached for the first time. If a value is set for time, it will count down from that value to 0, with alarm sounding for 1 minute.
3. If no value is set for timer (ie. 0), timer will start count up to 9999; in which case the time will not increase anymore, but temperature control will remain.
4. If user reset the target temperature value after timer starts counting, the timer will NOT reset the countdown.


Step 3 To start temperature control, press the  button. Once temperature control starts, the button will light in red. 

Tip for maintaining certain humidity during operation:

To prevent samples from drying out during incubation, user could put a small glass of water at the corner of inner chamber before starting operation.

4.3 System Setup


Step 1 Press  to enter system setup page. Input the password and press  button. (Password: 1111)



Step 2 Press   to move between the 2 columns. The selected one will display in yellow. Press  to start changing the value. Press  again to confirm. If users do not press  to confirm for more than 3 seconds, system will automatically confirm.

Temperature Offset: Ranging from -10.0 to +10.0

Touch Panel Sensitivity: Ranging from 1 to 100. The larger the value is, the more sensitive the touch keypad become.




Step 3 Press  to return to home page.

4.4 Temperature Calibration

The system allows users to calibrate the temperature parameter to their environment standard. Normally there is no need to perform this action unless there is a temperature discrepancy between the factory environment and your working environment. In such cases, follow the guidelines below to perform user temperature calibration.

Step 1 Place a certified thermometer through thermometer hole, which is on the top of the oven chamber. Make sure the front door and the behind spare port is closed during the entire calibration process.


Note: It is suggested to place thermometer in the middle of chamber.



Step 2 Switch on the main power. When the screen is blue, press . Don't release the button until figure (on the upper left corner) is between 4000 and 5000. Once the figure reaches between 4000 and 5000, release , and the oven is now under Calibration Mode.




Step 3 Under Calibration Mode, Low Point Calibration is required to perform first. The default low point is 40°C, user only needs to input the heating time. It is suggested to input at least 50 minutes to reach temperature balance. When finish setting, press  button to start heating.



Note: Use   button to select Hour/Minute/Second and   button to adjust the value.

Step 4 Once present value reaches 40°C, the timer starts count down. When countdown is over, a beep-beep sound will alert. Press  button, a blank column will appear. Input the temperature value you have measured, and press  button to confirm.




Step 5 Please proceed to High Point Calibration. The default high point is 50°C, user only needs to input the heating time. It is suggested to input at least 50 minutes to reach temperature balance. When finish setting, press  button to start heating.



Note: Use   button to select Hour/Minute/Second and   button to adjust the value.

Step 6 Once present value reaches 50°C, the timer starts count down.

When countdown is over, a beep-beep sound will alert. Press  button, a blank column will appear. Input the temperature value you have measured, and press  button to confirm.



Step 7 System will automatically adjust the temperature. Wait until both "Target" and "Temp." show 50°C and the calibration process is complete. After finishing calibration, **turn off the machine.**



Section 5 Function Control Software

NOTE: All rights reserved by Thistle Scientific co., Ltd.

5.1 Installation Instruction

1. Insert the CD into CD ROM and press the setup.exe in the Installer Folder for installation.
2. Follow up the instructions shown on the computer display screen to complete the installation.


5.2 Operation Instruction

After the installation completes, user can find a “MS” icon (as the image below) on user’s desktop. Double-click the icon to start the remote control operation.




***Notice:**

User will need to use the USB cable to connect the device with user's computer before starting the remote control operation. If the device is not connected or connected wrongly, the software will pop up a warning dialogue as the picture below. The dialogue will disappear once the device is connected. Click “OK” to close the dialogue.




Remote Control Interface Introduction

Start MO-MINI software program and then the below screen will be shown.




Error System Setup



A.1 Option

Option	Set up the file name format and record period.
--------	--



It is suggested select this check box to automatically save file every time users start running the machine. If not, system will pop up a file-saving dialogue.

Select this check box if it is needed to include device number into the file name.

Appoint desired string of words before/after the file name. You can preview the result at Example below.


The **default value** of how frequent the operation data is recorded. This value will apply on **ALL** machines and **EVERY** operation.

A.2 Rescan Device

Rescan Device	If the device has once disconnected, press "ReScan Device" to help the software find the reconnected device.
---------------	--


A.3 Language

Language	Set up system language.
----------	-------------------------



A.4 About


About	Display the software version and patent information.
-------	--



A.5 Device No.


1	System can control more than one Mini Incubator at the same time; click the blank to switch.
---	--

B Main Column:




B.1 Link Status

Indication whether the Mini Incubator is linked with computer or not.



Disconnected



Connected

B.2 RUN/STOP


When user finishes all the setting, Click "RUN" to start running the device and press the same button to stop operation.

B.3 Exit

Click to leave the remote control program.


(C) Setting & Record Column

There are 3 different tabs in this column: **Chart** (C.1), **Setting** (C.2), and **History** (C.3). We'll elaborate the function of each tab respectively in the paragraphs below.




③ .1 Chart:

The software will start recording the data as soon as Mini Incubator starts operating (whether start by software control or by button on machine.)



c.1-a

This area helps user distinguish the curves of different parameters in the running chart.




c.1-b

The program will display user's set values and the current condition in this section. User can compare the values with the curves in chart.


Note:

In this column, all the values are unchangeable. If user attempts to adjust the set values, user has to stop the device and enter the setting tab to edit the parameters.




Temp. SV	35.0
Temp. PV	34.9
Machine Time	00:05:00
Record Time	264

Compare this chart with the column above: at 264th second the temperature reaches 34.9°C



C.1-c

To set up how frequent the operation data is recorded. Press the button; it can be set either in minute or in second.

A dialog box titled "REC. Period" with a text input field containing the number "1" and a button labeled "Second".

*Note:

The frequency of data-recording here is only for **current** running. Next time when user starts the program, system will apply the default value in **Option** rather than the value here. Please refer to **A.1 Option** for setting up default value.


C.1-d

In Chart Option, user is able to edit the scale of chart.

A dialog box titled "Chart Option" with four rows of controls. Each row has a label, a text input field, and a blue dot with a line pointing to a description. The rows are: "Min Temp." with input "0" and description "Min Temp: The minimum of temperature scale in the chart"; "Max Temp." with input "70" and description "Max Temp: The maximum of temperature scale in the chart"; "Min Time" with input "0" and description "Min Time: The minimum of time scale in the chart"; and "Max Time" with input "0" and description "Max Time: The maximum of time scale in the chart".


C.1-e

Select either "Unit" or "All in standard" option for the time scale of the chart.

A dialog box titled "Time" with two radio button options: "Unit" (selected) and "All in standard".

☒ Unit	Select "Unit" and press "Apply" to show the part of trend in the range of time scale user set.
☐ All in standard	Select "All in standard" and press "Apply" to show the whole operation timespan in 1 sheet of the chart.

Image 1 (Unit)



Set desired time range. The max value is total record time.


Use the scroll bar to view the rest of the data.

The time scale is range from the user-defined minimum to maximum value. In this case, it is from 0-219 (unit: 2 sec.).

Note:

If the distance between minimum and maximum time is less than 20 sec/min, no scroll bar will display.

Image 2 (All in Standard)



The time range cannot be set.

The total recorded data will be concentrated into 1 sheet. In this case, it is the whole 438 seconds operation data.

C.1-f Reset Button

A rectangular button with a light gray background and a thin black border. The word "Reset" is centered in a black, sans-serif font.

Press "Reset" to restore to default settings.

C.1-f Apply Button

A rectangular button with a light gray background and a thin black border. The word "Apply" is centered in a black, sans-serif font.

Press "Apply" to confirm the new setting.

C.1-h Print Chart Button

A rectangular button with a light gray background and a thin black border. The words "Print Chart" are centered in a black, sans-serif font.

Press the button to print current-viewed chart.

c.2 Setting:


In this tab you can set desired temperature value for desired time.

The screenshot shows a software interface with three tabs: 'Chart', 'Setting', and 'History'. The 'Setting' tab is active. Inside the 'Setting' tab, there are two main sections. The first section is for temperature settings, labeled 'Temp. SV:', with a text input field containing '25.0' and a unit label '°C (0~70)'. The second section is for time settings, labeled 'Time:', with three stacked text input fields, each containing '0', and corresponding unit labels 'hour', 'minute', and 'second'. At the bottom of the 'Setting' tab, there are two buttons: 'Download' and 'Upload'.


Temp.SV	Set up desired temperature value.
Time	The countdown timer which starts count down once the present temperature value reaches the set value.
	Press Download button to receive setting data from machine.
	Press Upload button to launch above setting data to the machine.

C.3 History:

History page look very similar to Chart page. In History, you can view historic record data.



C.3-a Open File Button



Click to open desired historic record data (.log file)

C.3-b Chart Option

In Chart Option, user can edit the scale of chart.


A dialog box titled "Chart Option" with four input fields and their corresponding descriptions. The fields are: "Min Temp." with value "-10", "Max Temp." with value "100", "Min Time" with value "0", and "Max Time" with value "46". Each field has a blue dot to its right, connected by a line to its description. The descriptions are: "Min Temp: The minimum of temperature scale in the chart", "Max Temp: The maximum of temperature scale in the chart", "*Min Time: The minimum of time scale in the chart", and "*Max Time: The maximum of time scale in the chart".

Field	Value	Description
Min Temp.	-10	Min Temp: The minimum of temperature scale in the chart
Max Temp.	100	Max Temp: The maximum of temperature scale in the chart
Min Time	0	*Min Time: The minimum of time scale in the chart
Max Time	46	*Max Time: The maximum of time scale in the chart

*The distance between the minimum and maximum time scale should be more than 20 seconds/minutes, otherwise the resolution of the curves will be compromised.

C.3-c Time Scale Option Column

Select one of the option and press "Apply" to view the curves with different time scales.

A vertical column of three radio buttons under the heading "Time". The options are "Unit", "All in standard", and "All in scale". The "All in scale" option is selected, indicated by a blue dot.

Time


☐ Unit

☐ All in standard

☒ All in scale

Unit	Select "Unit" and press "Apply" to show the part of trend in the time scale user defined.
All in standard	Select "All in standard" and press "Apply" to show the whole operation data in 1 sheet of the chart.
All in scale	Select "All in scale" and press "Apply" to show the part of trend in the time scale user defined in 1 sheet.

Image 1 (Unit)



Set desired time range. The max value is total record time.


Use the scroll bar to view the rest of the data.

The time scale is range from the user-defined minimum to maximum value. In this case, it is from 0 to 25.

Note:

If the distance between minimum and maximum time is less than 20 sec/min, no scroll bar will display.


Image 2 (All in standard)



The time range cannot be set.

The total recorded data will be concentrated into 1 sheet. In this case, it is the whole 142 seconds operation data.

Image 3 (All in scale)



Set desired
time range
concentrated
in 1 sheet

The time scale is range from the user-defined minimum to maximum value, and data in the range will display in 1 sheet. In this case, it is from 0 to 25.

Section 6 Troubleshooting and Maintenance

6.1 Troubleshooting

Problem	Possible Causes and Solution
1. The system does not heat up	a. Door is not tightly closed. b. Temperature control is not activated.
2. The system cannot be powered up	Broken fuse. Replace the fuse by following the procedure on the next page
3. Temperature discrepancy between the display and thermometer	There might be some temperature discrepancy between factory setting and your environment. Adjust the temperature to your environment standard by following the Section 4.4 Temperature Calibration.

6.2 Warning Illustration

Alert Information	Description	Performance	Solution
 Finish!	The unit finishes temperature control. (Timer finishes count down.)	Warning pops up and beep-beep sound alerts.	Press  or  to go back to home page.
 Warning ! E0002	During operation, the heating device function abnormally hence cause great deviation between set temperature and measured temperature. E0001: indicates PV is 5°C above SV. E0002: (PV has reached SV) indicates PV has been 1°C lower than SV for more than 40 minutes. E0003: indicates PV is over 85°C. E0004: (PV has not reached SV) indicates PV has been 1°C lower than SV for more than 40 minutes.	Warning sign pops up with code and beep-beep sound alerts. System stops heating.	Press  or  to go back to home page.

--	--	--	--

***PV**: Present Value **SV**: Setting Value

Replacing the Fuse

For additional fuses, contact Thistle Scientific.


To replace the fuse:

1. Turn off the main power switch at the rear of image system and detach the power cord from the rear of image system.
2. Open the fuse compartment located inside the Power Entry Module by inserting a small flathead screwdriver into the slot below the ON/OFF switch. Turn the screwdriver to gently pry open the fuse compartment.

Note: the fuse compartment will not open with the power cord in place.

3. Pull the fuse holder out of the compartment and inspect the fuse. If the fuse is burned or there is a break in the fuse element, replace the fuse with an identical type of fuse (**2A/250V~**) as provided in the fuse holder (see figure below)

4. Place the fuse holder back into the compart
5. Snap the cover close.



6.3 Cleaning and Maintenance

The oven interior can be cleaned with a mild detergent or water. To clean the oven, please make sure the power cord is disconnected.

The system should NEVER come into contact with the following reagents, which could cause irreversible damage:

Acetone, Phenol, Chloroform, Carbon Tetrachloride, Ethanol, Methanol, Isopropyl alcohol, Alkalies

Section 7 Ordering Information

Cat. No.	Description
MINI-INCBASIC	Thistle Scientific Mini Incubator with one stainless steel shelf

Accessories

INCBASIC-SHELF	Stainless Steel shelf for MINI-INCBASIC. 250 x 240mm
INCBASIC-SW	Functional control software for MINI-INCBASIC, including USB cable

Section 8 Warranty

Thistle Scientific warrants apparatus of its manufacture against defects in materials and workmanship, under normal service, for **one year from the shipping date to purchaser**. This warranty excludes damages resulting from shipping, misuse, carelessness, or neglect. Consumable parts are not covered by our warranty. Thistle Scientific's liability under the warranty is limited to the receipt of reasonable proof by the customer that the defect is embraced within the terms of the warranty. All claims made under this warranty must be presented to Thistle Scientific within one year following the date of delivery of the product to the customer.

Contact information:

Address:

Thistle Scientific Ltd.
Unit 41, Somers Road Industrial Estate
Rugby
Warwickshire
UNITED KINGDOM
CV22 7DH

Tel: +44 (0)1788 565300

Email: enquiries@thistlescientific.co.uk

Or support@thistlescientific.co.uk