



# Hybrimmune™ - Hybridoma

**HYBRID YIELDS GENERATED BY ELECTROFUSION ARE 10-FOLD GREATER THAN PEG**

Experiment Number	Antigen Specific Clones	
	E-Fusion	PEG
1	20	0
2	10	0
3	400	24
4	151	21
<b>Mean</b>	<b>145</b>	<b>11</b>

*Table 1: Four different transgenic mice expressing Abs to human Ag were used to compare the efficiency of E-fusion to PEG fusion.*

*Source: Data courtesy of M. Coccia, PhD, Platform Development Group, Medarex Inc, Milpitas, CA)*

**Fusion Methods:** PEG fusion was performed with standard protocols. For E-Fusion, mouse and SP2/O cells were washed twice in Cytofusion medium then mixed in the fusion chamber and tri-phasic pulse applied. Cells were recovered after 30 min and cultured in 96 well plates at 5000 cells/ml. Antigen-specific clones were counted using ELISA or HTRF, normalized to 100M cells.

# Production System

## FEATURES & BENEFITS

**Hybrimmune Advantage** – Hybridoma production efficiency and cell viability are enhanced by specialized waveforms. The patented Ramp-K™ feature enhances cell compression, resulting in high fusion rates and excellent cell viability.

**Non-Uniform Waveform** – Provides rapid cell alignment and compression for increased fusion.

**Scale-up** – Direct scale-up from 2 ml to 9 ml to large-volume hybridoma production in the Hybrimmune system.

**Programmable** – Easy programmable user-friendly windows based software. Data logs are stored and retrieved easily.

## ELECTROFUSION APPLICATIONS

- Monoclonal Antibody Production
- Dendritic – Tumor Cell Fusions

## WAVEFORM GENERATOR SPECIFICATIONS

The Hybrimmune™ Waveform Generator is programmed using the Application Software.

The following parameters are available:

<b>Pulse Function</b>	Constant, linear, non-linear
<b>Pulse Amplitude</b>	100-1000 V
<b>Pulse Width Range</b>	20 – 1000 ms
<b>AC Start Peak Range</b>	5-75 V
<b>AC Stop Peak Range</b>	5-75 V
<b>AC Frequency</b>	0.2 to 2.0 MHz
<b>AC Duration</b>	0 to 60 sec

## FUSION CHAMBER SPECIFICATIONS

Both the optimization and production chambers have been engineered to have identical electrical characteristics to facilitate direct scale-up to production, once pulse parameters have been optimized. In addition, the small chamber has a transparent bottom to permit visualization of the cell alignment by inverted or regular microscope.

Parameter	Optimization Chamber	Production Chamber
<b>Volume</b>	2 ml	9 ml
<b>Outer ID</b>	45.72 mm	45.72 mm
<b>Inner OD</b>	38.10 mm	38.10 mm
<b>Gap</b>	3.81 mm	3.81 mm
<b>Well Height</b>	5 mm	18 mm
<b>Inner/Outer Radius</b>	0.8333	0.8333

For reuse, the fusion chamber can be cleaned with NaOH, sterilization by EtOH, or Spor-Klenz® for spores and mycoplasma.

## License Requirements

The use of the Hybrimmune as a commercial and therapeutic system requires a license from Collectis. Please contact BTX for more information.

## ORDERING INFORMATION

Catalog #	Description
<b>47-0300N</b>	Hybrimmune™ Electrofusion System <b>Includes:</b> Hybrimmune waveform generator, 2 ml and 9 ml coaxial chambers, BTXpress Cytofusio <sup>®</sup> Medium C, user interface software, cables and manual. Requires Windows based laptop or PC (not included).
<b>Accessories:</b>	
<b>47-0301</b>	User-Interface Application Software
<b>47-0030</b>	2 ml Optimization Chamber
<b>47-0020</b>	9 ml Production Chamber
<b>47-0001</b>	BTXpress Cytofusio <sup>®</sup> Medium C, 500 ml



Australian Distributors:  
fisher biotec australia  
free call: 1800 066 077  
email: [info@fisherbiotec.com](mailto:info@fisherbiotec.com)  
web: [www.fisherbiotec.com](http://www.fisherbiotec.com)