


# HOW TO SAFEGUARD YOUR CELL CULTURE FROM CONTAMINATION

## MOST COMMON CONTAMINANTS



## CONTAMINATION OF YOUR CELL CULTURE CAN BE THE RESULT OF ANY OF THESE CONDITIONS:

- Contaminated glass or plastic ware
- Contaminated media
- Contaminated cell line
- Improper technique
- Incubator & its components
- Airborne contaminants
- Lab environment & instruments
- Operators themselves!




CO<sub>2</sub> Incubator needs to be cleaned with a disinfectant or detergent first and then rinsed with sterile water followed by 70% alcohol.




Avoid chlorine-based detergents/disinfectants.



Change water (sterile water recommended) in water pan periodically.



Allow autoclaved water to cool before adding to water pan.




Change ULPA filter at least once a year.



Do not keep CO<sub>2</sub> Incubator door open for too long.



Avoid talking, laughing, coughing, when opening CO<sub>2</sub> Incubator door.



Always wear gloves.



Keep air clean: no eating, drinking, smoking.



Always wear clean lab coat.



Work surface of Biological Safety Cabinet / Laminar Flow Cabinets must be cleaned before and after use.



Avoid leaving medium and culture container on contaminated surfaces.

# ESCO

WORLD CLASS. WORLDWIDE.



Australian distributors:  
Fisher Biotec Australia  
free call: 1800 066 077  
email: [info@fisherbiotec.com](mailto:info@fisherbiotec.com)  
web: [www.fisherbiotec.com](http://www.fisherbiotec.com)