



## FavorPrep™ Fungi/ Yeast Genomic DNA Extraction Mini Kit

Cat. No.: FAFYG 000  
FAFYG 001  
FAFYG 001-1  
**(For Research Use Only)**

### Kit Contents:

Cat. No:	FAFYG 000 (4 preps_sample)	FAFYG 001 (50 preps)	FAFYG 001-1 (100 preps)
Beads Tube	4 pcs	50 pcs	100 pcs
FA Buffer	5 ml	60 ml	120 ml
FB Buffer	2.7 ml	32 ml	65 ml
TG1 Buffer	2 ml	22 ml	45 ml
TG2 Buffer	2 ml	15 ml	30 ml
W1 Buffer <sup>a</sup> (concentrate)	1.3 ml	22 ml	44 ml
Wash Buffer <sup>b</sup> (concentrate)	1 ml	10 ml	20 ml
Elution Buffer	0.5 ml	7 ml	15 ml
Lyticase solution <sup>■</sup>	250 µl	550 µl x 5	550 µl x 10
Proteinase K <sup>c</sup>	1 mg	11 mg	11 mg x 2
TG Mini Column	4 pcs	10 pcs x 5	10 pcs x 10
Collection Tube	8 pcs	100 pcs	100 pcs x 2
Elution Tube	4 pcs	50 pcs	100 pcs

Preparation of W1 Buffer and Wash Buffer by adding ethanol (96 ~ 100%) and **Store at RT.**

Ethanol volume for W1 Buffer <sup>a</sup>	0.5 ml	8 ml	16 ml
Ethanol volume for Wash Buffer <sup>b</sup>	4 ml	40 ml	80 ml

Preparation of Proteinase K solution by adding ddH<sub>2</sub>O and **Store the solution at 4 °C.**

ddH <sub>2</sub> O volume for Proteinase K <sup>c</sup>	0.1 ml	1.1 ml	1.1 ml x 2
---	--------	--------	------------

### Specification:

Principle:	mini spin column (silica matrix)
Sample size:	1~ 5 x10 <sup>6</sup>
Operation time:	~ 60 minutes
Binding capacity:	60 µg/ column
Column applicability:	centrifugation and vaccum

### Additional requirement to be provided by user

1. Microcentrifuge capable of speed at ~18,000 x g
2. 1.5 ml microcentrifuge tube
3. 96~100 % ethanol
4. Vortex
5. Heating block or water bath

### Important Notes:

1. Buffers provided in this system contain irritants. Wear gloves and lab coat when handling these buffers.
2. Store the Lyticase solution **at -20 °C on arrival.**
3. **Caution: Lyticase solution and FB Buffer containing 14 mM of β-mercaptoethanol is hazardous to human health. perform the procedures involving Lyticase solution and FB Buffer in a chemical fume hood.**
4. Add sterilized ddH<sub>2</sub>O to a Proteinase K tube to make a **10 mg/ml** stock solution. Vortex and make sure that Proteinase K has been completely dissolved. **Store the stock solution at 4 °C.**
5. Add required volume of ethanol (96- 100 %) to W1 Buffer and Wash Buffer when first open. **Store the solution at room temperature.**
6. Prepare a heating block or a water baths to 37 °C for step 4 and to 55 °C for step 9 before operation.



Australian distributors:  
Fisher Biotec Australia  
free call: 1800 066 077  
email: [info@fisherbiotec.com](mailto:info@fisherbiotec.com)  
web: [www.fisherbiotec.com](http://www.fisherbiotec.com)

## General Protocol:

### Please Read Important Notes Before Starting Following Steps.

1. Transfer  $1\sim 5 \times 10^6$  of cultures (fungal/ yeast cells) to a 1.5 ml microcentrifuge tube. (not provided)
2. Add 1 ml of **FA Buffer** to the cells and resuspend the cells by pipetting.
3. Descend the cells by centrifuging at 5,000 x g for 2 min and discard the supernatant completely.
4. Resuspend the cells in **550  $\mu$ l of FB buffer and add 50  $\mu$ l of lyticase solution, mix well by vortexing**. Incubate the sample at 37 °C for 30 min.  
**Caution: Lyticase solution and FB Buffer containing 14 mM of  $\beta$ -mercaptoethanol is hazardous to human health. perform the procedures involving Lyticase solution and FB Buffer in a chemical fume hood.**
5. **(Optional)** If RNA-free genomic DNA is required, add 8  $\mu$ l of 50 mg/ml **RNase A** (not provided) and incubate for 2 min at room temperature.
6. Descend the cells by centrifuging at 5,000 x g for 10 min. Remove the supernatant completely.
7. Add 350  $\mu$ l **TG1 Buffer** and mix well by pipetting. Transfer the sample mixture to a **bead tube**. (provided)
8. **Mix well by Plus-vortexing for 5 minutes.**
9. Add 20  $\mu$ l of **Proteinase K** (10 mg/ml) **and mix well by vortexing**. Incubate at 55 °C for 15 min, vortex 30 seconds for every 5 minutes incubation.
10. Descend the cells by centrifuging at 5,000 x g for 1 min and transfer 200  $\mu$ l of supernatant to a new 1.5 ml microcentrifuge tube. (not provided)
11. Add 200  $\mu$ l of **TG2 Buffer** and mix well by pipetting.
12. Add 200  $\mu$ l of **ethanol** (96-100%) and mix well by pulse-vortexing for 10 seconds.
13. Place a **TG Mini Column** in **Collection Tube**. Transfer the sample mixture (including any precipitate) carefully to **TG Mini Column**. Centrifuge at 11,000 x g for 30 second **then place the TG Mini Column to a new Collection Tube**.
14. Add 400  $\mu$ l of **W1 Buffer** to the **TG Mini Column**. Centrifuge at 11,000 x g for 30 seconds and discard the flow-through. Place the TG Mini Column back to the Collection Tube.  
**--Make sure ethanol has been added into W1 Buffer when first use.**
15. Add 750  $\mu$ l of **Wash Buffer** to the **TG Mini Column**. Centrifuge at 11,000 x g for 30 seconds and discard the flow-through. Place the TG Mini Column back to the Collection Tube.  
**--Make sure ethanol has been added into Wash Buffer when first use.**
16. Centrifuge at full speed ( $\sim 18,000$  x g) for an additional 3 min to dry the column.  
**Important Step!** This step will avoid the residual liquid to inhibit subsequent enzymatic reactions.
17. Place the **TG Mini Column** to a **Elution Tube**.
18. Add 50  $\sim$ 100  $\mu$ l of **Elution Buffer or ddH<sub>2</sub>O** to the membrane center of the **TG Mini Column**. Stand **TG Mini Column** for 3 min.  
**Important Step!** For effective elution, make sure that the elution solution is dispensed onto the membrane center and is absorbed completely.
19. Centrifuge at full speed ( $\sim 18,000$  x g) for 1 min to elute total DNA.
20. Store total DNA at 4°C or -20°C.