

Table 2. Catalog Numbers and Spectral Properties

Label/Dye	Ex/Em (nm)	Labeling size/Catalog number			
		5-20 ug	20-50 ug	50-100 ug	3x(5-20 ug)
CF@350	347/448	92270	92250	92230	---
CF@405S	404/431	92271	92251	92231	---
CF@405M	408/452	92272	92252	92232	---
CF@405L	395/545	92303	92304	92305	92445
CF@430	426/498	92316	92317	92318	---
CF@440	440/515	92319	92320	92321	---
CF@450	450/538	92322	92323	92324	---
CF@488A	490/515	92273	92253	92233	92446
CF@514	516/548	92331	92332	92333	---
CF@532	527/558	92289	92290	92291	---
CF@543	541/560	92287	92267	92247	---
CF@555	555/565	92274	92254	92234	---
CF@568	562/583	92275	92255	92235	92447
CF@570	568/591	92334	92335	92336	---
CF@583	583/606	92337	92338	92339	---
CF@583R	586/609	92455	92456	92457	---
CF@594	593/614	92276	92256	92236	92448
CF@594ST	593/620	92458	92459	92460	---
CF@597R	597/619	92461	92462	92463	---
CF@633	630/650	92277	92257	92237	---
CF@640R	642/662	92278	92258	92245	---
CF@647	650/665	92279	92259	92238	92449
CF@660C	667/685	92280	92260	92239	---
CF@660R	663/682	92281	92261	92243	---
CF@680	681/698	92282	92262	92240	---
CF@680R	680/701	92283	92263	92246	---
CF@700	695/720	92425	92426	92427	---
CF@750	755/777	92284	92264	92241	---
CF@770	770/797	92285	92265	92242	---
CF@790	784/806	92288	92268	92248	---
CF@800	797/816	92428	92429	92430	---
CF@820	822/835	92431	92432	92433	---
Biotin	N/A	92286	92266	92244	92444
DNP	N/A	92325	92326	92327	---
DIG	N/A	92328	92329	92330	92450
FITC	494/518	92294	92295	92296	---
Cyanine 555	555/565	92412	92413	92414	---
Cyanine 647	650/665	92416	92417	92418	---

Related Products

Catalog number	Product
22004	Ultrafiltration vial, 10K MWCO (pack of 5)
22018	Ultrafiltration vial, 3K MWCO (pack of 5)
30071	AccuOrange™ Protein Quantitation Kit
23012	TrueBlack® IF Background Suppressor System (Permeabilizing)
23013	TrueBlack® WB Blocking Buffer Kit
23007	TrueBlack® Lipofuscin Autofluorescence Quencher
23014	TrueBlack® Plus Lipofuscin Autofluorescence Quencher
40083	NucSpot® 470 Green Nuclear Counterstain
40061	RedDot™2 Far Red Nuclear Counterstain
23001	EverBrite™ Mounting Medium
23002	EverBrite™ Mounting Medium with DAPI
23003	EverBrite™ Hardset Mounting Medium
23004	EverBrite™ Hardset Mounting Medium with DAPI
23015	EverBrite™ Mounting Medium with NucSpot® 640
23016	EverBrite™ Hardset Mounting Medium with NucSpot® 640
23017	EverBrite TrueBlack® Hardset Mounting Medium
23018	EverBrite TrueBlack® Hardset Mounting Medium with DAPI
23019	EverBrite TrueBlack® Hardset Mounting Medium with NucSpot® 640
23020	EverBrite™ Mounting Medium with NucSpot® 680
23008	Drop-n-Stain EverBrite™ Mounting Medium
23009	Drop-n-Stain EverBrite™ Mounting Medium with DAPI
23005	CoverGrip™ Coverslip Sealant
22005	Mini Super ^{HT} Pap Pen 2.5 mm tip, ~400 uses
22006	Super ^{HT} Pap Pen 4 mm tip, ~800 uses
23006	Flow Cytometry Fixation/Permeabilization Kit
22023	Paraformaldehyde, 4% in PBS, Ready-to-Use Fixative
22016	Permeabilization Buffer
22017	Permeabilization and Blocking Buffer

Other Mix-n-Stain™ Antibody Labeling Kits

Catalog number	Product
92404-92424, 92454	Mix-n-Stain™ Maxi Labeling Kits
92298-92299	Mix-n-Stain™ R-PE Antibody Labeling Kits
92306-92307	Mix-n-Stain™ APC Antibody Labeling Kits
92340-92341, 92346	Mix-n-Stain™ RPE-CF@647T Antibody Labeling Kits
92310-92311	Mix-n-Stain™ APC-CF@750T Antibody Labeling Kits
92300-92302	Mix-n-Stain™ HRP Antibody Labeling Kits
92314-92315	Mix-n-Stain™ Alkaline Phosphatase Antibody Labeling Kits
92312-92313	Mix-n-Stain™ Glucose Oxidase Antibody Labeling Kits

Please [click here](#) for a full selection of products featuring bright and photostable CF® Dyes, including Mix-n-Stain™ Small Ligand Labeling Kits, primary & secondary antibodies, streptavidin, phalloidins, and much more.

Materials from Biotium are sold for research use only, and are not intended for food, drug, household, or cosmetic use.

CF Dye, Mix-n-Stain, and modified Mix-n-Stain labeling are covered by US and international patents.

Cubitainer is a registered trademark of Hedwin Corporation. Lightning Link is a registered trademark of Expedeon. Zenon and Image-iT are registered trademarks of Thermo Fisher Scientific.