




Designed in Singapore

ESCO
VACCIXCELL
ADHERENT BIOPROCESSING SPECIALIST



BioNOC™ II

The Heart of Tide Motion Bioreactors



fb fisher biotec australia

MORE THAN **20** YEARS IN BUSINESS Established 1997

Australian distributors:
Fisher Biotec Australia
free call: 1800 066 077
email: info@fisherbiotec.com
web: www.fisherbiotec.com

ESCO
Healthcare



"Discovery to Delivery"

Anchorage-dependent cells require surface to adhere in order to grow. Culture vessels such as T-flasks and roller bottles are typically used but the difficulty in scale up have led to the development of carriers. Carriers provide optimum environment for anchorage-dependent cells. The advantages of carriers are ease of scale up, cost reduction from serum and culture media, minimum risk of contamination, and reduction of handling process. Carriers are used to grow virus-generating or protein-producing adherent cells in a large-scale commercial production of vaccines and biologics.

Carriers come in different shapes and sizes. Microcarrier beads are spherical in shape and suitable in stirred tank bioreactors. Another form can be in sheets or fibers which are used in packed-bed or tide motion bioreactors. A wide variety of carriers are available in the market with different physical and chemical properties. The choice of carriers depend greatly on the type of cells to be cultured as different cells have varying anchorage requirements.

BioNOC™ II is a macroporous carrier that supports the growth of anchorage-dependent cells including animal, mammalian, and insect cells in either serum-containing or serum-free culture media. It is made of 100% pure polyethylene terephthalate (PET) nonwoven fabric manufactured according to cGMP guidelines. The special geometric design and surface treatment on the carrier enhance fluid mixing, immobilization efficiency, protection from shear forces, and nutrient transfer during cell culture.



Key Features

- Productivity: up to 5×10^8 cells/g of BioNOC™ II
- Growth area: up to 2,400 cm²/g of BioNOC™ II
- Low lint, convenient for adherent cell growth for protein or viral vaccine production
- Process control and quality assurance (cGMP guideline), full support of documentation


Note: Productivity and growth area varies on the cell line to be used.




BioNOC™ II features improved biocompatibility, long hydrophilicity, high porosity, low lint content and excellent mechanical strength. One gram of BioNOC™ II has a surface area of 2,400 cm² and provides space for up to 1x10⁹ cells to grow. Besides Tide Motion bioreactors, BioNOC™ II is compatible with other commercial packed-bed bioreactor systems. It has been used for a variety of cell strains in either research or production application.



BioNOC™ II macrocarriers



BioNOC™ II carriers are straight and uniform fibers that favor cell distribution and growth.





BioNOC™ II macrocarriers exhibit consistent wettability for over 36 months

Months after exposure to air

BioNOC™ II Structure

Rigid enough to support fixed bed:

- High load-bearing capability to support packed-bed structure
- Increased void fraction and reduced diffusion limitation
- 45° angle generates microscopic eddy for mixing.



Legend:


- BioNOC™ II
- Other Brand

Cell Morphology in BioNOC™ II

BioNOC™ II features improved biocompatibility, long hydrophilicity, high porosity, low lint content and excellent mechanical strength. One gram of BioNOC™ II has a surface area of 2,400 cm² and provides space for up to 1x10⁹ cells to grow. Besides Tide Motion bioreactors, BioNOC™ II is compatible with other commercial packed-bed bioreactor systems. It has been used for a variety of cell strains in either research or production application.

Interception – initial contacts of cells with fiber


Adhesion – cell attachment to fibers and cell-to-cell adhesion to form aggregates







High Porosity and Thickness

- High surface area to volume ratio: 160/cm
- Pore size mainly distributed between 50-200 µm
- Can yield up to 2-3 x 10⁸ VERO cells in 5.5 grams of carriers
- Fiber thickness ~0.4mm
- Allow easy exposure to culture media and aeration
- No oxygen limitation






Stem Cells Grown on and within the BioNOC™ II Fiber Network



Application

Growth of anchorage-dependent cells including animal, mammalian, and insect cells in packed-bed cell culture systems. Compatible with CelCradle™, TideXcell™, and other commercially available packed-bed culture equipment.



Anchorage-Dependent Cells	Insect Cells
Vero, BHK-21, rBHK, CHO-K1, rCHO, rC-127, HEK-293, HEK293A, HepG2, C3A, Hela, Huh-7, RK-13, L929, Human Foreskin Fibroblast Cells, Human Muscle Skeleton Cells, Human Mesenchymal Stem Cells, Human Embryonic Stem Cells	Sf9 Sf21 Hi-5



SEM figure Sf-9 cells in BioNOC™ II macrocarriers


Mesenchymal Stem Cells or MSCs cultured attaches to the 3D matrix and starts to grow in the fibers and the spaces between the fibers. This enables 3D growth and a gene switches on to secrete more extracellular matrix.

MSCs Morphology Alters When Cultured in 2D T-flask vs 3D BioNOC™ II Macrocarriers



MSCs exhibit polyhedral morphology in petri dishes, MSCs in BioNOC™ II are more spindle-like, and is unidirectional

CL-MSCs Secrete Extracellular Matrix (ECM) which Assembles into Fibril Networks on 3D BioNOC™ II




CL-MSCs Secrete Extracellular Matrix (ECM) which Assembles into Fibril Networks on 3D BioNOC™ II



Low Fibre Release

Growth of anchorage-dependent cells including animal, mammalian, and insect cells in packed-bed cell culture systems. Compatible with CelCradle™, TideXcell™, and other commercially available packed-bed culture equipment.




BioNOC™ II has low fibre release which means lower downstream processing costs


- This is favorable for vaccine production and cell therapy

Cells Grown on BioNoc™ II Macrocarriers


▶ Trypan Blue Stain of Vero Cells




▶ SEM Image of Vero Cells



▶ Trypan Blue Stain HEK Cells




▶ HE Stain of MSCs



▶ HE Stain CHO Cells



▶ Hematoxylin Stain of Insect Cells



Specification

Function	For Culture of Adherent Cells
Material	100% non-woven PET filament
Dimension	5 mm × 10 mm strip
Surface Area / g	Up to 2,400 cm ² /g
Capacity	up to 5X10 ⁸ cell/g of BioNOC™ II
Default Packed Volume	15 ml/g (5.5 g for CelCradle™, 110 g for TideXcell™ 2 L matrix vessel, 550 g for TideXcell™ 10 L matrix vessel, 1,100 g for TideXcell™ 20 L matrix vessel, 5,500 g for TideXcell™ 100 L matrix vessel)
Pore Size	50 ~200 μm
Porosity	90~94%
Endotoxin Test	< 0.25 EU/ml
Bioburden Test	< 1 CFU/g
Cytotoxicity Tested	Yes, pyrogen-free
Cell Line	CHO, CHO-K1, rCHO-hlgO, rC-127-TPA, HEK-293, VERO, SF-9, Hi-5, BHK-21, rBHK-Factor VIII, HepG2, HeLa, Huh 7, RK-13, ST, MDCK, MDBK, 3T3, MRC-5, CEF, Human foreskin fibroblast, human muscle skeleton cell, human mesenchymal cell, human embryonic stem cell, etc
Regulation Guideline	USP Class VI, USP <87>,<83>, ISO 10993-5
Storage	Room Temperature (20–22 °C), Keep away from direct sunlight
Shelf Life	2 years
Sterilization	Tolerate autoclave (121° C), gamma irradiation (25 kGy), and EO sterilization

Ordering Information

Product Name	Item Code	Package
BioNOC™ II Cell Culture Carriers (50 g)	1400018	50 gram per bottle
BioNOC™ II Cell Culture Carriers (250 g)	1400019	250 gram per bottle
BioNOC™ II Cell Culture Carriers (1000 g)	1400020	1000 gram per bottle

Literature Support

- Hu, Y.C., J.T. Lu and Y.C. Chung (2003). High-density cultivation of insect cells and production of recombinant baculovirus using a novel oscillating bioreactor. *Cytotechnology*, Volume 42, 145-153.
- Ho, L., C.L. Greene, A.W. Schmidt and L.H. Huang (2004). Cultivation of HEK 293 cell line and production of a member of the superfamily of G-protein coupled receptors for drug discovery applications using a highly efficient novel bioreactor. *Cytotechnology*, Volume 45, 117-123.
- Lu, J.T., Y.C. Chung, Z.R. Chan and Y.C. Hu (2005). A Novel Oscillating Bioreactor BelloCell: Implications for Insect Cell Culture and Recombinant Protein Production. *Biotechnol. Lett.*, Volume 27, 1059-1065. DOI 10.1007/s10529-005-8450-3.
- Mosbeux C., A.R. dos Santos Pedregal, D. Ribeiro de Sousa, V. Hendric, N. Joseph, M. Bensellam, D. Blankaert, T. Marique, C. Alloin, D. Parent, C. Liesnard, J.P. Van Vooren, S. Lowagie and J. Werenne (2006). Mesenchymal cells: metalloproteinases and adhesion on microcarriers. *Animal Cell Technology: Basic & Applied Aspects*, Volume 14, 1-7.
- Drugmand, J.-C., J.-F. Michiels, S.N. Agathos, and Y.-J. Schneider (2007). Growth of Mammalian and Lepidopteran Cells on BioNOC™ II Disks, a Novel Macroporous Microcarrier. *Cell Technology for Cell Products*, Volume 3, 781-784.
- Drugmand, J.-C., N. Havelange, S. Osumba, F. Bosco, F. Debras, F. Collignon, and J. Castillo (2010). New Disposable Fixed-Bed Bioreactor for Cell Culture and Virus Production Based on a Proprietary Agitation and Aeration System. *Cells and Culture*, Volume 4, 151-155.
- Weber, C., S. Pohl, R. Pertner, Pablo Pino-Grace, D. Freimark, C. Wallrapp, P. Geigle and P. Czermark (2010). Production Process for Stem Cell Based Therapeutic Implants: Expansion of the Production Cell Line and Cultivation of Encapsulated Cells. *Adv. Biochem. Eng. Biotechnol.*, Volume 123, 143-162.
- Drugmand J.-C., N. Havelange, F. Collignon, J. Castillo, P.-A. Girod (2012). 4g/l.day: Monoclonal Antibody Volumetric Productivity in the iCELLis™ Disposable Fixed-Bed Bioreactor Proceedings of the 21st Annual Meeting of the European Society for Animal Cell Technology (ESACT), Dublin, Ireland, June 7-10, 2009, Volume 5, 375-378.
- N. Havelange, M. Marigliano, M. Sainte-Marie, F. Debras, N. Tazir and J. Castillo (2012). Poxvirus Production on Chicken Embryo Fibroblasts in iCELLis™ Disposable Fixed-Bed Bioreactor. Proceedings of the 21st Annual Meeting of the European Society for Animal Cell Technology (ESACT), Dublin, Ireland, June 7-10, 2009. Volume 5, 719-722.
- Drugmand, J.-C., F. Collignon, S. Dubois, N. Havelange and J. Castillo (2012). Biomass Sensors in iCELLis™ Fixed-Bed Reactors: Data on CHO and Vero Cells. Proceedings of the 21st Annual Meeting of the European Society for Animal Cell Technology (ESACT), Dublin, Ireland, June 7-10, 2009. Volume 5, 391-394.
- Drugmand, J.-C., G. Esteban, N. Alaoui, N. Jafâr, N. Havelange, O. Berteau and J. Castillo (2012). On-line Monitoring: Animal Cell Cultivation in iCELLis™ Fixed-Bed Reactor Using Dielectric Measurements. Proceedings of the 21st Annual Meeting of the European Society for Animal Cell Technology (ESACT), Dublin, Ireland, June 7-10, 2009. Volume 5, 395-399.
- Sun, B., X.H. Yu, W. Kong, S. Sun, P. Yang, C. Zhu, H. Zhang, Y. Wu, Y. Chen, Y. Shi, X. Zhang and C. Jiang (2013). Production of influenza H1N1 vaccine from MDCK cells using a novel disposable packed-bed bioreactor. *Applied Microbiology and Biotechnology*, Volume 97 (3), 1063-1070.
- Portner, R. (2015). Bioreactors for Mammalian Cells. *Animal Cell Culture*, Volume 9, 89-135.

ESCO[®]
Healthcare



"Discovery to Delivery"

ESCO
VACCIXCELL
ADHERENT BIOPROCESSING SPECIALIST



Australian distributors:
Fisher Biotec Australia
free call: 1800 066 077
email: info@fisherbiotec.com
web: www.fisherbiotec.com