


//76

Axy RACKS



1.5ml and 2.0ml Microtube/Screw-cap tube Racks

80-well Axyracks




- Holds eighty 1.5/2.0ml micro and screw cap tubes
- Polypropylene; resistant to most lab chemicals and solvents
- Perfect for cryogenic storage of tubes
- Conical bottom wells allow superior stability of conical tubes
- Each well numbered for ease of sample identification
- Frosted side writing panel
- Available in a variety of fluorescent and non-fluorescent colors

| Color | Product | Axygen Cat # |
|----------------------|--------------|--------------|
| Assorted Fluorescent | 80-well Rack | R-80-AF |
| Fluorescent Blue | 80-well Rack | R-80-BF |
| Fluorescent Green | 80-well Rack | R-80-GF |
| Fluorescent Orange | 80-well Rack | R-80-OF |
| Fluorescent Red | 80-well Rack | R-80-RF |
| Fluorescent Pink | 80-well Rack | R-80-PF |
| Fluorescent Yellow | 80-well Rack | R-80-YF |

Packaging Description: 25 racks/case

50-well Axyracks




- Holds fifty 2.0/1.5ml micro and screw cap tubes
- polycarbonate; resistant to most lab chemicals and solvents
- Perfect for cryogenic storage to -80°C and below
- Alpha-numeric grid for easy sample identification

| | | |
|----------|--------------|---------|
| Black | 50-well Rack | R-50-BK |
| Blue | 50-well Rack | R-50-B |
| Green | 50-well Rack | R-50-G |
| Red | 50-well Rack | R-50-R |
| Yellow | 50-well Rack | R-50-Y |
| Assorted | 50-well Rack | R-50-A |

Packaging Description: 10 racks/pack, 10 packs/case

81-well Axyracks



- Holds eighty-one 2.0/1.5ml micro and screw cap tubes
- Polycarbonate; resistant to most chemical and organic solvents
- Perfect for cryogenic storage to -80°C and below

| | | |
|--------|--------------|-----------|
| Green | 81-well Rack | R-81-G-PC |
| Red | 81-well Rack | R-81-R-PC |
| Yellow | 81-well Rack | R-81-Y-PC |

Packaging Description: 10 racks/case

48-place aluminum Axyracks



- Holds 48 screw cap tubes of all sizes
- Long lasting aluminum construction
- Unique ribbed-port design allows one-hand operation
- Non-skid benchtop pads for a firm platform
- Center support for added strength and stability


Metallic Anodized Aluminum Rack RA-48

Packaging Description: 4 racks/pack, 4 packs/case

This unique storage rack is made of anodized aluminum and is intended for lab bench work. It features notched holes which 'lock' with Axygen screw cap tubes for one-handed operation.

Screwcap Tube AxyRack
For use with Axygen screw cap microtubes





Axyracks, Caps and CLOSURES

77 //



Color

Assorted

Packaging Description:

Product

80 Well 4 Sided Tube Rack

5 racks/pack, 4 packs/case

Axygen Cat #

R-480-A

- Holds up to 4x 50ml, 16x 15ml, 13x100 or 12x75, 32x 2.0ml or 1.5ml, and 32x 0.5ml Tubes.
- Assorted pack of Pink, Orange, Yellow, Blue, and Green
- Racks that can be positioned on 4-different sides and interlocked to help manage tubes in your lab.

4-sided Tube Rack

