


Deep Well PLATES



384-well Diamond Plates

- Excellent chemical resistance
- Standard SBS footprint
- Centrifugation to 4,000 RCF
- Diamond shaped wells allow complete sample recovery
- Ultra-flat surface ensures proper sealing of heat film

50µl RigiPlate™




Volume to Brim of Well	Working Volume	Volume with AxyMAT	Color	Product Description	Packaging	Axygen Cat #
50µl	49.5µl	-	Clear	384-well RigiPlate™	Bulk	PCR-384-RGD-C
-	-	30µl	Natural	Sealing Mat for 384-well PCR Microplates	Bulk	AM-384-PCR-RD
-	-	6.5µl	Natural	Volume Reducing Sealing Mat	Bulk	AM-384-PCR-RV

Packaging Description: Plates = 10 plates/pack, 2 packs/sleeve, 5 sleeves/case; Mats = 10 mats/pack, 5 packs/case

Colors available by special order. Contact customer service for details.

120µl deep well plate




130µl	120µl	-	Clear	384-well Deep well MicroPlate	Bulk	P-384-120SQ-C
-	-	90µl	Natural	Sealing Mat for 384-well MicroPlates with Square Wells	Bulk	AM-384-DW-SQ
-	-	90µl	Natural	Chemically Resistant Sealing Mat for 384-well MicroPlates with Square Wells	Bulk	AM-384-SQ-IMP

Packaging Description: Plates = 5 plates/pack, 10 packs/case; Mats = 10 mats/pack, 5 packs/case

Black and white plates are available by special order.

240µl deep well plate



250µl	240µl	-	Clear	384-well Deep Well Plate	Bulk	P-384-240SQ-C
-	-	210µl	Natural	Sealing Mat for 384-well Deep Well Plates with Square Wells	Bulk	AM-384-DW-SQ
-	-	210µl	Natural	Chemically Resistant Sealing Mat for 384 Deep Well Plates with Square Wells	Bulk	AM-384-SQ-IMP

Packaging Description: Plates = 5 plates/pack, 10 packs/case; Mats = 10 mats/pack, 5 packs/case

Black and white plates are available by special order.

Maximum RCF for plates - 4,000

EasyPress - an economical applicator for efficiently applying mats to plates (see page 67).





Deep Well PLATES

96-well Round and Rectangular Bottoms



0.6ml deep well plate

Volume to Brim of Well	Working Volume	Volume with Axymat	Color	Product Description	Packaging	Axygen Cat #
620µl	560µl	-	Clear	96-well Deep Well Plate, 600µl	Bulk	P-DW-500-C
620µl	560µl	-	Clear	96-well Deep Well Plate, 600µl Pre-Sterilized	Bulk	P-DW-500-C-S
-	-	440µl	Natural	Sealing Mat for 600µl 96-well Deep Well Plates	Bulk	AM-500UL-RD



Packaging Description: Plates = 5 plates/pack, 10 packs/case; Mats = 10 mats/pack, 5 packs/case



1.1ml deep well plate

Volume to Brim of Well	Working Volume	Volume with Axymat	Color	Product Description	Packaging	Axygen Cat #
1.2ml	1.1ml	-	Clear	96-well Deep Well Plate, 1.1ml	Bulk	P-DW-11-C
1.2ml	1.1ml	-	Clear	96-well Deep Well Plate, 1.1ml Pre-Sterilized	Bulk	P-DW-11-C-S
-	-	1.0ml	Natural	Sealing Mat for 1.1ml 96-well Deep Well Plates	Bulk	AM-2ML-RD
-	-	1.0ml	Natural	Chemically Resistant Sealing Mat for 1.1ml 96-well Deep Well Plates	Bulk	AM-2ML-RD-IMP


Packaging Description: Plates = 5 plates/pack, 10 packs/case; Mats = 10 mats/pack, 5 packs/case



1.6ml deep well plate - square wells

Volume to Brim of Well	Working Volume	Volume with Axymat	Color	Product Description	Packaging	Axygen Cat #
1.7ml	1.6ml	-	Clear	96-well Deep Well Plate, 2.0ml Round Wells	Bulk	P-DW-20-C
1.7ml	1.6ml	-	Clear	96-well Deep Well Plate, 2.0ml Round Wells, Pre-Sterilized	Bulk	P-DW-20-C-S
-	-	1.5ml	Natural	Sealing Mat for 2.0ml 96-well Round Deep Well Plates	Bulk	AM-2ML-RD
-	-	1.5ml	Natural	Chemically Resistant Sealing Mat for 2.0ml 96-well Deep Well Plates with Round Wells	Bulk	AM-2ML-RD-IMP



Packaging Description: Plates = 5 plates/pack, 10 packs/case; Mats = 10 mats/pack, 5 packs/cas



2.2ml deep well plate - round wells

Volume to Brim of Well	Working Volume	Volume with Axymat	Color	Product Description	Packaging	Axygen Cat #
2.3ml	2.2ml	-	Clear	96-well Deep Well Plate, 2.0ml Rectangular Wells	Bulk	P-2ML-SQ-C
2.3ml	2.2ml	-	Clear	96-well Deep Well Plate, 2.0ml Rectangular Wells, Pre-Sterilized	Bulk	P-2ML-SQ-C-S
-	-	1.9ml	Natural	Sealing Mat for 2.0ml 96-well Deep Well Plates with Rectangular Wells	Bulk	AM-2ML-SQ
-	-	1.9ml	Natural	Chemically Resistant Sealing Mat for 2.0ml 96-well Deep Well Plates with Rectangular Wells	Bulk	AM-2ML-SQ-IMP

Packaging Description: Plates = 5 plates/pack, 5 packs/case; Mats = 10 mats/pack, 5 packs/case





// 46 Deep Well/Assay PLATES

Round and 'V' Well options

5.0ml deep well plate - rectangular wells



Volume to Brim of Well	Working Volume	Volume with AxyMat	Color	Product Description	Packaging	Axygen Cat #
4.9ml	4.6ml	-	Clear	48-well Deep Well Plate, 5.0ml Rectangular Wells	Bulk	P-5ML-48-C
4.9ml	4.6ml	-	Clear	48-well Deep Well Plate, 5.0ml Rectangular Wells, Pre-Sterilized	Bulk	P-5ML-48-C-S
-	-	4.0ml	Natural	Chemically Resistant Sealing Mat for 5.0ml 48 Well Deep Well Plates	Bulk	AM-48-IMP

Packaging Description: Plates = 25 plates/case; Mats = 10 mats/pack, 5 packs/case


Round well assay plate



550µl	500µl	-	Clear	96-well Assay Round Bottom Plate	Bulk	P-96-450R-C
550µl	500µl	-	Clear	96-well Assay Round Bottom Plate Pre-Sterilized	Bulk	P-96-450R-C-S
-	-	350µl	Natural	Sealing Mat for 96-well Assay Plates with Round Wells	Bulk	AM-2ML-RD
-	-	350µl	Natural	Chemically Resistant Sealing Mat for 96-well Assay Plates with Round Wells	Bulk	AM-2ML-RD-IMP

Packaging Description: Plates = 10 plates/pack, 5 packs/case; Mats = 10 mats/pack, 5 packs/case


'V' well assay plate



520µl	480µl	-	Clear	96-well Assay 'V' Bottom Plate	Bulk	P-96-450V-C
520µl	480µl	-	Clear	96-well Assay 'V' Bottom Plate Pre-Sterilized	Bulk	P-96-450V-C-S
-	-	300µl	Natural	Sealing Mat for 96-well Assay Plates with 'V' Wells	Bulk	AM-2ML-RD
-	-	300µl	Natural	Chemically Resistant Sealing Mat for 96-well Assay Plates with Round Wells	Bulk	AM-2ML-RD-IMP

Packaging Description: Plates = 10 plates/pack, 5 packs/case; Mats = 10 mats/pack, 5 packs/case

Centrifuge adapter for Deep well plates



-	-	-	Aluminum	Centrifuge adapter for Deep well plates		P-ADAPTER-1
---	---	---	----------	---	--	-------------

Packaging Description: 2 Adapters/case; (includes one of each of Axygen's popular deep well plates)