



plate centrifuge

Patent Pending

- Quickly spin down droplets and condensation
- Use before and after thermal cycling to improve PCR yield
- Accepts skirted, non-skirted and all standard PCR plates
- Less than 1/4 the size of most plate centrifuges



ORDERING INFORMATION

Catalog No.	Description
C1000-AXY*	Axygen plate centrifuge with two position rotor

* Add (-2) for 230V EU or (-3) for 230V UK

SPECIFICATIONS

Capacity:	2 PCR Plates
Rotor Speed:	2500 rpm
G-Force:	500 xg
Rotor:	Fixed Vertically
Dimensions (WxDxH):	7.5 x 8.2 x 7.2 in. 19 x 21 x 18.3 cm
Weight:	6 lbs / 3 kg
Electrical:	120V or 230V, 50-60 Hz

