

# AxySpin™ Chromatography Columns

MINI, MIDI AND MAXI COLUMNS

FOR SLURRY-TYPE CHROMATOGRAPHIC RESINS



- 1.5 ml capacity
- Fits both 1.5 ml and 2 ml microfuge tubes
- Will fit most microcentrifuges



- 2.5 ml capacity
- Fits both 1.5 ml and 2 ml microfuge tubes
- Fits 15 ml conical centrifuge tubes with optional adapter
- Recommended for use in low-speed centrifuges with swinging bucket rotors



- 20 ml capacity
- 3 ml resin bed volume
- Fits 50 ml conical centrifuge tubes
- Recommended for use in low-speed centrifuges with swinging bucket rotors
- **COMING SOON**

## FEATURES:

- Three column sizes for a wide range of purification scales
- Snap-off ends for packaging of prehydrated resin slurries
- Columns fit common, disposable centrifuge tubes
- End caps allow columns to be resealed

AxySpin Columns are designed to transport prehydrated chromatographic media, such as gel filtration and affinity resins. After filling, the columns are prepared for use by simply snapping off the bottom tab. Each is designed to fit commonly available centrifuge tubes. End caps are available for resealing the columns. The three column sizes will accommodate a wide range of different resin volumes. AxySpin Columns are sold either fritted or unfritted by the case. A special 15 ml tube adapter allows the Midi column to be used in disposable 15 ml conical low-speed centrifuge tubes. End caps and polyethylene frits are available separately. AxySpin columns are available either fritted\* or unfritted.

\* FRIT MATERIAL CONSISTS OF 1.5MM THICK HYDROPHOBIC POLYETHYLENE WITH A 20µM PORE SIZE.



# AxySpin™ Chromatography Columns

MINI, MIDI AND MAXI COLUMNS FOR SLURRY-TYPE CHROMATOGRAPHIC RESINS



*AxySpin Mini*



*15 ml tube adapters*



*AxySpin Midi*



*End caps*

## AXYPREP CHROMATOGRAPHY COLUMNS

AxySpin Mini Columns, fritted  
 AxySpin Mini Columns, unfritted  
 AxySpin Midi Columns, fritted  
 AxySpin Midi Columns, unfritted  
 AxySpin Maxi Columns, fritted  
 AxySpin Maxi Columns, unfritted  
 AxySpin Midi 15 ml Tube Adapters  
 AxySpin end caps

## PACK SIZES

4,000/case  
 4,000/case  
 2,000/case  
 2,000/case  
 500/case  
 500/case  
 1,000/case  
 1,000/bag

## CATALOG NO.

AS-MNF  
 Please inquire  
 AS-MDF  
 Please inquire  
 AS-MXF  
 Please inquire  
 AS-MD-ADP  
 AS-EC

