

AxyPrep™ SPE Syringe Column

The AxySep Syringe Column is a fritted, snap-together device which is intended to be filled with a sorbent. It is an exceptionally useful and versatile device for carrying out syringe-mediated SPE (solid phase extraction) chromatography. This column is molded in medical-grade polypropylene and comes disassembled into two components. The internal volume is approximately 1 ml. After filling the lower part of the column, the top component snaps into place. The top component has a complementary fitting which is designed to fit common disposable syringe tips. SPE columns are sold fritted*.



FEATURES:

- *Versatile design*
- *Easy-to-use*
- *Fits all commonly available disposable syringes*

* FRIT MATERIAL CONSISTS OF 1.5MM THICK HYDROPHOBIC POLYETHYLENE WITH A 20µM PORE SIZE.

AXYPREP SPE SYRINGE COLUMNS

SPE Syringe Columns, fritted

PACK SIZES CATALOG NO.

500/bag SPE-SCF