



SERVICE FEATURES AND BENEFITS

Fisher Biotec offers a comprehensive range of quality products at competitive prices and specializes in the manufacture and distribution of leading-edge products for Molecular Biology, Genomics, and Proteomics research.

The company is fully certified to conform with the requirements of AS/NZS ISO 9001:2000 for the manufacture of molecular biology reagents and kits for research and diagnostic applications, and for the distribution of molecular biology products and equipment.

INFORMATION

For more information on any of our products or services please visit us on the web at:

www.fisherbiotec.com

CONTACT US

Freecall 1800 066 077 or

email: info@fisherbiotec.com

Limited time – special offer

Special offer valid until 29 February 2012

Gel Tank Specials – Reduced Prices

Until 29 February 2012 order any Biocom Gel System and make huge savings! Check out our special prices listed below. Send your orders via email to info@fisherbiotec.com or free call us on 1800 066 077.

The Biocom range of gel tanks includes:

The multiSUB series of [Horizontal Gel Units](#) offering the most versatile solution for DNA and RNA agarose gel electrophoresis currently available.

All 5 units offer a combination of economy of gel and buffer volume, with gel size and sample number versatility. Gel size and sample number requirements can be exactly matched in each unit, with the option of additional gel tray sizes. All units feature removable UV transparent trays. Systems are available with 1, 2 or 3 tray options (dependent on model). Easy to use, leak proof 'plug and go' gel casting dams are included as standard to allow gels to be rapidly cast whilst the multiSUB unit is in use for gel running. Lid connectors are compatible with most major power supplies, with adapters available to provide complete compatibility.

The omniPAGE range of [Vertical Gel Units](#) combining ease of use with high resolution separations.

Three sizes, Mini 10 x 10cm, Mini Wide 20 x 10cm and Maxi 20 x 20cm, share a host of common features including a guaranteed leak proof seal required for trouble free, rapid and uncomplicated gel casting. Utilising a built in gel running module eliminates time consuming transfer of glass plates during casting. Glass plates with permanently bonded spacers guarantee perfect spacer alignment. Fast set up times while even pressure bars and ultra soft seals guarantee lead proof casting. Once the gel has polymerised, the gel running module is simply inserted into the gel tank for electrophoresis.

See the Biocom [Catalog](#) for a complete range of gel systems.



Listed below are examples of savings that can be made. Special prices apply to the complete range of Biocom Gel Systems – contact us for more information.

| Catalog No.* | Model | Description | List Price | Special Price |
|-------------------|-----------------------|--|------------|---------------|
| BCMSMINIDUO-WS | Multi Sub Mini Duo | 7 x 7 & 7 x 10cm UV Tray, 2 x 8-sample 1mm thick combs, casting dams, and loading guides | \$367 | \$295 |
| BCMSMIDIDUO-WS | Multi Sub Midi Duo | 10 x 7 & 10 x 10cm UV Tray, 2 x 16-sample 1mm thick combs, casting dams, and loading guides | \$397 | \$315 |
| BCMSCHOICETRIO-WS | Multi Sub Choice Trio | 15 x 7, 15 x 10 & 15 x 15cm UV Tray, 2 x 20-sample 1mm thick combs, casting dams, and loading guides | \$453 | \$360 |
| BCMSMAXIDUO-WS | Multi Sub Maxi Duo | 20 x 10 & 20 x 20cm UV Tray, 2 x 20-sample 1mm thick combs, casting dams, and loading guides | \$777 | \$620 |
| BCMSSCREENTRIO-WS | Multi Sub Screen Trio | 26 x 16, 26 x 24 & 26 x 32cm UV Tray, 6 x 28-sample 1mm thick combs, casting dams, and loading guides | \$1,013 | \$810 |
| BCVS10DSYS-WS | omniPAGE Mini | (Sliding Gates & Screw Clamps) Vertical, 10 x 10cm Dual, 2 sets of Glass Plates, 1mm thick bonded Spacers, 2 x 12-sample 1mm thick combs, and caster | \$712 | \$570 |
| BCVS10CBS-V-WS | omniPAGE Mini Blot | (Sliding Gates and Screw Clamps) complete system for Mini Vertical Electrophoresis and Blotting including: Vertical unit, blotting insert, and accessories | \$992 | \$795 |
| BCVS10WDSYS-WS | omniPAGE Mini Wide | 20 x 10cm Dual, 2 sets of Glass Plates with 1mm bonded Spacers, and 2 x 24-sample 1mm thick combs, and caster | \$932 | \$745 |
| BCVS20DSYS-WS | omniPAGE Maxi | 20 x 20cm Dual, 2 sets of Glass Plates, 1mm thick bonded Spacers, 24-sample 1mm thick combs, and Caster | \$1,346 | \$1,075 |

* Quote these Catalog Nos. to receive offer.

Visit us on the web at: www.fisherbiotec.com