

# Buffers

## 10x TBE (Tris-Borate-EDTA) Buffer

Each pack contains 4 sachets of buffer mix. Each Sachet is sufficient to make 1 litre of 10x buffer.

Catalog No.	Pack Size
TBE-1	4 sachets

## 10x Tris- Glycine (Laemmli) Buffer

Each pack contains 4 sachets of buffer mix. Each Sachet is sufficient to make 1 litre of 10x buffer.

Catalog No.	Pack Size
TGB-1	4 sachets

## 10x PBS (Phosphate Buffered Saline)

This high purity starch has been specifically prepared and tested for starch gel electrophoresis applications. Suitable for applications such as serum protein analysis and isozyme analysis.

Catalog No.	Pack Size
PBS-1	4 sachets

## 2x Gel Loading Buffer

Contains 5% Ficoll 400, 0.8% bromophenol blue, and 0.8% xylene cyanol in TE Buffer pH 8.0. The bromophenol blue and xylene cyanol dyes are used for tracking migration during electrophoresis. In 0.5 - 1.2% agarose gels using 0.5x TBE bromophenol blue migrates at approximately 300bp and xylene cyanol at approximately 4kb.

Catalog No.	Pack Size
GLB2-5	5 ml

## 6x Gel Loading Buffer

Contains 15% Ficoll 400, 0.25% bromophenol blue, and 0.25% xylene cyanol in TE Buffer pH 8.0. The bromophenol blue and xylene cyanol dyes are used for tracking migration during electrophoresis. In 0.5 - 1.2% agarose gels using 0.5x TBE bromophenol blue migrates at approximately 300bp and xylene cyanol at approximately 4kb.

Catalog No.	Pack Size
GLB6-2	1 ml
GLB6-5	5 ml
GLB6-10	10 ml

### SERVICE FEATURES AND BENEFITS

Fisher Biotec offers a comprehensive range of quality products at competitive prices and specialises in the manufacture and distribution of leading-edge products for Molecular Biology, Genomics, and Proteomics research.

Our management system has been certified to JAS-ANZ ISO 9001:2008 for the sale and supply of equipment and consumables to the pharmaceutical and laboratory health care industry.

### INFORMATION

For more information on any of our products or services please visit us on the Web at:

[www.fisherbiotec.com](http://www.fisherbiotec.com)

### CONTACT US

Freecall 1800 066 077 or

email: [info@fisherbiotec.com](mailto:info@fisherbiotec.com)



fisher biotec  
australia



Quality  
ISO 9001

SAI GLOBAL