

Limited time – special offer

Special offer valid until 30 April 2012

Electroporation Cuvette Specials

BTX manufactures universal fit **Cuvettes Plus™** for use in electroporation of bacteria, yeast, mammalian and plant cells. Cuvettes are compatible with most electroporation systems.

Pay **only \$190** if you order any **BTX Electroporation Cuvettes Plus** pack by 30 April.

45-0124 (610) WS*
45-0125 (620) WS*
45-0126 (640) WS*

Cuvette Plus 1mm 90µl, Sterile Pkg/50 Grey
Cuvette Plus 2mm 400µl, Sterile Pkg/50 Blue
Cuvette Plus 4mm 800µl, Sterile Pkg/50 Yellow

*Quote Catalog Nos. to receive offer



SERVICE FEATURES AND BENEFITS

Fisher Biotec offers a comprehensive range of quality products at competitive prices and specializes in the manufacture and distribution of leading-edge products for Molecular Biology, Genomics, and Proteomics research.

Our management system has been certified to JAS-ANZ ISO 9001:2008 for the sale and supply of equipment and consumables to the pharmaceutical and laboratory health care industry.

INFORMATION

For more information on any of our products or services please visit us on the Web at:

www.fisherbiotec.com

CONTACT US

Freecall 1800 066 077 or

email: info@fisherbiotec.com

fb

fisher biotec
australia

BTX[®]
HARVARD APPARATUS



KEY FEATURES

- Sterile Transfer Pipette included with every Cuvette Plus package, allowing quick and easy removal of each sample
- Three cuvette sizes: 1 mm, 2 mm and 4 mm, are available to suite all applications
- Round cuvette caps allow for single-hand removal and are colour coded for easy recognition

PRODUCT DESCRIPTION

BTX® Cuvettes Plus are designed for use in electroporation of bacteria, yeast, mammalian and plant cells. Each sterilized Cuvettes Plus package includes a disposable cuvette and a transfer pipette. The cuvettes are molded with embedded aluminum electrodes, washed prior to packaging and gamma irradiated for guaranteed sterility. Electrode gap sizes are available in 1mm, 2 mm and 4 mm for high field strengths up to 25.0 kV/cm. Colour-coded round caps are leak resistant and allow for quick and easy removal. BTX® disposable Cuvettes Plus require no maintenance. Simply discard carefully after use. Cuvettes Plus are also compatible with most commercially available electroporators.

Send your orders via email to: info@fisherbiotec.com
or FREECALL us on **1800 066 077**



Quality
ISO 9001
SAI GLOBAL