



Gel Cutting Tips

A BETTER WAY TO EXCISE SPOTS & BANDS FROM GELS

- designed for removing bands from agarose gels
- one handed operation
- push button gel release
- ejector button tip release
- avoid cross contamination
- tips fit standard 1000 μ l pipettors



compatible with
the new Axygen
AXYPET

Ordering details:

Cat.No.

TGL-1140

TGL-1165

TGL-1140-R

TGL-1165-R

Description

1.1 mm x 4 mm Gel Cutting Tips

1.1 mm x 6.5 mm Gel Cutting Tips

1.1 mm x 4 mm Gel Cutting Tips

1.1 mm x 6.5 mm Gel Loading Tips

Packaging

Bag/250 tips

Bag/250 tips

Rack/48 tips

Rack/48 tips

Packaging Description: Bulk = 250 tips/pack, 10 packs/case;
Racked = 48 tips/rack, 10 racks/unit, 5 units/case

www.axygen.com

AXYGEN[®]

Innovation and Integrity