

H1 High Grade Membrane

CelluSep H1 individual-use 50 cm lengths are packed wet in hard plastic vials and are available in quantities of 10 or 5 vials per pack. H1 is also available in 10m and 5m continuous rolls.

CelluSep H1

PREPARED DIALYSIS MEMBRANES



High grade, ultra-clean, ready to use

Heavy-metal free, packed wet, EDTA treated

H1 10m Rolls

Part No.	Flat Width	Nominal MWCO	Vol/cm
1510-45	45 mm	15,000	6.42
1510-34	34 mm	15,000	3.46
1510-12	12 mm	15,000	0.46
1510-25	25 mm	15,000	1.98
1010-50	50 mm	10,000	7.94
1010-45	45 mm	10,000	6.42
1010-32	32 mm	10,000	3.27
1010-24	24 mm	10,000	1.67
1010-12	12 mm	10,000	0.46
0810-40	40 mm	8,000	5.10
0810-50	50 mm	8,000	7.94
0810-24	24 mm	8,000	1.67
0810-32	32 mm	8,000	3.27
0810-12	12 mm	8,000	0.46
0310-19	19 mm	3,500	1.15
0310-46	45 mm	3,500	6.74

H1 50cm Vials - 10 Pack

Part No.	Flat Width	Nominal MWCO	Vol/cm
1-5050-34	34 mm	50,000	3.46
1-2550-40	40 mm	25,000	5.10
1-5050-28	28 mm	50,000	2.50
1-2550-35	35 mm	25,000	3.91
1-2550-12	12 mm	25,000	.46
1-2550-25	25 mm	25,000	1.98
1-1550-32	32 mm	15,000	3.27
1-1550-45	45 mm	15,000	6.42
1-1550-25	25 mm	15,000	1.98
1-1050-40	40 mm	10,000	5.10
1-1550-10	10 mm	15,000	.32
1-1050-32	32 mm	10,000	3.27
1-1050-10	10 mm	10,000	.32
1-1050-25	25 mm	10,000	1.98
1-0550-40	40 mm	5,000	5.10
1-0550-30	30 mm	5,000	2.86
1-0550-20	20 mm	5,000	1.27
1-0250-45	45 mm	2,000	6.42

Part No.	Flat Width	Nominal MWCO	Vol/cm
1-0550-18	18 mm	5,000	1.04
1-0250-18	18 mm	2,000	1.04
1-0250-38	38 mm	2,000	4.59
1-0150-45	45 mm	1,000	6.42
1-0150-38	38 mm	1,000	4.59
1-0150-18	18 mm	1,000	1.04

H1 50cm Vials - 5 Pack

Part No.	Flat Width	Nominal MWCO	Vol/cm
5-5050-34	34 mm	50,000	3.46
5-2550-40	40 mm	25,000	5.10
5-5050-28	28 mm	50,000	2.50
5-2550-35	35 mm	25,000	3.91
5-2550-12	12 mm	25,000	0.46
5-2550-25	25 mm	25,000	1.98
5-1550-32	32 mm	15,000	3.27
5-1550-45	45 mm	15,000	6.42
5-1550-25	25 mm	15,000	1.98
5-1550-10	10 mm	15,000	0.32
5-1050-40	40 mm	10,000	5.10
5-1050-32	32 mm	10,000	3.27
5-1050-10	10 mm	10,000	0.32
5-1050-25	25 mm	10,000	1.98
5-0550-40	40 mm	5,000	5.10
5-0550-30	30 mm	5,000	2.86
5-0550-20	20 mm	5,000	1.27
5-0250-45	45 mm	2,000	6.42
5-0550-18	18 mm	5,000	1.04
5-0250-38	38 mm	2,000	4.59
5-0250-18	18 mm	2,000	1.04
5-0150-45	45 mm	1,000	6.42
5-0150-18	18 mm	1,000	1.04
5-0150-38	38 mm	1,000	4.59

H1 5m Rolls

Part No.	Flat Width	Nominal MWCO	Vol/cm
1505-45	45 mm	15,000	6.42
1505-34	34 mm	15,000	3.46
1505-25	25 mm	15,000	1.98
1505-12	12 mm	15,000	0.46
1005-45	45 mm	10,000	6.42
1005-50	50 mm	10,000	7.94
1005-32	32 mm	10,000	3.27
1005-12	12 mm	10,000	0.46
1005-24	24 mm	10,000	1.67
0805-40	40 mm	8,000	5.10
0805-50	50 mm	8,000	7.94
0805-24	24 mm	8,000	1.67
0805-32	32 mm	8,000	3.27
0805-12	12 mm	8,000	0.46