

Protein Pre-Cast Gel Complete Kit

Protein Pre-Cast Gel Complete Kit is a convenient combination package of pre-cast gels, buffers, markers and staining solution. It is sufficient for 10 mini-gel assays.



■ Kit Contents

- Pre-Cast Gel : Choice of EzWay™ Gel (1 year shelf life) or Tris-Glycine Gel (standard type of gel)
- Sample Buffer
- Running Buffer
- Protein Marker for MW estimation
- Prestained Protein Marker for transfer efficiency confirmation
- Staining Solution

■ KG0001 Protein Pre-Cast Gel Complete Kit I

Cat.No.	Description	Size
KG XXXX	Tris-Glycine-PAG, gel (any combination)	10/pk
KTG030	Tris-Glycine SDS Running Buffer (10X)	500 ml
KTG020-5	Tris-Glycine SDS Sample Buffer (5X)	1.5 ml
K18000-S	EzWay™ Protein-PreBlue MW Marker, 10 bands	10 applications
K08000-S	EzWay™ Protein-Blue MW Marker, 11 bands	10 applications
K14020-S	EzWay™ Protein-Blue II Staining Solution	125 ml x 2 ea

■ KG0002 Protein Pre-Cast Gel Complete Kit II

Cat.No.	Description	Size
KG XXXX	EzWay-PAG, gel (any combination)	10/pk
KTG030	Tris-Glycine SDS Running Buffer (10X)	500 ml
KTG020-5	Tris-Glycine SDS Sample Buffer (5X)	1.5 ml
K18000-S	EzWay™ Protein-PreBlue MW Marker, 10 bands	10 applications
K08000-S	EzWay™ Protein-Blue MW Marker, 11 bands	10 applications
K14020-S	EzWay™ Protein-Blue II Staining Solution	125 ml x 2 ea

■ KG0003 Protein Pre-Cast Gel Complete Kit III

Cat.No.	Description	Size
KG XXXX	EzWay-PAG, gel (any combination)	10/pk
KTR030	Tricine Running Buffer (10X)	500 ml
KTR020-5	Tricine Sample Buffer (5X)	1.5 ml
K18000-S	EzWay™ Protein-PreBlue MW Marker, 10 bands	10 applications
K08000-S	EzWay™ Protein-Blue MW Marker, 11 bands	10 applications
K14020-S	EzWay™ Protein-Blue II Staining Solution	125 ml x 2 ea

■ KG0004 Protein Pre-Cast Gel Complete Kit IV

Cat.No.	Description	Size
KG XXXX	EzWay-PAG, gel (any combination)	10/pk
KAS030	Aspartate Running Buffer (10X)	500 ml
KAS020-5	Aspartate Sample Buffer (5X)	1.5 ml
K18000-S	EzWay™ Protein-PreBlue MW Marker, 10 bands	10 applications
K08000-S	EzWay™ Protein-Blue MW Marker, 11 bands	10 applications
K14020-S	EzWay™ Protein-Blue II Staining Solution	125 ml x 2 ea