

Microtubes




for the ultimate in sample integrity

THE DIFFERENCE IS CLEAR

AXYGEN'S MICROCENTRIFUGE TUBES ARE PRODUCED FROM MEDICAL GRADE POLYPROPYLENE

- **MAXIMUM RECOVERY & Standard Snaplock tubes** – 0.5, 1.5, 1.7 and 2.0ml
- **MAXIMUM RECOVERY & Standard Screw-cap tubes** – 0.5, 1.5 and 2.0ml conical and free-standing
- **Durable Plastic Tube Racks** – 50-, 80- and 81-position
- **Aluminum Tube Racks** – 48-position for screw cap tubes

proteomics • drug discovery • genomics • high throughput screening



The closure on Oxygen Screw Cap Tubes feature an 'O' ring that is formulated and manufactured from a specially selected material that is resistant to many commonly-used laboratory chemicals. The Alignment Orientation Tab (US patent # 5891327) allows the microtube to be aligned with the centrifuge rotor for multiple applications, including desalting. This ensures that the microtube can be re-spun without disturbing the partially pelletized sediment.


- **MAXYclear™** for the ultimate in sample visibility - found in all of Axygen's SnapLock tubes
- A wide range of colored tubes for easy sample identification
- **SnapLock tubes are guaranteed to stay closed** - even during boiling or autoclaving!

MAXYclear SnapLock Microtubes

Axygen 1.5ml, 1.7ml and 2.0ml Snaplock microtubes are engineered to close with a "snap" and stay closed during centrifugation up to 20,000 RCF. The 0.6ml and 1.7ml MAXYclear tubes are manufactured of a copolymer material to offer exceptional strength and chemical resistance. All Axygen SnapLock microtubes have frosted cap surfaces for labelling, frosted panels on the side for writing, a 0.28 mm piercing port in the center of the cap, are tested to 55 psi and are ready for autoclaving.

MAXYclear Screw Cap Microtubes

Axygen screw cap microtubes serve a broad range of laboratory needs including collection, processing, centrifugation, long-term storage and packaging/shipping of restriction enzymes and other reagents. Color-coded caps provide a syringe port which is easily pierced for radioisotope labelling and fraction collecting. The threaded 'O' ring design ensures a uniform closure allowing centrifugation up to 20,000 RCF. Ribs are designed to fit racks and facilitate single-handed opening and closing when used with our aluminum AxyRack, found on page 76. See page 79 for information on RCF, or "g-force."



- Will withstand up to 20,000 RCF
- SnapLock cap assures efficient seal, yet is still easy to 'flip' open
- Thin membrane in cap center allows easy access to sample by syringe needle
- Certified RNase-, DNase- and endotoxin-safe for total sample integrity
- **MAXIMUM RECOVERY** - available in most styles for minimal sample retention

Integrity



// 70

Snaplock MICROTUBES



Maxyclear

- Co-polymer tubes offer exceptional strength and chemical resistance
- Homopolymer tubes are more resistant to autoclaving
- 0.6ml, 1.5ml and 2.0ml microtubes are certified boil proof

0.6ml - Maxyclear



Color	Volume	Microtube Material	Packaging	Axygen Cat #
MaxyClear	0.6ml	Homopolymer	Bulk	MCT-060-C
Assorted*	0.6ml	Homopolymer	Bulk	MCT-060-A
Blue	0.6ml	Homopolymer	Bulk	MCT-060-B
Green	0.6ml	Homopolymer	Bulk	MCT-060-G
Orange	0.6ml	Homopolymer	Bulk	MCT-060-O
Red	0.6ml	Homopolymer	Bulk	MCT-060-R
Violet	0.6ml	Homopolymer	Bulk	MCT-060-V
Amber	0.6ml	Homopolymer	Bulk	MCT-060-X
Yellow	0.6ml	Homopolymer	Bulk	MCT-060-Y

Packaging Description:
Maximum G-force:

1000 tubes/pack, 10 packs/case
14,000 RCF

1.7ml - Maxyclear



MaxyClear	1.7ml	Co-polymer	Bulk	MCT-175-C
Assorted*	1.7ml	Co-polymer	Bulk	MCT-175-A
Blue	1.7ml	Co-polymer	Bulk	MCT-175-B
Green	1.7ml	Co-polymer	Bulk	MCT-175-G
Orange	1.7ml	Co-polymer	Bulk	MCT-175-O
Red	1.7ml	Co-polymer	Bulk	MCT-175-R
Violet	1.7ml	Co-polymer	Bulk	MCT-175-V
Amber	1.7ml	Co-polymer	Bulk	MCT-175-X
Yellow	1.7ml	Co-polymer	Bulk	MCT-175-Y

Packaging Description:
Maximum G-force:

500 tubes/pack, 10 packs/case
20,000 RCF

1.5ml - Maxyclear



MaxyClear	1.5ml	Homopolymer	Bulk	MCT-150-C
Assorted*	1.5ml	Homopolymer	Bulk	MCT-150-A
Blue	1.5ml	Homopolymer	Bulk	MCT-150-B
Green	1.5ml	Homopolymer	Bulk	MCT-150-G
Orange	1.5ml	Homopolymer	Bulk	MCT-150-O
Red	1.5ml	Homopolymer	Bulk	MCT-150-R
Violet	1.5ml	Homopolymer	Bulk	MCT-150-V
Amber	1.5ml	Homopolymer	Bulk	MCT-150-X
Yellow	1.5ml	Homopolymer	Bulk	MCT-150-Y

Packaging Description:
Maximum G-force:

500 tubes/pack, 10 packs/case
20,000 RCF

2.0ml - Maxyclear




MaxyClear	2.0ml	Homopolymer	Bulk	MCT-200-C
Assorted*	2.0ml	Homopolymer	Bulk	MCT-200-A
Blue	2.0ml	Homopolymer	Bulk	MCT-200-B
Green	2.0ml	Homopolymer	Bulk	MCT-200-G
Orange	2.0ml	Homopolymer	Bulk	MCT-200-O
Red	2.0ml	Homopolymer	Bulk	MCT-200-R
Violet	2.0ml	Homopolymer	Bulk	MCT-200-V
Amber	2.0ml	Homopolymer	Bulk	MCT-200-X
Yellow	2.0ml	Homopolymer	Bulk	MCT-200-Y

Packaging Description:
Maximum G-force:

500 tubes/pack, 10 packs/case
20,000 RCF

* colors included: blue, green, violet, orange, red and yellow in approximately equal quantities



Snaplock MICROTUBES

71 //

Maxyclear

Maxyclear



Color	Volume	Microtube Material / Description	Packaging	Axygen Cat #
Spectrum	0.6ml	Homopolymer	Bulk	MCT-060-SP
Spectrum	1.5ml	Homopolymer	Bulk	MCT-150-SP
Spectrum	1.7ml	Co-polymer	Bulk	MCT-175-SP
Spectrum	2.0ml	Homopolymer	Bulk	MCT-200-SP

Packaging Description:

0.6ml = 1000 tubes/pack, 10 packs/case; other sizes = 500 tubes/pack, 10 packs/case. 2 packs of blue, green, orange, red and yellow tubes per case

Maximum G-force:

0.6ml = 14,000 RCF; other size tubes = 20,000 RCF



Maxymum Recovery				
Color	Volume	Microtube Material / Description	Packaging	Axygen Cat #
MaxyClear	0.6ml	Homopolymer, Maxymum Recovery	Bulk	MCT-060-L-C
MaxyClear	1.5ml	Homopolymer, Maxymum Recovery	Bulk	MCT-150-L-C
MaxyClear	1.7ml	Co-polymer, Maxymum Recovery	Bulk	MCT-175-L-C
MaxyClear	2.0ml	Homopolymer, Maxymum Recovery	Bulk	MCT-200-L-C

Packaging Description:

0.6ml = 500 tubes/pack, 10 packs/case; other sizes = 250 tubes/pack, 10 packs/case

Maximum G-force:

0.6ml = 14,000 RCF; other size tubes = 20,000 RCF



Pre-sterilized				
Color	Volume	Microtube Material / Description	Packaging	Axygen Cat #
MaxyClear	0.6ml	Homopolymer, Pre-Sterilized	Bulk	MCT-060-C-S
MaxyClear	1.5ml	Homopolymer, Pre-Sterilized	Bulk	MCT-150-C-S
MaxyClear	1.7ml	Co-polymer, Pre-Sterilized	Bulk	MCT-175-C-S
MaxyClear	2.0ml	Homopolymer, Pre-Sterilized	Bulk	MCT-200-C-S

Packaging Description:

0.6ml = 500 tubes/pack, 10 packs/case; other sizes = 250 tubes/pack, 10 packs/case

Maximum G-force:

0.6ml = 14,000 RCF; other size tubes = 20,000 RCF

Capless				
Color	Volume	Microtube Material / Description	Packaging	Axygen Cat #
MaxyClear	1.5ml	Homopolymer, No Cap	Bulk	MCT-150-NC
MaxyClear	2.0ml	Homopolymer, No Cap	Bulk	MCT-200-NC


Packaging Description:

500 tubes/pack, 10 packs/case


Maximum G-force:

20,000 RCF

1.5 and 2.0ml



for 1.5ml and 2.0ml tubes



Tissue Grinder				
Color	Volume	Microtube Material / Description	Packaging	Axygen Cat #
Blue	1.5 & 2.0ml	Pestle, Pre-Sterilized	Individually Wrapped	PES-15-B-SI

Packaging Description:

Individually Wrapped, 100 /pack, 4 packs/case

Micro Tube Packaging

All Microtubes are shipped in re-sealable plastic bags packed in individual boxes – keeping your supply safe from contamination and easily stored.

Autoclaving Protocol

Autoclave at 121°C/250°F at 15 psi (1 atm.) for 15 minutes. Do not use drying cycle. Remove from autoclave once temperature has reached at least 80°C.

Autoclaving Caution

Excessive exposure to autoclaving heat may cause warping or decreased structural integrity. Do not exceed 121°C / 250°F, 15 psi (1 atm.), 15 minutes.



- Durable and inert polypropylene construction
- Designed to fit 1.5ml & 2.0ml centrifuge tubes
- Individually wrapped and pre-sterilized for convenience

www.axygen.com

// 72

Screw Cap TUBES



With Caps

0.5ml, conical



All non-sterile cap tubes are packaged with unattached, non-looped caps with 'o'-rings

Color	Volume	Style	Tube Type	Oxygen Cat #
Clear	0.5ml	Conical	Screw Cap Tubes and Caps	SCT-45-C
Clear	0.5ml	Conical	Screw Cap Tubes and Caps	SCT-050-C
Assorted†	0.5ml	Conical	Screw Cap Tubes and Caps	SCT-050-A
Grey	0.5ml	Conical	Screw Cap Tubes and Caps	SCT-050-GR
White	0.5ml	Conical	Screw Cap Tubes and Caps	SCT-050-W
Pink	0.5ml	Conical	Screw Cap Tubes and Caps	SCT-050-P
Red	0.5ml	Conical	Screw Cap Tubes and Caps	SCT-050-R
Orange	0.5ml	Conical	Screw Cap Tubes and Caps	SCT-050-O
Yellow	0.5ml	Conical	Screw Cap Tubes and Caps	SCT-050-Y
Green	0.5ml	Conical	Screw Cap Tubes and Caps	SCT-050-G
Blue	0.5ml	Conical	Screw Cap Tubes and Caps	SCT-050-B
Violet	0.5ml	Conical	Screw Cap Tubes and Caps	SCT-050-V
Brown	0.5ml	Conical	Screw Cap Tubes and Caps	SCT-050-BR

Packaging Description: 500 tubes/pack, 8 packs/case: 150/bag, 10 bags/pack, 8 packs/case
Maximum G-force: 20,000 RCF

0.5ml, self standing



Color	Volume	Style	Tube Type	Oxygen Cat #
Clear	0.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-050-SS-C
Assorted†	0.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-050-SS-A
Grey	0.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-050-SS-GR
White	0.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-050-SS-W
Pink	0.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-050-SS-P
Red	0.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-050-SS-R
Orange	0.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-050-SS-O
Yellow	0.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-050-SS-Y
Green	0.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-050-SS-G
Blue	0.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-050-SS-B
Violet	0.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-050-SS-V
Brown	0.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-050-SS-BR

Packaging Description: 500 tubes/pack, 8 packs/case: 150/bag, 10 bags/pack, 8 packs/case
Maximum G-force: 20,000 RCF

1.5ml, conical



Color	Volume	Style	Tube Type	Oxygen Cat #
Clear	1.5ml	Conical	Screw Cap Tubes and Caps	SCT-150-C
Assorted†	1.5ml	Conical	Screw Cap Tubes and Caps	SCT-150-A
Grey	1.5ml	Conical	Screw Cap Tubes and Caps	SCT-150-GR
White	1.5ml	Conical	Screw Cap Tubes and Caps	SCT-150-W
Pink	1.5ml	Conical	Screw Cap Tubes and Caps	SCT-150-P
Red	1.5ml	Conical	Screw Cap Tubes and Caps	SCT-150-R
Orange	1.5ml	Conical	Screw Cap Tubes and Caps	SCT-150-O
Yellow	1.5ml	Conical	Screw Cap Tubes and Caps	SCT-150-Y
Green	1.5ml	Conical	Screw Cap Tubes and Caps	SCT-150-G
Blue	1.5ml	Conical	Screw Cap Tubes and Caps	SCT-150-B
Violet	1.5ml	Conical	Screw Cap Tubes and Caps	SCT-150-V
Brown	1.5ml	Conical	Screw Cap Tubes and Caps	SCT-150-BR

Packaging Description: 500 tubes/pack, 8 packs/case: 150/bag, 10 bags/pack, 8 packs/case
Maximum G-force: 20,000 RCF

1.5ml, self standing



Color	Volume	Style	Tube Type	Oxygen Cat #
Clear	1.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-150-SS-C
Assorted†	1.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-150-SS-A
Grey	1.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-150-SS-GR
White	1.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-150-SS-W
Pink	1.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-150-SS-P
Red	1.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-150-SS-R
Orange	1.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-150-SS-O
Yellow	1.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-150-SS-Y
Green	1.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-150-SS-G
Blue	1.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-150-SS-B
Violet	1.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-150-SS-V
Brown	1.5ml	Self Standing	Screw Cap Tubes and Caps	SCT-150-SS-BR

Packaging Description: 500 tubes/pack, 8 packs/case: 150/bag, 10 bags/pack, 8 packs/case
Maximum G-force: 20,000 RCF

* Maximum Recovery Screw Cap Tubes available as special order items. Please call our customer service department for more information.



Screw Cap TUBES

With and Without Caps



2.0ml, conical

Color	Volume	Style	Tube Type	Axygen Cat #
Clear	2.0ml	Conical	Screw Cap Tubes and Caps	SCT-200-C
Assorted [†]	2.0ml	Conical	Screw Cap Tubes and Caps	SCT-200-A
Grey	2.0ml	Conical	Screw Cap Tubes and Caps	SCT-200-GR
White	2.0ml	Conical	Screw Cap Tubes and Caps	SCT-200-W
Pink	2.0ml	Conical	Screw Cap Tubes and Caps	SCT-200-P
Red	2.0ml	Conical	Screw Cap Tubes and Caps	SCT-200-R
Orange	2.0ml	Conical	Screw Cap Tubes and Caps	SCT-200-O
Yellow	2.0ml	Conical	Screw Cap Tubes and Caps	SCT-200-Y
Green	2.0ml	Conical	Screw Cap Tubes and Caps	SCT-200-G
Blue	2.0ml	Conical	Screw Cap Tubes and Caps	SCT-200-B
Violet	2.0ml	Conical	Screw Cap Tubes and Caps	SCT-200-V
Brown	2.0ml	Conical	Screw Cap Tubes and Caps	SCT-200-BR



Packaging Description: 500 tubes/pack, 8 packs/case; [†]50/bag, 10 bags/pack, 8 packs/case
Maximum G-force: 20,000 RCF

2.0ml, self standing

Clear	2.0ml	Self Standing	Screw Cap Tubes and Caps	SCT-200-SS-C
Assorted [†]	2.0ml	Self Standing	Screw Cap Tubes and Caps	SCT-200-SS-A
Grey	2.0ml	Self Standing	Screw Cap Tubes and Caps	SCT-200-SS-GR
White	2.0ml	Self Standing	Screw Cap Tubes and Caps	SCT-200-SS-W
Pink	2.0ml	Self Standing	Screw Cap Tubes and Caps	SCT-200-SS-P
Red	2.0ml	Self Standing	Screw Cap Tubes and Caps	SCT-200-SS-R
Orange	2.0ml	Self Standing	Screw Cap Tubes and Caps	SCT-200-SS-O
Yellow	2.0ml	Self Standing	Screw Cap Tubes and Caps	SCT-200-SS-Y
Green	2.0ml	Self Standing	Screw Cap Tubes and Caps	SCT-200-SS-G
Blue	2.0ml	Self Standing	Screw Cap Tubes and Caps	SCT-200-SS-B
Violet	2.0ml	Self Standing	Screw Cap Tubes and Caps	SCT-200-SS-V
Brown	2.0ml	Self Standing	Screw Cap Tubes and Caps	SCT-200-SS-BR



Packaging Description: 500 tubes/pack, 8 packs/case; [†]50/bag, 10 bags/pack, 8 packs/case
Maximum G-force: 20,000 RCF

Pre-sterilized

Clear	0.5ml	Conical	Pre-Sterilized [‡]	SCT-45-C-S
Clear	0.5ml	Conical	Pre-Sterilized [‡]	SCT-050-C-S
Assorted [†]	0.5ml	Conical	Pre-Sterilized [‡]	SCT-050-A-S
Clear	0.5ml	Self Standing	Pre-Sterilized [‡]	SCT-050-SS-C-S
Assorted [†]	0.5ml	Self Standing	Pre-Sterilized [‡]	SCT-050-SS-A-S
Clear	1.5ml	Conical	Pre-Sterilized [‡]	SCT-150-C-S
Assorted [†]	1.5ml	Conical	Pre-Sterilized [‡]	SCT-150-A-S
Clear	1.5ml	Self Standing	Pre-Sterilized [‡]	SCT-150-SS-C-S
Assorted [†]	1.5ml	Self Standing	Pre-Sterilized [‡]	SCT-150-SS-A-S
Clear	2.0ml	Conical	Pre-Sterilized [‡]	SCT-200-C-S
Assorted [†]	2.0ml	Conical	Pre-Sterilized [‡]	SCT-200-A-S
Clear	2.0ml	Self Standing	Pre-Sterilized [‡]	SCT-200-SS-C-S
Assorted [†]	2.0ml	Self Standing	Pre-Sterilized [‡]	SCT-200-SS-A-S




All pre-sterilized screw cap tubes are packaged with non-looped caps attached and include 'O' rings

Packaging Description: 500 tubes/pack, 8 packs/case; [†]50/bag, 10 bags/pack, 8 packs/case
Maximum G-force: 20,000 RCF

Tubes only

Clear	0.5ml	Conical	Screw Cap Tubes	ST-45
Clear	0.5ml	Conical	Screw Cap Tubes	ST-050
Clear	0.5ml	Self Standing	Screw Cap Tubes	ST-050-SS
Clear	1.5ml	Conical	Screw Cap Tubes	ST-150
Clear	1.5ml	Self Standing	Screw Cap Tubes	ST-150-SS
Clear	2.0ml	Conical	Screw Cap Tubes	ST-200
Clear	2.0ml	Self Standing	Screw Cap Tubes	ST-200-SS
Amber	0.5ml	Conical	Screw Cap Tubes	ST-050-X
Amber	0.5ml	Self Standing	Screw Cap Tubes	ST-050-SS-X
Amber	1.5ml	Conical	Screw Cap Tubes	ST-150-X
Amber	1.5ml	Self Standing	Screw Cap Tubes	ST-150-SS-X
Amber	2.0ml	Conical	Screw Cap Tubes	ST-200-X
Amber	2.0ml	Self Standing	Screw Cap Tubes	ST-200-SS-X



Screw Cap Tubes (Tubes Only)
 Amber Screw Cap Tubes (Tubes Only)

Packaging Description: 500 tubes/pack, 8 packs/case
Maximum G-force: 20,000 RCF

[†] Maximum Recovery Screw Cap Tubes available as special order items. [‡] Individual Colors Available as Special Order Items.
 Please call our Customer Service department for more information.



// 74

Microtube CAPS

Regular and Looped

Caps without 'O' rings



Cap Color

- Clear
- Assorted¹
- Grey
- White
- Pink
- Red
- Orange
- Yellow
- Green
- Blue
- Violet
- Brown

Description

- Non-Looped Screw Caps
- Non-Looped Screw Caps
- Non-Looped Screw Caps
- Non-Looped Screw Caps
- Non-Looped Screw Caps
- Non-Looped Screw Caps
- Non-Looped Screw Caps
- Non-Looped Screw Caps
- Non-Looped Screw Caps
- Non-Looped Screw Caps
- Non-Looped Screw Caps
- Non-Looped Screw Caps
- Non-Looped Screw Caps

Axygen Cat #

- SC-C
- SC-A
- SC-GR
- SC-W
- SC-P
- SC-R
- SC-O
- SC-Y
- SC-G
- SC-B
- SC-V
- SC-BR

Packaging Description:

500 caps/pack, 8 packs/case; ¹50/bag, 10 bags/pack, 8 packs/case

Looped caps without 'O' rings



- Clear
- Assorted¹
- Grey
- White
- Pink
- Red
- Orange
- Yellow
- Green
- Blue
- Violet
- Brown

- Looped Screw Caps
- Looped Screw Caps
- Looped Screw Caps
- Looped Screw Caps
- Looped Screw Caps
- Looped Screw Caps
- Looped Screw Caps
- Looped Screw Caps
- Looped Screw Caps
- Looped Screw Caps
- Looped Screw Caps
- Looped Screw Caps
- Looped Screw Caps

- SC-LP-C
- SC-LP-A
- SC-LP-GR
- SC-LP-W
- SC-LP-P
- SC-LP-R
- SC-LP-O
- SC-LP-Y
- SC-LP-G
- SC-LP-B
- SC-LP-V
- SC-LP-BR

Packaging Description:

500 caps/pack, 8 packs/case; ¹50/bag, 10 bags/pack, 8 packs/case

Caps with 'O' rings



- Clear
- Assorted¹
- Grey
- White
- Pink
- Red
- Orange
- Yellow
- Green
- Blue
- Violet
- Brown

- Non-Looped Screw Caps with 'O' rings
- Non-Looped Screw Caps with 'O' rings
- Non-Looped Screw Caps with 'O' rings
- Non-Looped Screw Caps with 'O' rings
- Non-Looped Screw Caps with 'O' rings
- Non-Looped Screw Caps with 'O' rings
- Non-Looped Screw Caps with 'O' rings
- Non-Looped Screw Caps with 'O' rings
- Non-Looped Screw Caps with 'O' rings
- Non-Looped Screw Caps with 'O' rings
- Non-Looped Screw Caps with 'O' rings
- Non-Looped Screw Caps with 'O' rings
- Non-Looped Screw Caps with 'O' rings

- SCO-C
- SCO-A
- SCO-GR
- SCO-W
- SCO-P
- SCO-R
- SCO-O
- SCO-Y
- SCO-G
- SCO-B
- SCO-V
- SCO-BR

Packaging Description:

500 caps/pack, 8 packs/case; ¹50/bag, 10 bags/pack, 8 packs/case

Looped caps with 'O' rings



- Clear
- Assorted¹
- Grey
- White
- Pink
- Red
- Orange
- Yellow
- Green
- Blue
- Violet
- Brown

- Looped Screw Caps with 'O' rings
- Looped Screw Caps with 'O' rings
- Looped Screw Caps with 'O' rings
- Looped Screw Caps with 'O' rings
- Looped Screw Caps with 'O' rings
- Looped Screw Caps with 'O' rings
- Looped Screw Caps with 'O' rings
- Looped Screw Caps with 'O' rings
- Looped Screw Caps with 'O' rings
- Looped Screw Caps with 'O' rings
- Looped Screw Caps with 'O' rings
- Looped Screw Caps with 'O' rings
- Looped Screw Caps with 'O' rings

- SCO-LP-C
- SCO-LP-A
- SCO-LP-GR
- SCO-LP-W
- SCO-LP-P
- SCO-LP-R
- SCO-LP-O
- SCO-LP-Y
- SCO-LP-G
- SCO-LP-B
- SCO-LP-V
- SCO-LP-BR

Packaging Description:

500 caps/pack, 8 packs/case; ¹50/bag, 10 bags/pack, 8 packs/case