



# Pipette Tip COMPATIBILITY

Axygen Axyjet (single channel)				(multi-channel)	Axygen Axyjet-16			Biohit Proline (single channel)				(multi-channel)	Biohit e-Line (single channel)			(multi-channel)	Biohit m-Line (s/ch)										
0.1-2µl	0.5-10µl	2-20µl	5-50µl	10-100µl	20-200µl	50-250µl	100-1000µl	1-10µl	5-50µl	20-200µl	50-300µl	0.2-2µl	0.5-10µl	5-50µl	0.2-10µl	5-100µl	25-250µl	0.2-10µl	5-120µl	10-300µl	50-1000µl	0.2-10µl	5-120µl	10-300µl	0.5-10µl	2-20µl	10-100µl

< Pipettor Brand/make

< Volume Range

## NON FILTER TIPS

																		# ebed	<b>Volume/Axygen Tip</b>
																		14	0.5-60µl, T-384
																		14	0.5-30µl, T-1536
																		10	0.5-10µl, T-300
																		10	0.5-10µl, TXL-10
																		10	0.5-20µl, T-400
																		11	1-200µl, T-200
																		11	1-200µl, TE-204
																		12	1-200µl, T-205-WB
																		11	1-200µl, T-210
																		11	1-200µl, TR-222
																		12	1-250µl, TR-333
																		12	1-300µl, T-350
																		12	100-1000µl, T-1000
13	100-1000µl, TE-1004																		
13	100-1000µl, T-1005-WB																		
13	100-1200µl, T-1200																		
13	15-1250µl, MTX-1250																		

## AEROSOL BARRIER TIPS

																		# ebed	<b>Volume/Axygen Tip</b>
																		16	0.5-30µl, TF-384-30
																		16	0.5-10µl, TF-300
																		16	0.5-10µl, TXLF-10
																		16	0.5-10µl, TF-400
																		16	0.5-20µl, TF-420
																		17	0.5-20µl, TF-20
																		17	0.5-50µl, TF-50
																		17	1-100µl, TF-100
																		17	1-150µl, TF-150
																		17	1-200µl, TF-200
																		17	1-200µl, TF-205-WB
																		18	1-300µl, TF-350
																		18	10-500µl, TF-500
18	100-1000µl, TF-1000																		
18	100-1000µl, TF-1005-WB																		
18	100-1100µl, TF-1200																		
18	15-1250µl, MTXF-1250																		

FOR TECHNICAL INFORMATION, VISIT [WWW.AXYGEN.COM](http://WWW.AXYGEN.COM)

# Pipette Tip COMPATIBILITY



Pipettor Brand/make >

Volume Range >

Volume/Axygen Tip

page #

Biohit m-Line (single) (multi-ch)		BrandTech Transferpette (single channel)							BrandTech Transferpette (multi-channel)						Capp Aero (single/multi-channel)					(384 format multi-channel)			Eppendorf Electronic								
20-200µl	100-1000µl	0.1-1µl	0.5-10µl	2-20µl 'Ultra Micro'	2-20µl	5-50µl	10-100µl	20-200µl	25-250µl	100-1000µl	0.5-10µl	2-20µl	2.5-25µl	5-50µl	10-100µl	20-200µl	30-300µl	0.2-2µl	0.5-10µl	5-50µl	10-100µl	25-200µl	50-300µl	100-1000µl	0.2-2µl	0.5-10µl	5-50µl	0.5-10µl	5-100µl	2-200µl	50-1000µl


## NON FILTER TIPS

0.5-60µl, T-384	14
0.5-30µl, T-1536	14
0.5-10µl, T-300	10
0.5-10µl, TXL-10	10
0.5-20µl, T-400	10
1-200µl, T-200	11
1-200µl, TE-204	11
1-200µl, T-205-WB	12
1-200µl, T-210	11
1-200µl, TR-222	11
1-250µl, TR-333	12
1-300µl, T-350	12
100-1000µl, T-1000	12
100-1000µl, TE-1004	13
100-1000µl, T-1005-WB	13
100-1200µl, T-1200	13
15-1250µl, MTX-1250	13

## AEROSOL BARRIER TIPS

0.5-30µl, TF-384-30	16
0.5-10µl, TF-300	16
0.5-10µl, TXLF-10	16
0.5-10µl, TF-400	16
0.5-20µl, TF-420	16
0.5-20µl, TF-20	17
0.5-50µl, TF-50	17
1-100µl, TF-100	17
1-150µl, TF-150	17
1-200µl, TF-200	17
1-200µl, TF-205-WB	17
1-300µl, TF-350	18
10-500µl, TF-500	18
100-1000µl, TF-1000	18
100-1000µl, TF-1005-WB	18
100-1100µl, TF-1200	18
15-1250µl, MTFX-1250	18





# Pipette Tip COMPATIBILITY



Pipettor Brand/make >

Volume Range >

Volume/Axygen Tip

page #

Labnet Labpette (single channel)						Labnet Labpette (multi-channel)				Labnet BioPette E (single channel)			Labnet BioPette E (multi-channel)			Matrix Impact/Impact <sup>2</sup> (s/mc)			Nichiryo/Oxford Benchmate				Rainin EDP3 (single channel)								
0.1-2µl	0.5-10µl	2-20µl	10-100µl	20-200µl	100-1000µl	1-10µl	5-50µl	20-200µl	50-300µl	0.5-10µl	2-20µl	10-200µl	100-1000µl	0.5-10µl	2-20µl	10-200µl	100-1200µl	2-125µl	5-250µl	15-1250µl	0.5-10µl	2-20µl	10-50µl	10-100µl	20-200µl	100-1000µl	E-MC-10, 0.5-10µl	EZ-25, 2.5-25µl	EZ-100, 10-100µl	EZ-250, 25-250µl	EZ-1000, 100-1000µl

## NON FILTER TIPS

0.5-60µl, T-384	14
0.5-30µl, T-1536	14
0.5-10µl, T-300	10
0.5-10µl, TXL-10	10
0.5-20µl, T-400	10
1-200µl, T-200	11
1-200µl, TE-204	11
1-200µl, T-205-WB	12
1-200µl, T-210	11
1-200µl, TR-222	11
1-250µl, TR-333	12
1-300µl, T-350	12
100-1000µl, T-1000	12
100-1000µl, TE-1004	13
100-1000µl, T-1005-WB	13
100-1200µl, T-1200	13
15-1250µl, MTX-1250	13

Labnet Labpette (single channel)						Labnet Labpette (multi-channel)				Labnet BioPette E (single channel)			Labnet BioPette E (multi-channel)			Matrix Impact/Impact <sup>2</sup> (s/mc)			Nichiryo/Oxford Benchmate				Rainin EDP3 (single channel)							
0.1-2µl						0.5-10µl				2-20µl			10-100µl			20-200µl			100-1000µl											
1-10µl						5-50µl				20-200µl			50-300µl			0.5-10µl			2-20µl			10-50µl				10-100µl				
20-200µl						100-1000µl				0.5-10µl			2-20µl			10-200µl			100-1200µl			2-125µl				5-250µl				
15-1250µl						0.5-10µl				2-20µl			10-50µl			10-100µl			20-200µl			100-1000µl								
E-MC-10, 0.5-10µl						EZ-25, 2.5-25µl				EZ-100, 10-100µl			EZ-250, 25-250µl			EZ-1000, 100-1000µl														

## AEROSOL BARRIER TIPS

0.5-30µl, TF-384-30	16
0.5-10µl, TF-300	16
0.5-10µl, TXLF-10	16
0.5-10µl, TF-400	16
0.5-20µl, TF-420	16
0.5-20µl, TF-20	17
0.5-50µl, TF-50	17
1-100µl, TF-100	17
1-150µl, TF-150	17
1-200µl, TF-200	17
1-200µl, TF-205-WB	17
1-300µl, TF-350	18
10-500µl, TF-500	18
100-1000µl, TF-1000	18
100-1000µl, TF-1005-WB	18
100-1100µl, TF-1200	18
15-1250µl, MTXF-1250	18

Labnet Labpette (single channel)						Labnet Labpette (multi-channel)				Labnet BioPette E (single channel)			Labnet BioPette E (multi-channel)			Matrix Impact/Impact <sup>2</sup> (s/mc)			Nichiryo/Oxford Benchmate				Rainin EDP3 (single channel)							
0.1-2µl						0.5-10µl				2-20µl			10-100µl			20-200µl			100-1000µl											
1-10µl						5-50µl				20-200µl			50-300µl			0.5-10µl			2-20µl			10-50µl				10-100µl				
20-200µl						100-1000µl				0.5-10µl			2-20µl			10-200µl			100-1200µl			2-125µl				5-250µl				
15-1250µl						0.5-10µl				2-20µl			10-50µl			10-100µl			20-200µl			100-1000µl								
E-MC-10, 0.5-10µl						EZ-25, 2.5-25µl				EZ-100, 10-100µl			EZ-250, 25-250µl			EZ-1000, 100-1000µl														



