

## Mini Incubator and Mini Water Bath

The Mini Incubator and Mini Water Bath are small units designed for personal use, small laboratories and classrooms. Each features a broad temperature range to meet a variety of small sample applications.

Wrap-around heating elements ensure even heating and uniformity within the Mini Incubator. A hydraulic thermostat controls heating and temperature stability. A thermometer port in the top of the unit allows users to monitor temperature. The large window in the door gives a clear view of the contents which can include up to 48 petri dishes or 28 standard microplates. Shelves can be placed at two points within the corrosion resistant chamber. A cord port allows small mixers to be placed in the incubator.

The Mini Water Bath has a 6 liter capacity to accommodate test tube racks or up to three 250 ml flasks. A powerful 550W heater allows for quick temperature adjustment. A hydraulic thermostat ensures accurate temperature control. The seamless, stainless steel chamber is corrosion resistant and easy to clean. A stainless steel gable cover helps to retain heat in the bath and control evaporation. A clip is included for suspending a thermometer in the bath.



### Specifications

	Mini Incubator	Mini Water Bath
Temperature range	Ambient +5° to 60°C	Ambient +5° to 100°C
Capacity	9.2 liters/0.375 cu. ft.	6 liter
Controls	Analog	Analog
Exterior dimensions (W x D x H)	28.5 x 28 x 33.5 cm 11.2 x 11 x 13.2 in.	36 x 28 x 32 cm 14 x 11 x 13 in.
Interior dimensions (W x D x H)	23 x 20 x 20 cm 9 x 7.9 x 7.9 in.	32 x 17 x 18 cm 12.6 x 6.7 x 6.9 in.
Weight	8.3 kg/ 18.4 lb	7.7 kg/ 17 lb
Electrical	120V or 230V, 50/60 Hz	120V or 230V, 50/60 Hz

### Ordering Information

SBI510	Mini Incubator, 120V
SBI510-2	Mini Incubator, 230V, EU power cord
SBI510-3	Mini Incubator, 230V, UK power cord
SBI10	Shelf for Mini Incubator, each
SBW110	Mini Water Bath with stainless steel gable cover and thermometer clip, 120V
SBW110-2	Mini Water Bath with stainless steel gable cover and thermometer clip, 230V, EU power cord
SBW110-3	Mini Water Bath with stainless steel gable cover and thermometer clip, 230V, UK power cord
W0020-110C	Non-toxic, Teflon coated thermometer for -20° to 110°C, 1° resolution, 305 mm long, 76 mm immersion