



SERANA (WA) PTY LTD
ACN: 128 305 635

Lot 535, 6 Juniper Way
Davenport, Bunbury WA 6230
P.O. Box 418, Bunbury WA 6231

Phone: (08) 9725 6868
Fax: (08) 9725 6867
Email: info@serana-australia.com
Website: www.serana-australia.com

Sample Certificate of Analysis

Product: **Secoll lymphocyte separating solution**
Density: 1.077 g/ml
Batch no.:
Production date:
Expiry date:
Storage: +2°C to RT

Specifications:

Visual: Clear colorless liquid free of particulate matter
pH: 6.5 to 7.5
Osmolality: 280 to 310 mOsm/kg
Density: 1.076 to 1.078 g/ml @ 20°C
Endotoxin: < 10 EU/ml
Sterility: Incubation 37°C: negative
Incubation 20°C: negative

SAMPLE