

T-Series Products

Cellu-Sep® Regenerated Cellulose Tubular Membranes have proven time after time to be the smartest value in dialysis tubing available to laboratories anywhere around the globe.

Cellu-Sep is low in price and includes extra value items that don't cost you extra. Each box of T1, T2 & T3 includes a technical manual with instructions, analytical data sheet and disposable scissors.

Cellu-Sep®
T-SERIES
DIALYSIS MEMBRANES ▶



with free dialysis accessories in every box



You also get two free Locking Dialysis Membrane Clamps with every box - just to help make your dialysis procedures run a little bit easier. And to protect your order of Cellu-Sep, we always ship each box in a sealed bag with a desiccant to control excess humidity.

Cellu-Sep T1/Nominal MWCO: 3,500

Part No.	Flat Width	Ø Dry	Vol/cm	Length
5015-19	19	12.1 mm	1.15	15 M
5015-46	46	29.3 mm	6.74	15 M
5015-55	55	35.0 mm	9.62	15 M
5030-19	19	12.1 mm	1.15	30 M
5030-46	46	29.3 mm	6.74	30 M

Cellu-Sep T2/Nominal MWCO: 6,000 - 8,000

Part No.	Flat Width	Ø Dry	Vol/cm	Length
8015-10	10	6.37 mm	.32	15 M
8015-100	100	63.7 mm	31.8	15 M
8015-23	23	14.6 mm	1.67	15 M
8015-25	25	15.9 mm	1.98	15 M
8015-32	32	20.4 mm	3.27	15 M
8015-40	40	25.5 mm	5.10	15 M
8015-50	50	31.8 mm	7.94	15 M
8030-10	10	6.37 mm	.32	30 M
8030-23	23	14.6 mm	1.67	30 M
8030-25	25	15.9 mm	1.98	30 M
8030-32	32	20.4 mm	3.27	30 M
8030-40	40	25.5 mm	5.10	30 M
8030-50	50	31.8 mm	7.94	30 M

Cellu-Sep T3/Nominal MWCO: 12,000 - 14,000

Part No.	Flat Width	Ø Dry	Vol/cm	Length
1205-76	76	48.4 mm	18.4	5 M
1210-76	76	48.4 mm	18.4	10 M
1215-10	10	6.37 mm	.32	15 M
1215-25	25	15.9 mm	1.98	15 M
1215-33	25	15.9 mm	1.98	15 M
1215-45	45	28.6 mm	6.42	15 M
1230-10	10	6.37 mm	.32	30 M
1230-25	25	15.9 mm	1.98	30 M
1230-33	33	21.0 mm	3.46	30 M
1230-45	45	28.6 mm	6.42	30 M

Cellu-Sep T4/Nominal MWCO: 12,000 - 14,000

Standard Grade, includes all extra value items except data sheet

Part No.	Flat Width	Ø Dry	Vol/cm	Length
1415-75	75	47.7 mm	17.9	15 M
1430-10	10	6.37 mm	.32	30 M
1430-25	25	15.9 mm	1.98	30 M
1430-45	45	28.6 mm	6.42	30 M
1430-33	33	21.0 mm	3.46	30 M