

HUGE SAVINGS on Electrophoresis Systems

MultiSUB™ Mini Duo System

A small lab bench footprint

- Run up to 64 samples
- Economical low gel and buffer volumes
- Supplied with 7x7cm and 7x10cm trays

\$520 (ex-GST)




MultiSUB™ Midi Duo System

Ideal for rapid electrophoresis

- Run up to 100 samples
- Economical low buffer volumes
- Supplied with 10x7cm and 10x10cm trays

\$565 (ex-GST)



omniPAGE™ Mini Vertical Protein System

Rapid gel casting and loading

- Run up to 4 gels / 80 samples
- Low buffer volumes
- Compatible with all 8x10cm and 10x10cm precast gels

\$945 (ex-GST)



CATALOG No	DESCRIPTION	QTY.	PRICE (ex-GST)
OLMSMiniDuo-CP	MultiSUB Mini Duo Horizontal System, 7x7cm and 7x10cm UV Tray, 2x8-sample, 1mm thick combs, casting dams & loading guides	1	\$520
OLMSMidiDuo-CP	MultiSUB Midi Duo Horizontal System, 10x7cm and 10x10cm UV Tray, 2x16-sample, 1mm thick combs casting dams & loading guides	1	\$565
OLVS10DSYS-CP	omniPAGE Mini Vertical Protein System, complete system for mini vertical electrophoresis	1	\$945

Terms & Conditions: All prices quoted in AUD and exclude GST. Offers exclude servicing, spare parts and freight. Unless otherwise stated, accessories shown are extra. Offer valid for stock on hand. Fisher Biotech reserves the right to withdraw offers at any time.