

## Material Safety Data Sheet

### FavorPrep™ Tissue Genomic DNA Extraction Mini Kit

<b>Cat. No.:</b> <b>Contents</b>	<b>FATGK 000</b> <b>(4 preps)</b>	<b>FATGK 001</b> <b>(50 preps)</b>	<b>FATGK 001-1</b> <b>(100 preps)</b>	<b>FATGK 001-2</b> <b>(300 preps)</b>
FATG1 Buffer	1.5 ml	15 ml	30 ml	70 ml
FATG2 Buffer	1.5 ml	15 ml	30 ml	70 ml
W1 Buffer (Concentrate)	1.3 ml	22 ml	44 ml	124 ml
Wash Buffer (Concentrate)	1 ml	10 ml	20 ml	55 ml
Elution Buffer	1 ml	15 ml	30 ml	90 ml
Proteinase K (Liquid)	100 µl	1050 µl	1050 µl × 2	1600 µl × 4

**Favorgen Biotech Corp.**

Address: Ping-Tung Agricultural Biotechnology Park  
No. 8-1, Yuanxi 2<sup>nd</sup> Road, Ping Tung 908126,  
TAIWAN

Phone Number: +886-8-762-1829

technical@favorgen.com

www.favorgen.com

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



## Section 1. IDENTIFICATION

### 1.1 Product Identifier

Commercial Product Name: **FATG1 Buffer**

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

Product form: Laboratory chemical mixture/Lysis Buffer

Relevant identified uses: For research use only, not for diagnostic use.

Uses advised against:

-----

### 1.3 Details of the supplier of the safety data sheet

Company: Favorgen Biotech Corp.

Address: Ping-Tung Agricultural Biotechnology Park No. 8-1, Yuanxi 2<sup>nd</sup> Road, Ping Tung 908126, TAIWAN

Phone Number: +886-8-762-1829

### 1.4 Emergency telephone number

+886-8-762-1829

## Section 2. HAZARDS IDENTIFICATION

### 2.1 Classification of the substance or mixture

Hazard identification Hazard classes/categories

Skin corrosion/irritation	Category 2
Serious eye damage/eye irritation	Category 2A

GHS hazard pictograms:



GHS07

Signal word: **WARNING**

Hazard statement(s)

H315	Causes skin irritation
H319	Causes serious eye irritation

Precautionary statement(s)

P264	Wash hands thoroughly after handling
P280	Wear protective gloves/protective clothing/eye protection/face protection
P302 + P352	IF ON SKIN: Wash with plenty of soap and water
P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
P332 + P313	If skin irritation occurs: Get medical advice/attention
P337 + P313	If eye irritation persists: Get medical advice/attention
P362 + P364	Take off contaminated clothing and wash it before reuse.

### 2.2 Label elements

According to section 1.5.2 of Annex I to CLP

GHS symbol



GHS07

WARNING

### 2.3 Other hazards

Not applicable

## Section 3. COMPOSITION/INFORMATION ON INGREDIENTS

### 3.1 Substances or 3.2 Mixtures

Chemical: **Sodium dodecyl sulfate**  
Concentration: 1~5%

CAS No.: 151-21-3

Always wear recommended Personal Protective Equipment.

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



## Section 4. FIRST AID MEASURES

### 4.1 Description of first-aid measures

#### If inhaled

After inhalation: fresh air. If breathing stops: mouth-to-mouth breathing or artificial respiration. Oxygen if necessary. Immediately call in physician.

#### In case of skin contact

After skin contact: Take off immediately all contaminated clothing. Rinse skin with water/shower.

#### In case of eye contact

After eye contact: rinse out with plenty of water. Call in ophthalmologist. Remove contact lenses.

#### If swallowed

After swallowing: immediately make victim drink water (two glasses at most). Consult a physician.

#### If in eyes

Rinse cautiously with water for several minutes

### 4.2 Most important symptoms and effects, both acute and delayed

H315 - Causes skin irritation

H319 - Causes serious eye irritation

### 4.3 Indication of any immediate medical attention and special treatment needed

Wear Personal Protective Equipment. If skin irritation occurs: Get medical advice

## Section 5. FIRE-FIGHTING MEASURES

### 5.1 Extinguishing media

All extinguishers can be used.

### 5.2 Special hazards arising from the substance or mixture

None

### 5.3 Advice for firefighters

Wear self-contained breathing equipment and protective suit.

### 5.4 Additional Information

----

## Section 6. ACCIDENTAL RELEASE MEASURE

### 6.1 Personal precautions, protective equipment, and emergency procedures

Do not breathe vapours/Ensure adequate ventilation/Avoid contact with skin, eyes or clothing/Use personal protection equipment/Regular staff training is necessary.

### 6.2 Environmental precautions

Do not allow mixture to enter ground water system.

### 6.3 Methods and material for containment and cleaning up

Soak up with inert absorbent material.

### 6.4 Reference to other sections

---

## Section 7. HANDLING AND STORAGE

### 7.1 Precautions for safe handling

Always wear recommended Personal Protective Equipment.

### 7.2 Conditions for safe storage, including any incompatibilities

Keep in a dry, cool, and well-ventilated place.

### 7.3 Specific end use(s)

For research use only.

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



## Section 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

### 8.1 Control parameters: None

Chemical: Sodium dodecyl sulfate  
Concentration: 1~5%

CAS No.: 151-21-3

### 8.2 Exposure controls

Good ventilation and exhaust system in the workplace, and the floor has drainage and washing facilities, which can resist the erosion of chemicals. Must maintain the highest level of cleanliness in the workplace.

#### Personal protective equipment

**Eye/face protection:** Safety glasses. Used to prevent splash hazards. Government standards such as NIOSH (US) or EN 166(EU).

**Skin protection:** Recommended to avoid contamination with these hazards.

**Hand protection:** Protective gloves, use for short times chemical resistant latex gloves with code EN 374-3 level 1.

**Respiratory protection:** Not necessary, if the workplace is well-ventilated.

**Personal Hygiene:** Wash hands before breaks and at the end of work. Handle in accordance with good industrial hygiene and safety practice.

## Section 9. PHYSICAL AND CHEMICAL PROPERTIES

### 9.1 Information on basic physical and chemical properties

Appearance: liquid/Color: clear/Odor: characteristic /pH: 8.0~8.6

### 9.2 Other information

---

## Section 10. STABILITY AND REACTIVITY

### 10.1 Reactivity

No data available.

### 10.2 Chemical stability

Stable under normal condition.

### 10.3 Possibility of hazardous reactions

None.

### 10.4 Conditions to avoid

No data available.

### 10.5 Incompatible materials

Not dangerous reaction known under normal use.

### 10.6 Hazardous decomposition products

No decomposition if stored and applied as directed.

## Section 11. TOXICOLOGICAL INFORMATION

### 11.1 Information on toxicological effects

#### Acute toxicity

Not Classified

Chemical Name	Oral LD50 (Rat)	Dermal LD50 (Rabbit)	Inhalation LC50 (Rat)
Sodium dodecyl sulfate CAS No.: 151-21-3	1,228 mg/kg	200 mg/kg	3,900 mg/m <sup>3</sup>

**Skin corrosion/irritation**

Not Classified

**Serious eye damage/irritation**

Not Classified

**Respiratory or skin sensitization**

Not Classified

**Germ cell mutagenicity**

Not Classified

**Carcinogenicity**

Not Classified

**Reproductive toxicity**

Not Classified

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



Specific target organ toxicity (STOT) – single exposure	Not Classified
Specific target organ toxicity (STOT) – repeated exposure	Not Classified
Aspiration hazard	Not Classified

## 11.2 Information on other hazards

No data available

## Section 12. ECOLOGICAL INFORMATION

### 12.1 Toxicity

Following information is valid for pure substances.

Sodium dodecyl sulfate

CAS No.: 151-21-3

EC50 Daphnia magna /48h: 1.8 mg/L

### 12.2 Persistence and degradability

No data available.

### 12.3 Bioaccumulative potential

No data available.

### 12.4 Mobility in soil

No data available.

### 12.5 Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be a PBT or a vPvB.

### 12.6 Endocrine disrupting properties

No data available.

### 12.7 Other adverse effects

No data available.

## Section 13. DISPOSAL CONSIDERATIONS

Refer to local regulations.

### 13.1 Waste treatment methods

Dispose of contents/container in accordance with local regulation.

Empty containers should be taken to an approved waste handling site for recycling or disposal.

## Section 14. TRANSPORT INFORMATION

### 14.1 UN number

ADR/RID: -

IMDG: -

IATA: -

### 14.2 UN proper shipping name

ADR/RID: Not dangerous goods

IMDG: Not dangerous goods

IMDG: Not dangerous goods

### 14.3 Transport hazard class(es)

ADR/RID: -

IMDG: -

IATA: -

### 14.4 Packaging group

ADR/RID: -

IMDG: -

IATA: -

### 14.5 Environmental hazards

ADR/RID: no

IMDG Marine pollutant: no

IATA-DGR: no

### 14.6 Special precautions for user

None

### 14.7 Maritime transport in bulk according to IMO instruments

Not Applicable.

Further information

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



Not classified as dangerous in transport regulations.

## Section 15. REGULATIONS INFORMATION

### 15.1 Safety, health, and environmental regulations/legislation specific for the substance or mixture

Look for your country-specific regulations.

### 15.2 Chemical safety assessment

Not necessary

## Section 16. OTHER INFORMATION

### 16.1 List of R, H and P phrases

Full text of H-Statements referred to under sections 2 and 3.

### 16.2 Training Advice

Regular safety training.

### 16.3 Recommended Restriction on Use

Only for professional/research user.

## Section 1. IDENTIFICATION

### 1.1 Product Identifier

Commercial Product Name: **FATG2 Buffer**

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

Product form: Laboratory chemical mixture/Lysis Buffer

Relevant identified uses: For research use only, not for diagnostic use.

Uses advised against:

-----

### 1.3 Details of the supplier of the safety data sheet

Company: Favorgen Biotech Corp.

Address: Ping-Tung Agricultural Biotechnology Park No. 8-1, Yuanxi 2<sup>nd</sup> Road, Ping Tung 908126, TAIWAN

Phone Number: +886-8-762-1829

### 1.4 Emergency telephone number

+886-8-762-1829

## Section 2. HAZARDS IDENTIFICATION

### 2.1 Classification of the substance or mixture

Hazard identification Hazard classes/categories

Acute oral toxicity	Category 4
Skin corrosion/irritation	Category 2
Serious eye damage/eye irritation	Category 2

GHS hazard pictograms:



GHS07

Signal word: **WARNING**

Hazard statement(s)

H302

Harmful if swallowed

H315

Causes skin irritation

H319

Causes serious eye irritation.

H412

Harmful to aquatic life with long lasting effects

Precautionary statement(s)

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



P264	Wash hands thoroughly after handling
P280	Wear protective gloves/protective clothing/eye protection/face protection
P270	Do not eat, drink or smoke when using this product
P273	Avoid release to the environment
P301 + P310	IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician
P302 + P352	IF ON SKIN: Wash with plenty of soap and water
P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
P330	Rinse mouth
P332 + P313	If skin irritation occurs: Get medical advice/attention

## 2.2 Label elements

According to section 1.5.2 of Annex I to CLP

GHS symbol



GHS07

WARNING

## 2.3 Other hazards

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

## Section 3. COMPOSITION/INFORMATION ON INGREDIENTS

### 3.1 Substances or 3.2 Mixtures

Chemical:	<b>Guanidine hydrochloride</b>	CAS No.: 50-01-1
Concentration:	36 - <50%	
Classification:	H302, Acute oral Tox. 4 H315, Skin Irritation 2 H319, Eye Irritation 2	

Chemical:	<b>Polyoxyethylene sorbitan monolaurate (Tween® 20)</b>	CAS No.: 9005-64-5
Concentration:	10 - <33%	
Classification:	No criteria for classification or naming of chemical not required.	

## Section 4. FIRST AID MEASURES

### 4.1 Description of first-aid measures

#### If inhaled

After inhalation: fresh air. If breathing stops: mouth-to-mouth breathing or artificial respiration. Oxygen if necessary. Immediately call in physician.

#### In case of skin contact

After skin contact: Take off immediately all contaminated clothing. Rinse skin with water/shower.

#### In case of eye contact

After eye contact: rinse out with plenty of water. Call in ophthalmologist. Remove contact lenses.

#### If swallowed

After swallowing: immediately make victim drink water (two glasses at most). Consult a physician.

### 4.2 Most important symptoms and effects, both acute and delayed

H302 - Harmful if swallowed  
H315 - Causes skin irritation  
H319 - Causes serious eye irritation  
H412 - Harmful to aquatic life

### 4.3 Indication of any immediate medical attention and special treatment needed

No data available

## Section 5. FIRE-FIGHTING MEASURES

### 5.1 Extinguishing media

All extinguishers can be used.

### 5.2 Special hazards arising from the substance or mixture

None known.

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



## 5.3 Advice for firefighters

No specific fire or explosion hazard.

## 5.4 Additional Information

---

## Section 6. ACCIDENTAL RELEASE MEASURE

### 6.1 Personal precautions, protective equipment, and emergency procedures

Do not breathe vapours/Ensure adequate ventilation/Avoid contact with skin, eyes or clothing/Use personal protection equipment/Regular staff training is necessary.

### 6.2 Environmental precautions

Do not allow mixture to enter ground water system.

### 6.3 Methods and material for containment and cleaning up

Soak up with inert absorbent material.

### 6.4 Reference to other sections

---

## Section 7. HANDLING AND STORAGE

### 7.1 Precautions for safe handling

Always wear recommended Personal Protective Equipment.

### 7.2 Conditions for safe storage, including any incompatibilities

Keep in a dry, cool, and well-ventilated place.

### 7.3 Specific end use(s)

For research use only.

## Section 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

### 8.1 Control parameters: None

Chemical: guanidine hydrochloride

CAS No.: 50-01-1

NIOSH: not listed

OSHA: not listed

Chemical: polyoxyethylene sorbitan monolaurate

CAS No.: 9005-64-5

### 8.2 Exposure controls

Good ventilation and exhaust system in the workplace, and the floor has drainage and washing facilities, which can resist the erosion of chemicals. Must maintain the highest level of cleanliness in the workplace.

#### Personal protective equipment

**Eye/face protection:** Safety glasses. Used to prevent splash hazards. Government standards such as NIOSH (US) or EN 166(EU).

**Skin protection:** Recommended to avoid contamination with these hazards.

**Hand protection:** Protective gloves, use for short times chemical resistant latex gloves with code EN 374-3 level 1.

**Respiratory protection:** Not necessary, if the workplace is well-ventilated.

**Personal Hygiene:** Wash hands before breaks and at the end of work. Handle in accordance with good industrial hygiene and safety practice.

## Section 9. PHYSICAL AND CHEMICAL PROPERTIES

### 9.1 Information on basic physical and chemical properties

Appearance: liquid/Color: slightly yellow/Odor: odorless/pH: 3.4~5.6/Specific gravity: 1.15 g/cm<sup>3</sup>

### 9.2 Other information

---



# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



## Section 10. STABILITY AND REACTIVITY

### 10.1 Reactivity

No data available.

### 10.2 Chemical stability

No data available.

### 10.3 Possibility of hazardous reactions

**Note:** Can form very reactive substances with oxidizing agents.

### 10.4 Conditions to avoid

Not necessary.

### 10.5 Incompatible materials

**Note:** Avoid contact with strong acids or alkalines.

### 10.6 Hazardous decomposition products

Decompositions are not observed during the expiration period under recommended conditions.

## Section 11. TOXICOLOGICAL INFORMATION

### 11.1 Information on toxicological effects

#### Acute toxicity

Not Classified

Chemical Name	Oral LD50 (Rat)	Dermal LD50 (Rabbit)	Inhalation LC50 (Rat)
Guanidine hydrochloride CAS No.: 50-01-1	475-907 mg/kg	2000 mg/kg	[4h] 3181-7655 µg/m <sup>3</sup>
Polyoxyethylene sorbitan monolaurate CAS No.: 9005-64-5	No data available	No data available	No data available

#### Skin corrosion/irritation

Skin irritation/ Irritating to eyes

#### Serious eye damage/irritation

Not Classified

#### Respiratory or skin sensitization

Not Classified

#### Germ cell mutagenicity

Not Classified

#### Carcinogenicity

Not Classified

#### Reproductive toxicity

Not Classified

#### Specific target organ toxicity (STOT) – single exposure

Not Classified

#### Specific target organ toxicity (STOT) – repeated exposure

Not Classified

#### Aspiration hazard

Not Classified

### 11.2 Information on other hazards

No data available.

## Section 12. ECOLOGICAL INFORMATION

### 12.1 Toxicity

**Following information is valid for pure substances.**

Chemical: guanidine hydrochloride

CAS No.: 50-01-1

LC50leuciscus idus/96h: 1759 mg/L

WGK (DE): 1 WGK No.: 0788

Storage class (VCI): 12

Chemical: polyoxyethylene sorbitan monolaurate (Tween® 20)

CAS No.: 9005-64-5

WGK (DE): 1

Storage class (VCI): 10-11

### 12.2 Persistence and degradability

No data available.

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



## 12.3 Bioaccumulative potential

No data available.

## 12.4 Mobility in soil

No data available.

## 12.5 Results of PBT and vPvB assessment

No data available.

## 12.6 Endocrine disrupting properties

No data available.

## 12.7 Other adverse effects

No data available.

## Section 13. DISPOSAL CONSIDERATIONS

Refer to local regulations.

### 13.1 Waste treatment methods

Dispose of contents/container in accordance with local regulation.

Empty containers should be taken to an approved waste handling site for recycling or disposal.

## Section 14. TRANSPORT INFORMATION

### 14.1 UN number

ADR/RID: -

IMDG: -

IATA: -

### 14.2 UN proper shipping name

ADR/RID: Not dangerous goods

IMDG: Not dangerous goods

IMDG: Not dangerous goods

### 14.3 Transport hazard class(es)

ADR/RID: -

IMDG: -

IATA: -

### 14.4 Packaging group

ADR/RID: -

IMDG: -

IATA: -

### 14.5 Environmental hazards

ADR/RID: no

IMDG Marine pollutant: no

IATA-DGR: no

### 14.6 Special precautions for user

None.

### 14.7 Maritime transport in bulk according to IMO instruments

Not applicable.

#### Further information

Not classified as dangerous in transport regulations.

## Section 15. REGULATIONS INFORMATION

### 15.1 Safety, health, and environmental regulations/legislation specific for the substance or mixture

Look for your country-specific regulations.

### 15.2 Chemical safety assessment

Not necessary.

## Section 16. OTHER INFORMATION

### 16.1 List of R, H and P phrases

Full text of H-Statements referred to under sections 2 and 3.

### 16.2 Training Advice

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



Regular safety training.

## 16.3 Recommended Restriction on Use

Only for professional/research user.

## Section 1. IDENTIFICATION

### 1.1 Product Identifier

Commercial Product Name: **W1 Buffer**

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

Product form: Laboratory chemical mixture

Relevant identified uses: For research use only, not for diagnostic use.

Uses advised against:

-----

### 1.3 Details of the supplier of the safety data sheet

Company: Favorgen Biotech Corp.

Address: Ping-Tung Agricultural Biotechnology Park No. 8-1, Yuanxi 2<sup>nd</sup> Road, Ping Tung 908126, TAIWAN

Phone Number: +886-8-762-1829

### 1.4 Emergency telephone number

+886-8-762-1829

## Section 2. HAZARDS IDENTIFICATION

### 2.1 Classification of the substance or mixture

Hazard identification Hazard classes/categories

Acute oral toxicity	Category 4
Skin corrosion/irritation	Category 2
Serious eye damage/eye irritation	Category 2

GHS hazard pictograms:



GHS07

Signal word: **WARNING**

Hazard statement(s)

H302	Harmful if swallowed
H315	Causes skin irritation
H319	Causes serious eye irritation
H412	Harmful to aquatic life with long lasting effects

Precautionary statement(s)

P264	Wash hands thoroughly after handling
P270	Do not eat, drink or smoke when using this product
P273	Avoid release to the environment
P280	Wear protective gloves/protective clothing/ eye protection/ face protection
P301 + P310	IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician
P302 + P352	IF ON SKIN: Wash with plenty of soap and water
P305 + P351 + P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
P330	Rinse mouth
P332 + P313	If skin irritation occurs: Get medical advice/attention

### 2.2 Label elements

According to section 1.5.2 of Annex I to CLP

GHS symbol



GHS07

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



## WARNING

### Hazard contents

Guanidine hydrochloride, 20~50%, CAS-No. 50-01-1

### 2.3 Other hazards

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

## Section 3. COMPOSITION/INFORMATION ON INGREDIENTS

### 3.1 Substances or 3.2 Mixtures

Chemical:	<b>Guanidine hydrochloride</b>	CAS No.: 50-01-1
Concentration:	20 - <50 %	
Classification:	H302, Acute oral Tox. 4, H315, Skin Irrit. 2, H319, Eye Irrit. 2	

## Section 4. FIRST AID MEASURES

### 4.1 Description of first-aid measures

#### If inhaled

After inhalation: fresh air. If breathing stops: mouth-to-mouth breathing or artificial respiration. Oxygen if necessary. Immediately call in physician.

#### In case of skin contact

After skin contact: Take off immediately all contaminated clothing. Rinse skin with water/shower.

#### In case of eye contact

After eye contact: rinse out with plenty of water. Call in ophthalmologist. Remove contact lenses.

#### If swallowed

After swallowing: immediately make victim drink water (two glasses at most). Consult a physician.

### 4.2 Most important symptoms and effects, both acute and delayed

H302 - Harmful if swallowed  
H315 - Causes skin irritation  
H319 - Causes serious eye irritation  
H412 - Harmful to aquatic life

### 4.3 Indication of any immediate medical attention and special treatment needed

No data available

## Section 5. FIRE-FIGHTING MEASURES

### 5.1 Extinguishing media

All extinguishers can be used.

### 5.2 Special hazards arising from the substance or mixture

None known.

### 5.3 Advice for firefighters

No specific fire or explosion hazard.

### 5.4 Additional Information

---

## Section 6. ACCIDENTAL RELEASE MEASURE

### 6.1 Personal precautions, protective equipment, and emergency procedures

Do not breathe vapours/Ensure adequate ventilation/Avoid contact with skin, eyes or clothing/Use personal protection equipment/Regular staff training is necessary.

### 6.2 Environmental precautions

Do not allow mixture to enter ground water system.

### 6.3 Methods and material for containment and cleaning up

Soak up with inert absorbent material.

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



## 6.4 Reference to other sections

---

## Section 7. HANDLING AND STORAGE

### 7.1 Precautions for safe handling

Always wear recommended Personal Protective Equipment.

### 7.2 Conditions for safe storage, including any incompatibilities

Keep in a dry, cool, and well-ventilated place.

### 7.3 Specific end use(s)

For research use only.

## Section 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

### 8.1 Control parameters: None

Chemical:

Guanidine hydrochloride

CAS No.: 50-01-1

NIOSH: not listed

OSHA: not listed

### 8.2 Exposure controls

Good ventilation and exhaust system in the workplace, and the floor has drainage and washing facilities, which can resist the erosion of chemicals. Must maintain the highest level of cleanliness in the workplace.

#### Personal protective equipment

**Eye/face protection:** Safety glasses. Used to prevent splash hazards. Government standards such as NIOSH (US) or EN 166(EU).

**Skin protection:** Recommended to avoid contamination with these hazards.

**Hand protection:** Protective gloves, use for short times chemical resistant latex gloves with code EN 374-3 level 1.

**Respiratory protection:** Not necessary, if the workplace is well-ventilated.

**Personal Hygiene:** Wash hands before breaks and at the end of work. Handle in accordance with good industrial hygiene and safety practice.

## Section 9. PHYSICAL AND CHEMICAL PROPERTIES

### 9.1 Information on basic physical and chemical properties

Appearance: liquid/Color: clear/Odor: odorless/pH: 6.0~7.0

### 9.2 Other information

---

## Section 10. STABILITY AND REACTIVITY

### 10.1 Reactivity

No data available.

### 10.2 Chemical stability

No data available.

### 10.3 Possibility of hazardous reactions

**Note:** Can form very reactive substances with oxidizing agents.

### 10.4 Conditions to avoid

Not necessary.

### 10.5 Incompatible materials

**Note:** Avoid contact with strong acids or alkalines.

### 10.6 Hazardous decomposition products

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



Decompositions are not observed during the expiration period under recommended conditions.

## Section 11. TOXICOLOGICAL INFORMATION

### 11.1 Information on toxicological effects

#### Acute toxicity

Not Classified

Chemical Name	Oral LD50 (Rat)	Dermal LD50 (Rabbit)	Inhalation LC50 (Rat)
Guanidine hydrochloride CAS No.: 50-01-1	475-907 mg/kg	2000 mg/kg	[4h] 3181-7655 µg/m <sup>3</sup>

#### Skin corrosion/irritation

Skin irritation/Irritating to eyes

#### Serious eye damage/irritation

Not Classified

#### Respiratory or skin sensitization

Not Classified

#### Germ cell mutagenicity

Not Classified

#### Carcinogenicity

Not Classified

#### Reproductive toxicity

Not Classified

#### Specific target organ toxicity (STOT) – single exposure

Not Classified

#### Specific target organ toxicity (STOT) – repeated exposure

Not Classified

#### Aspiration hazard

Not Classified

### 11.2 Information on other hazards

No data available.

## Section 12. ECOLOGICAL INFORMATION

### 12.1 Toxicity

Following information is valid for pure substances.

Guanidine hydrochloride

CAS No.: 50-01-1

LC50leuciscus idus/96h: 1759 mg/L

WGK (DE): 1 WGK No.: 0788

Storage class (VCI): 12

### 12.2 Persistence and degradability

No data available.

### 12.3 Bioaccumulative potential

No data available.

### 12.4 Mobility in soil

No data available.

### 12.5 Results of PBT and vPvB assessment

No data available.

### 12.6 Endocrine disrupting properties

No data available.

### 12.7 Other adverse effects

No data available.

## Section 13. DISPOSAL CONSIDERATIONS

Refer to local regulations.

### 13.1 Waste treatment methods

Dispose of contents/container in accordance with local regulation.

Empty containers should be taken to an approved waste handling site for recycling or disposal.

## Section 14. TRANSPORT INFORMATION

### 14.1 UN number

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



ADR/RID: -

IMDG: -

IATA: -

## 14.2 UN proper shipping name

ADR/RID: Not dangerous goods

IMDG: Not dangerous goods

IMDG: Not dangerous goods

## 14.3 Transport hazard class(es)

ADR/RID: -

IMDG: -

IATA: -

## 14.4 Packaging group

ADR/RID: -

IMDG: -

IATA: -

## 14.5 Environmental hazards

ADR/RID: no

IMDG Marine pollutant: no

IATA-DGR: no

## 14.6 Special precautions for user

None

## 14.7 Maritime transport in bulk according to IMO instruments

Not Applicable.

### Further information

Not classified as dangerous in transport regulations.

## Section 15. REGULATIONS INFORMATION

### 15.1 Safety, health, and environmental regulations/legislation specific for the substance or mixture

Look for your country-specific regulations.

### 15.2 Chemical safety assessment

Not necessary.

## Section 16. OTHER INFORMATION

### 16.1 List of R, H and P phrases

Full text of H-Statements referred to under sections 2 and 3.

### 16.2 Training Advice

Regular safety training.

### 16.3 Recommended Restriction on Use

Only for professional/research user.

## Section 1. IDENTIFICATION

### 1.1 Product Identifier

Commercial Product Name: Wash Buffer

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

Product form: Laboratory chemical mixture

Relevant identified uses: For research use only, not for diagnostic use.

Uses advised against:

-----

### 1.3 Details of the supplier of the safety data sheet

Company: FAVORGEN Biotech Corp.

Address: Ping-Tung Agricultural Biotechnology Park No. 8-1, Yuanxi 2<sup>nd</sup> Road, Ping Tung 908126, TAIWAN

Phone Number: +886-8-762-1829

### 1.4 Emergency telephone number

+886-8-762-1829

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



## Section 2. HAZARDS IDENTIFICATION

### 2.1 Classification of the substance or mixture

No hazardous substance as specified in Regulation (EC) No 1272/2008.

### 2.2 Label elements

Does not need labelling as hazardous.

### 2.3 Other hazards

Not applicable.

## Section 3. COMPOSITION/INFORMATION ON INGREDIENTS

### 3.1 Substances or 3.2 Mixtures

Chemical:	chemicals/mixture < 1%,	no declaration necessary	CAS No.: -----
Concentration:	0.1 - <1%		

## Section 4. FIRST AID MEASURES

### 4.1 Description of first-aid measures

#### If inhaled

Supply fresh air.

#### In case of skin contact

Remove contaminated clothing. Rinse skin with water/shower.

#### In case of eye contact

Rinse out with plenty of water.

#### If swallowed

Consult a doctor if feeling unwell.

### 4.2 Most important symptoms and effects, both acute and delayed

No information available.

### 4.3 Indication of any immediate medical attention and special treatment needed

No information available.

## Section 5. FIRE-FIGHTING MEASURES

### 5.1 Extinguishing media

All extinguishers can be used.

### 5.2 Special hazards arising from the substance or mixture

None.

### 5.3 Advice for firefighters

No explosion hazard.

### 5.4 Additional Information

---

## Section 6. ACCIDENTAL RELEASE MEASURE

### 6.1 Personal precautions, protective equipment, and emergency procedures

Ensure adequate ventilation.

### 6.2 Environmental precautions

Do not allow mixture to enter ground water system.

### 6.3 Methods and material for containment and cleaning up

Soak up with inert absorbent material.

### 6.4 Reference to other sections

---



# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



## Section 7. HANDLING AND STORAGE

### 7.1 Precautions for safe handling

Always wear recommended Personal Protective Equipment.

### 7.2 Conditions for safe storage, including any incompatibilities

Keep container tightly closed in a dry and well-ventilated place.

### 7.3 Specific end use(s)

For research use only.

## Section 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

### 8.1 Control parameters:

Contains no substances with occupational exposure limit.

### 8.2 Exposure controls

#### Personal protective equipment

**Personal Hygiene:** Wash hands before breaks and at the end of work. Handle in accordance with good industrial hygiene and safety practice.

**Eye/face protection:** Tight sealing safety glasses.

**Hand protection:** Protective gloves.

**Respiratory protection:** Not necessary, if the workplace is well-ventilated.

## Section 9. PHYSICAL AND CHEMICAL PROPERTIES

### 9.1 Information on basic physical and chemical properties

Appearance: liquid/Color: Clear/Odor: odorless/pH: 7.5~8.5

### 9.2 Other information

---

## Section 10. STABILITY AND REACTIVITY

### 10.1 Reactivity

None.

### 10.2 Chemical stability

No data available.

### 10.3 Possibility of hazardous reactions

None.

### 10.4 Conditions to avoid

Not necessary.

### 10.5 Incompatible materials

None.

### 10.6 Hazardous decomposition products

Decompositions are not observed during the expiration period under recommended conditions.

## Section 11. TOXICOLOGICAL INFORMATION

### 11.1 Information on toxicological effects

Acute toxicity	Not Classified
Skin corrosion/irritation	Not Classified
Serious eye damage/irritation	Not Classified
Respiratory or skin sensitization	Not Classified
Germ cell mutagenicity	Not Classified
Carcinogenicity	Not Classified

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



<b>Reproductive toxicity</b>	Not Classified
<b>Specific target organ toxicity (STOT) – single exposure</b>	Not Classified
<b>Specific target organ toxicity (STOT) – repeated exposure</b>	Not Classified
<b>Aspiration hazard</b>	Not Classified

## 11.2 Information on other hazards

No data available.

## Section 12. ECOLOGICAL INFORMATION

### 12.1 Toxicity

Contains no substances known to be hazardous to the ecosystem.

### 12.2 Persistence and degradability

No data available.

### 12.3 Bioaccumulative potential

No data available.

### 12.4 Mobility in soil

No data available.

### 12.5 Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be a PBT or a vPvB.

### 12.6 Endocrine disrupting properties

No data available.

### 12.7 Other adverse effects

No data available.

## Section 13. DISPOSAL CONSIDERATIONS

Refer to local regulations.

### 13.1 Waste treatment methods

Dispose of contents/container in accordance with local regulation.

Empty containers should be taken to an approved waste handling site for recycling or disposal.

## Section 14. TRANSPORT INFORMATION

### 14.1 UN number

ADR/RID: -

IMDG: -

IATA: -

### 14.2 UN proper shipping name

ADR/RID: Not dangerous goods

IMDG: Not dangerous goods

IMDG: Not dangerous goods

### 14.3 Transport hazard class(es)

ADR/RID: -

IMDG: -

IATA: -

### 14.4 Packaging group

ADR/RID: -

IMDG: -

IATA: -

### 14.5 Environmental hazards

ADR/RID: no

IMDG Marine pollutant: no

IATA-DGR: no

### 14.6 Special precautions for user

None.

### 14.7 Maritime transport in bulk according to IMO instruments

Not Applicable.

### Further information

Not classified as dangerous in transport regulations.

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



## Section 15. REGULATIONS INFORMATION

### 15.1 Safety, health, and environmental regulations/legislation specific for the substance or mixture

Look for your country-specific regulations.

### 15.2 Chemical safety assessment

Not necessary.

## Section 16. REGULATIONS INFORMATION

### 16.1 List of R, H and P phrases

Full text of H-Statements referred to under sections 2 and 3.

### 16.2 Training Advice

Regular safety training.

### 16.3 Recommended Restriction on Use

Only for professional/research user.

## Section 1. IDENTIFICATION

### 1.1 Product Identifier

Commercial Product Name: Elution Buffer

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

Product form: Laboratory chemical mixture

Relevant identified uses: For research use only, not for diagnostic use.

Uses advised against:

-----

### 1.3 Details of the supplier of the safety data sheet

Company: FAVORGEN Biotech Corp.

Address: Ping-Tung Agricultural Biotechnology Park No. 8-1, Yuanxi 2<sup>nd</sup> Road, Ping Tung 908126, TAIWAN

Phone Number: +886-8-762-1829

### 1.4 Emergency telephone number

+886-8-762-1829

## Section 2. HAZARDS IDENTIFICATION

### 2.1 Classification of the substance or mixture

No hazardous substance as specified in Regulation (EC) No 1272/2008.

### 2.2 Label elements

Does not need labelling as hazardous.

### 2.3 Other hazards

Not applicable.

## Section 3. COMPOSITION/INFORMATION ON INGREDIENTS

### 3.1 Substances or 3.2 Mixtures

Chemical:	chemicals/mixture < 1%,	no declaration necessary	CAS No.: -----
Concentration:	0.1 - <1%		

## Section 4. FIRST AID MEASURES

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



## 4.1 Description of first-aid measures

### If inhaled

Supply fresh air.

### In case of skin contact

Remove contaminated clothing. Rinse skin with water/shower.

### In case of eye contact

Rinse out with plenty of water.

### If swallowed

Consult a doctor if feeling unwell.

## 4.2 Most important symptoms and effects, both acute and delayed

No information available.

## 4.3 Indication of any immediate medical attention and special treatment needed

No information available.

## Section 5. FIRE-FIGHTING MEASURES

### 5.1 Extinguishing media

All extinguishers can be used.

### 5.2 Special hazards arising from the substance or mixture

None.

### 5.3 Advice for firefighters

No explosion hazard.

### 5.4 Additional Information

---

## Section 6. ACCIDENTAL RELEASE MEASURE

### 6.1 Personal precautions, protective equipment, and emergency procedures

Ensure adequate ventilation.

### 6.2 Environmental precautions

Do not allow mixture to enter ground water system.

### 6.3 Methods and material for containment and cleaning up

Soak up with inert absorbent material.

### 6.4 Reference to other sections

---

## Section 7. HANDLING AND STORAGE

### 7.1 Precautions for safe handling

Always wear recommended Personal Protective Equipment.

### 7.2 Conditions for safe storage, including any incompatibilities

Keep container tightly closed in a dry and well-ventilated place.

### 7.3 Specific end use(s)

For research use only.

## Section 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

### 8.1 Control parameters:

Contains no substances with occupational exposure limit.

### 8.2 Exposure controls

#### Personal protective equipment

**Personal Hygiene:** Wash hands before breaks and at the end of work. Handle in accordance with good industrial hygiene and safety practice.

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



**Eye/face protection:** Tight sealing safety glasses.

**Hand protection:** Protective gloves.

**Respiratory protection:** Not necessary, if the workplace is well-ventilated.

## Section 9. PHYSICAL AND CHEMICAL PROPERTIES

### 9.1 Information on basic physical and chemical properties

Appearance: liquid/Color: Clear/Odor: odorless/pH: 8.0~9.0

### 9.2 Other information

---

## Section 10. STABILITY AND REACTIVITY

### 10.1 Reactivity

None.

### 10.2 Chemical stability

No data available.

### 10.3 Possibility of hazardous reactions

None.

### 10.4 Conditions to avoid

Not necessary.

### 10.5 Incompatible materials

None.

### 10.6 Hazardous decomposition products

Decompositions are not observed during the expiration period under recommended conditions.

## Section 11. TOXICOLOGICAL INFORMATION

### 11.1 Information on toxicological effects

Acute toxicity	Not Classified
Skin corrosion/irritation	Not Classified
Serious eye damage/irritation	Not Classified
Respiratory or skin sensitization	Not Classified
Germ cell mutagenicity	Not Classified
Carcinogenicity	Not Classified
Reproductive toxicity	Not Classified
Specific target organ toxicity (STOT) – single exposure	Not Classified
Specific target organ toxicity (STOT) – repeated exposure	Not Classified
Aspiration hazard	Not Classified

### 11.2 Information on other hazards

No data available.

## Section 12. ECOLOGICAL INFORMATION

### 12.1 Toxicity

Contains no substances known to be hazardous to the ecosystem.

### 12.2 Persistence and degradability

No data available.

### 12.3 Bioaccumulative potential

No data available.

### 12.4 Mobility in soil

No data available.

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



## 12.5 Results of PBT and vPvB assessment

This mixture does not contain any substances that are assessed to be a PBT or a vPvB.

## 12.6 Endocrine disrupting properties

No data available.

## 12.7 Other adverse effects

No data available.

## Section 13. DISPOSAL CONSIDERATIONS

Refer to local regulations.

### 13.1 Waste treatment methods

Dispose of contents/container in accordance with local regulation.

Empty containers should be taken to an approved waste handling site for recycling or disposal.

## Section 14. TRANSPORT INFORMATION

### 14.1 UN number

ADR/RID: -

IMDG: -

IATA: -

### 14.2 UN proper shipping name

ADR/RID: Not dangerous goods

IMDG: Not dangerous goods

IMDG: Not dangerous goods

### 14.3 Transport hazard class(es)

ADR/RID: -

IMDG: -

IATA: -

### 14.4 Packaging group

ADR/RID: -

IMDG: -

IATA: -

### 14.5 Environmental hazards

ADR/RID: no

IMDG Marine pollutant: no

IATA-DGR: no

### 14.6 Special precautions for user

None

### 14.7 Maritime transport in bulk according to IMO instruments

Not Applicable.

#### Further information

Not classified as dangerous in transport regulations.

## Section 15. REGULATIONS INFORMATION

### 15.1 Safety, health, and environmental regulations/legislation specific for the substance or mixture

Look for your country-specific regulations.

### 15.2 Chemical safety assessment

Not necessary.

## Section 16. REGULATIONS INFORMATION

### 16.1 List of R, H and P phrases

Full text of H-Statements referred to under sections 2 and 3.

### 16.2 Training Advice

Regular safety training.

### 16.3 Recommended Restriction on Use

Only for professional/research user.

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



## Section 1. IDENTIFICATION

### 1.1 Product Identifier

Commercial Product Name: **Proteinase K (Liquid)**

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

Product form: Laboratory chemical mixture

Relevant identified uses: For research use only, not for diagnostic use.

Uses advised against:

-----

### 1.3 Details of the supplier of the safety data sheet

Company: FAVORGEN Biotech Corp.

Address: Ping-Tung Agricultural Biotechnology Park No. 8-1, Yuanxi 2<sup>nd</sup> Road, Ping Tung 908126, TAIWAN

Phone Number: +886-8-762-1829

### 1.4 Emergency telephone number

+886-8-762-1829

## Section 2. HAZARDS IDENTIFICATION

### 2.1 Classification of the substance or mixture

Hazard identification Hazard classes/categories

Resp. Sens.	Category 1
Sensitization/skin	Category 1

GHS hazard pictograms:



GHS08

Signal word: **DANGER**

Hazard statement(s)

H317

May cause an allergic skin reaction.

H334

May cause allergy or asthma symptoms or breathing difficulties if inhaled.

Precautionary statement(s)

P101

If medical advice is needed, have product container or label at hand.

P102

Keep out of reach of children.

P103

Read label before use.

P261

Avoid breathing dust/fume/gas/mist/vapours/spray.

P280

Wear protective gloves/protective clothing/eye protection/face protection/hearing protection/...

P284

[In case of inadequate ventilation] wear respiratory protection.

P304+P340

IF INHALED: Remove person to fresh air and keep comfortable for breathing.

P342+P311

If experiencing respiratory symptoms: Call a POISON CENTER/doctor.

P501

Dispose of contents/container in accordance with local/regional/national/international regulations.

Supplemental Hazard Statements

None

### 2.2 Label elements

According to section 1.5.2 of Annex I to CLP

GHS symbol



GHS08

DANGER

### 2.3 Other hazards

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



## Section 3. COMPOSITION/INFORMATION ON INGREDIENTS

### 3.1 Substances or 3.2 Mixtures

Chemical: **Proteinase K**  
Concentration: 1- ≤2.5%  
Classification: H334, Resp. Sens. 1,  
H315, Skin Irrit. 2,  
H317, Skin Sens.1  
H319, Eye Irrit. 2  
H335, STOT SE 3

CAS No.: 39450-01-6

## Section 4. FIRST AID MEASURES

### 4.1 Description of first-aid measures

#### If inhaled

After inhalation: Supply fresh air; consult doctor in case of complaints.

#### In case of skin contact

After skin contact: Immediately wash with water and soap and rinse thoroughly. Seek medical treatment.

#### In case of eye contact

After eye contact: Rinse out with plenty of water. Call in ophthalmologist. Remove contact lenses. Call a doctor immediately.

#### If swallowed

After swallowing: Call a doctor immediately.

### 4.2 Most important symptoms and effects, both acute and delayed

The most important known symptoms and effects are described in the labelling (see section 2.2) and/or in section 11.

### 4.3 Indication of any immediate medical attention and special treatment needed

No data available

## Section 5. FIRE-FIGHTING MEASURES

### 5.1 Extinguishing media

All extinguishers can be used.

### 5.2 Special hazards arising from the substance or mixture

None known.

### 5.3 Advice for firefighters

No specific fire or explosion hazard.

### 5.4 Additional Information

---

## Section 6. ACCIDENTAL RELEASE MEASURE

### 6.1 Personal precautions, protective equipment, and emergency procedures

Do not breathe vapours/Ensure adequate ventilation/Avoid contact with skin, eyes or clothing/Use personal protection equipment/Regular staff training is necessary.

### 6.2 Environmental precautions

Do not allow mixture to enter ground water system.

### 6.3 Methods and material for containment and cleaning up

Ensure adequate ventilation.  
Soak up the dissolved solution with inert absorbent material.

### 6.4 Reference to other sections

---



# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



## Section 7. HANDLING AND STORAGE

### 7.1 Precautions for safe handling

Always wear recommended Personal Protective Equipment.

### 7.2 Conditions for safe storage, including any incompatibilities

Storage conditions: Tightly closed. Dry. Keep locked up or in an area accessible only to qualified or authorized persons.

### 7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated.

## Section 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

### 8.1 Control parameters: None

Ingredients with workplace control parameters.

### 8.2 Exposure controls

Good ventilation and exhaust system in the workplace, and the floor has drainage and washing facilities, which can resist the erosion of chemicals. Must maintain the highest level of cleanliness in the workplace.

#### Personal protective equipment

**Eye/face protection:** Safety glasses. Used to prevent splash hazards. Government standards such as NIOSH (US) or EN 166(EU).

**Skin protection:** Recommended to avoid contamination with these hazards.

**Hand protection:** Protective gloves, use for short times chemical resistant latex gloves with code EN 374-3 level 1.

**Respiratory protection:** Not necessary, if the workplace is well-ventilated.

**Personal Hygiene:** Wash hands before breaks and at the end of work. Handle in accordance with good industrial hygiene and safety practice.

## Section 9. PHYSICAL AND CHEMICAL PROPERTIES

### 9.1 Information on basic physical and chemical properties

Appearance: liquid/Color: Colorless/Odor: Odorless/pH: 7.5

### 9.2 Other information

---

## Section 10. STABILITY AND REACTIVITY

### 10.1 Reactivity

No further relevant information available.

### 10.2 Chemical stability

The product is chemically stable under standard ambient conditions (room temperature). Stable under recommended storage conditions.

### 10.3 Possibility of hazardous reactions

No data available.

### 10.4 Conditions to avoid

No information available.

### 10.5 Incompatible materials

Strong oxidizing agents

### 10.6 Hazardous decomposition products

In the event of fire: see section 5.

## Section 11. TOXICOLOGICAL INFORMATION

# Material Safety Data Sheet

According to Regulation (EC) No 1907/2006



## 11.1 Information on toxicological effects

Acute toxicity	Not Classified
Skin corrosion/irritation	Causes skin irritation.
Serious eye damage/irritation	Causes serious eye irritation.
Respiratory or skin sensitization	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
Germ cell mutagenicity	Not Classified
Carcinogenicity	Based on available data, the classification criteria are not met.
Reproductive toxicity	Not Classified
Specific target organ toxicity (STOT) – single exposure	Not Classified
Specific target organ toxicity (STOT) – repeated exposure	Not Classified
Aspiration hazard	Not Classified

## 11.2 Information on other hazards

No data available.

## Section 12. ECOLOGICAL INFORMATION

### 12.1 Toxicity

No data available.

### 12.2 Persistence and degradability

No data available.

### 12.3 Bioaccumulative potential

No data available.

### 12.4 Mobility in soil

No data available.

### 12.5 Results of PBT and vPvB assessment

No data available.

### 12.6 Endocrine disrupting properties

No data available.

### 12.7 Other adverse effects

No additional data available.

## Section 13. DISPOSAL CONSIDERATIONS

Refer to local regulations.

### 13.1 Waste treatment methods

Dispose of contents/container in accordance with local regulation.

Empty containers should be taken to an approved waste handling site for recycling or disposal.

## Section 14. TRANSPORT INFORMATION

### 14.1 UN number

ADR/RID: -

IMDG: -

IATA: -

### 14.2 UN proper shipping name

ADR/RID: Not dangerous goods

IMDG: Not dangerous goods

IMDG: Not dangerous goods

### 14.3 Transport hazard class(es)

ADR/RID: -

IMDG: -

IATA: -

### 14.4 Packaging group

ADR/RID: -

IMDG: -

IATA: -

### 14.5 Environmental hazards

ADR/RID: no

IMDG Marine pollutant: no

IATA-DGR: no

**14.6 Special precautions for user**

None

**14.7 Maritime transport in bulk according to IMO instruments**

Not Applicable.

**Further information**

Not classified as dangerous in transport regulations.

**Section 15. REGULATIONS INFORMATION**

**15.1 Safety, health, and environmental regulations/legislation specific for the substance or mixture**

This material safety data sheet complies with the requirements of Regulation (EC) No.1907/2006.

**15.2 Chemical safety assessment**

For this product, a chemical safety assessment was not carried out

**Section 16. OTHER INFORMATION**

**16.1 List of R, H and P phrases**

Full text of H-Statements referred to under sections 2 and 3.

**16.2 Training Advice**

Regular safety training.

**16.3 Recommended Restriction on Use**

Only for professional/research user.