

Macrovolume Tips

5 mL, Macrovolume

Fits Axygen® Axypet pipettors, Corning® Lambda™ Plus pipettors, and other Gilson-style pipettors.

- ▶ RNase-/DNase-free
- ▶ Nonpyrogenic



Graduation marks at 0.5 mL and 1 mL increments

Cat. No.	Volume (μL)	Pack Type	Packaging
T-5000-C	1 - 5	Bulk	250/pk, 10 pk/cs

10 mL, Macrovolume

Fits most research grade pipettors.



Cat. No.	Volume (μL)	Pack Type	Packaging
T-10ML-C	1 - 10	Bulk	200/pk, 10 pk/cs