



PCR Tubes

Axygen polypropylene, thin wall PCR tubes with and without caps are designed to fit standard 96-well blocks and are available in an assortment of colors. The consistent and ultra-thin wall allows precise thermal transfer. An even insertion cap ensures a tight fit to reduce evaporation. Available with Axxygen Maxymum Recovery® surface.

Standard Surface Cat. No.	Description	Volume (μL)	Color	Packaging
PCR-02-C	Thin wall, flat cap	0.2	Clear	1,000/pk, 10 pk/cs
PCR-02-A	Thin wall, flat cap	0.2	Assorted	1,000/pk, 10 pk/cs
PCR-02-NC	Thin wall, without cap	0.2	Clear	1,000/pk, 10 pk/cs
PCR-05-C	Thin wall, flat cap	0.5	Clear	1,000/pk, 10 pk/cs
PCR-05-A	Thin wall, flat cap	0.5	Assorted	1,000/pk, 10 pk/cs