

# Axygen PCR Consumables

to assist with your Molecular Biology requirements

**AXYGEN**

A Corning Brand

## Axygen PCR Microplates and Sealing Films



**SPECIAL  
DISCOUNTED  
PRICING**

### Features and Benefits

- ▶ 96 well and 384 well
- ▶ Available in clear and white for standard and real-time PCR
- ▶ Available in a wide variety of sizes and configurations
- ▶ Certified RNase-/DNase-free

### Ordering Information

#### 96 & 384 Well PCR Microplates

Cat. No.	Description	Color	Volume	Qty/Cs	Price
AX-PCR-96M2-HS-C	96 well PCR half skirt microplate with a single notch	Clear	200 µL	50	<b>\$226</b>
AX-PCR-96-AB-C	96 well PCR microplate with half skirt, compatible with ABI	Clear	200 µL	50	<b>\$341</b>
AX-PCR-96-FLT-C	96 well PCR flat top microplate with no skirt	Clear	200 µL	100	<b>\$428</b>
AX-PCR-96-LC480-W-NF	96 well PCR microplate with half skirt, compatible with Roche 480	White	100 µL	50	<b>\$226</b>
AX-PCR-384-LC480-W-NF	384 well PCR microplate, compatible with Roche Lightcycler 480	White	20 µL	50	<b>\$296</b>
AX-PCR-384M2-C	384 well PCR microplate with single notch	Clear	55 µL	50	<b>\$231</b>

#### Axygen® Adhesive Sealing Films

Cat. No.	Description	Qty/Pk	Application	
AX-PCR-TS-1	CyclerSeal sealing film, polyester	100	Non-optical PCR, storage	<b>\$62</b>
AX-PCR-SP-1	AxySeal sealing film, polyester	100	Non-optical PCR, ELISA	<b>\$66</b>
AX-UC-500-1	UltraClear, pressure sensitive sealing film, polyester	100	Real-time PCR	<b>\$217</b>

\*Sealing films are not autoclavable.

**FB**  
AUSTRALIA

#### Australian Distributors:

**FB AUSTRALIA**

Free Call: 1800 066 077

Email: [info@fisherbiotec.com](mailto:info@fisherbiotec.com)

Web: [fisherbiotec.com.au](http://fisherbiotec.com.au)