

StrepMAB-Immo

Cat. No.: 2-1517-001

Description	StrepMAB-Immo is a murine, high-affinity monoclonal antibody especially suited for capturing and detection of Strep-tag®II and Twin-Strep-tag® fusion proteins.
Subclass	IgG1
Specificity	The nearly irreversible binding is validated for fusion proteins carrying a C- or N-terminal Strep tag®II with an N-terminal SerAla linker (protein-SA-WSHQPFEK or SA-WSHQPFEK-protein). Strep tag®II proteins without the N-terminal SerAla are bound with reduced affinity in the cases tested so far. The Twin-Strep-tag® always contains the SerAla linker (WSHQPFEK-GGGG-GGGG-GGSA-WSHQPFEK).
Form	Lyophilized from PBS solution (2.5 mM KCl, 8 mM Na ₂ HPO ₄ , 200 mM NaCl, 1.5 mM KH ₂ PO ₄ , pH 7.4)
Amount	100 µg
Reconstitution	Reconstitute antibody in 200 µl water to get a concentration of 0.5 mg/ml in PBS. Aliquot and store at -20 °C.
Stability	24 months after shipping
Storage	Recommended: -25 °C to -15 °C (both lyophilized and after reconstitution)
Shipping	Room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.
Application	StrepMAB-Immo is suited for capturing of Strep-tag®II and Twin-Strep-tag® fusion proteins on solid phases (e.g., microplates or biochips) or for detection of those proteins in applications like ELISA, FACS, immunohistochemistry, and -cytochemistry. For detection in western blot, the use of StrepMAB-Classic or Strep-Tactin® conjugates is recommended.
Manufacturer	IBA Lifesciences

For research use only

Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale IBA Lifesciences grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at info@iba-lifesciences.com or at IBA Lifesciences GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.