

FB

Sleek, Reliable, and Intuitive

Esco CelCulture[®]

CO2 Incubator

To Grow & Maintain Cell Cultures

170L Stainless-Steel Chamber

Special Price of \$6,400 (ex GST)

KEY FEATURES

- Direct heat and air jacket for fast and uniform heating
- Heat-resistant, maintenance-free, in-chamber IR CO2 sensor
- ULPA continuous, active air filtration system
- Inlet filter for gas inputs, and air-exchange
- ISO Class 5 work zone with easy clean rounded corners with cable access port
- Exterior surfaces with Isocide™ antimicrobial powder coating
- Validated 90° C moist heat decontamination cycle
- On board data-logging, and with 2-year warranty



35% DISCOUNT (when purchased with incubator)



COA-2001-F (170L) Roller Base

COA-2005-F 2-Stage Gas Regulator for CO2

Catalog No.	Description	Volume	QTY	Price (ex GST)
2170002	Esco CCL-170B-8, 170L CO2 Incubator, IR sensor, with stainless steel chamber, ULPA filter, moist heat decon	170L	1	\$6,400
5170478	COA-2001-F (170L) Roller Base for CCL-170		1	\$289
5170481	COA-2005-F 2-Stage Gas Regulator for CO2		1	\$299



FB

CERTIFIED
ISO 9001



Australian Distributors

Freecall: 1800 066 077 | Email: info@fisherbiotec.com
fisherbiotec.com

Terms and Conditions

All prices are quoted in AUD and exclude GST.
Offers exclude servicing, spare parts and freight.
Valid until 31st Dec 2025