

FB

Sorenson Bioscience™

Barrier Tips

UNIVERSAL FIT

Universal Fit, Sterile Filter Tips

10 - 1250 UL

- Fits standard universal pipettors
- Features an inert barrier that prevents pipettor contamination
- Barrier will not trap liquid or contaminate samples in the event of over-pipetting
- Available in low binding formats



Catalog No.	Description	Qty
SOR-30340-1	Sorenson 10ul Sharp Tip, Barrier, Sterile	960
SOR-30360T-1	Sorenson 10ul LongReach Tip, Barrier, Sterile, Low Binding	960
SOR-30370T-1	Sorenson 20ul Standard Tip, Barrier, Sterile, Low Binding	960
SOR-36060T-1	Sorenson 100ul MultiGuard Barrier Tip, Sterile, Low Binding	960
SOR-30380T-1	Sorenson 100ul V Tip, Barrier, Sterile, Low Binding	960
SOR-30550T-1	Sorenson 200ul MultiGuard Barrier Tip, Sterile, Low Binding	960
SOR-30390T-1	Sorenson 200ul V Tip, Barrier, Sterile, Low Binding	960
SOR-14200T-1	Sorenson 1000µl Barrier Tip, Sterile, Low Binding	1000
SOR-34000-1	Sorenson 1000ul XL Barrier Tip, Sterile	768

