


# Calibration kit, d12mm

## DESCRIPTION

McFarland Turbidity Standards, Ø12mm

## SPECIFICATIONS

Diameter	12 mm
----------	-------



## CAT. NUMBER

21255	ø12mm
-------	-------

## ACCESSORIES



### DEN-1

BS-050102-AAF

Densitometer (suspension  
turbidity detector)

Densitometer is designed for measurement of cell suspension's turbidity in the range of 0.0–6.0 McFarland units ( $0 - 180 \times 10^7$  cells/ml).

Densitometer provides the opportunity to measure solution ...

[read more](#)



### DEN-1B

BS-050104-AAF

Densitometer (suspension  
turbidity detector)

Densitometer is designed for measurement of cell suspension's turbidity in the range of 0.0–6.0 McFarland units ( $0 - 180 \times 10^7$  cells/ml).

Densitometers provide the opportunity to measure solution turbidity in ...

[read more](#)



### DEN-600

BS-050109-AAA

Photometer

DEN-600 is a compact, portable, rechargeable battery powered photometer. DEN-600 comprises of 600 nm wavelength optical system, which enables to apply - 1) OD600 method that estimated total number of ...

[read more](#)