

# Glass test tubes 16mm

## DESCRIPTION

Glass Test Tubes 16x100mm, high borosilicate, PP Cap with silicone pad.  
Packing - 100 pcs/box

## SPECIFICATIONS

Height	100 mm
Diameter	16 mm



## CAT. NUMBER

**BS-050102-MK** Glass sample tubes with lids 16x100, 100 pcs.

## ACCESSORIES



### DEN-1

BS-050102-AAF

Densitometer (suspension  
turbidity detector)

Densitometer is designed for measurement of cell suspension's turbidity in the range of 0.0–6.0 McFarland units ( $0 - 180 \times 10^7$  cells/ml).

Densitometer provides the opportunity to measure solution ...

[read more](#)



### DEN-1B

BS-050104-AAF

Densitometer (suspension  
turbidity detector)

Densitometer is designed for measurement of cell suspension's turbidity in the range of 0.0–6.0 McFarland units ( $0 - 180 \times 10^7$  cells/ml).

Densitometers provide the opportunity to measure solution turbidity in ...

[read more](#)



### DEN-600

BS-050109-AAA

Photometer

**DEN-600** is a compact, portable, rechargeable battery powered photometer. DEN-600 comprises of 600 nm wavelength optical system, which enables to apply - 1) OD600 method that estimated total number of ...

[read more](#)