

# PDM, dimpled mat

## DESCRIPTION

Dimpled mat PDM prevents different size tubes from rolling around the platform.

## SPECIFICATIONS

Dimension

215 × 215 mm



## CAT. NUMBER

PDM

PDM, dimpled mat

## ACCESSORIES



**MR-1**  
 BS-010152-AAG  
 Mini-Rocker Shaker

Mini Rocker-Shaker **MR-1** provides regulated gentle rocking motion of the platform and is ideal for mini gel destaining after electrophoresis, conducting Northern, Southern and Western blot analysis.

Shaker is a ...

[read more](#)



**3D**  
 BS-010151-AAG  
 Sunflower Mini-Shaker

"Sunflower" **3D** Mini-Shaker provides adjustable three-dimensional smooth rotation of the platform and is designed for mixing blood samples, for minigel staining and destaining, sample washing, blot hybridization reactions.

Mini-Shaker is ...

[read more](#)



**Multi Bio 3D**  
 BS-010125-AAG  
 Programmable mini-shaker

Programmable mini-shaker **Multi Bio 3D** is designed for a variety of applications: hybridization reactions, cell growing, gel washing, soft extraction and homogenisation of biological components in solutions.

Multi Bio 3D ...

[read more](#)