


PV-6/10, platform

DESCRIPTION

Platform for 6 × 10 ml tubes (max diameter 15 mm).

SPECIFICATIONS

Capacity	6
Tube volume	10 ml
Hole diameters	13 mm



CAT. NUMBER

BS-010207-BK

PV-6/10



V-32
Including two platforms PV-32, PL-1
Multi-Vortex

Multi-Vortex V-32 is intended for intensive stirring of bacterial and yeast cell, washing from the culture medium and extraction of metabolites and enzymes from cells and cell cultures.

It is ...

[read more](#)