

SC-18/02, block

DESCRIPTION

20 × 0.2 ml microtubes + 12 × 1.5 ml microtubes

SPECIFICATIONS

Amount of 0.2 ml tubes	20
Amount of 1.5 ml tubes	12
Max. RPM	1400
Max. temperature	100°C



CAT. NUMBER

BS-010120-CK

SC-18/02



TS-100

Without thermoblock

Thermo-Shaker for microtubes and PCR plates

Thermo-Shaker **TS-100** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. Functions of heating (up to +100°C) and mixing can be performed both simultaneously and ...

[read more](#)