

SC-24NC, block

DESCRIPTION

24 x 1.5 ml microtubes

SPECIFICATIONS

Amount of 1.5 ml tubes	24
Max. RPM	1400
Max. temperature	100°C



CAT. NUMBER

BS-010143-GK

SC-24NC

ACCESSORIES



TS-100C

Without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker **TS-100C** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples down ...

[read more](#)



TS-100C Smart

Software included, without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker **TS-100C Smart** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples ...

[read more](#)