

## SC-96AC, block

### DESCRIPTION

96-well unskirted or semi-skirted microplate (0.2 ml) for PCR or 12 × 8 - 0.2ml strips or 96 tubes of 0.2 ml.

### SPECIFICATIONS

|   |       |
|---|-------|
| 96-well unskirted or semi-skirted microplate (0.2 ml) | 1     |
| Max. RPM  | 1400  |
| Max. temperature                                      | 100°C |



### CAT. NUMBER

BS-010143-FK

SC-96AC

### ACCESSORIES



**TS-100C**

Without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker **TS-100C** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples down ...

[read more](#)



**TS-100C Smart**

Software included, without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker **TS-100C Smart** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples ...

[read more](#)