

SV-16/8, platform

DESCRIPTION

16/8/8 sockets for 2*-1.5/0.5/0.2 ml microtest tubes (Ø 11/8/6 mm).

*at 75% volume, max 3400 rpm

SPECIFICATIONS

Diameter	11/8/6 mm
Tube volume	2*-1.5/0.5/0.2 ml
Capacity	16/8/8



CAT. NUMBER

BS-010210-CK

SV-16/8



MSV-3500

With all 4 platforms included
Multi Speed Vortex

Multi Speed Vortex **MSV-3500** is designed for soft or intensive mixing of reagents in different size and type plastic tubes (0.2 to 50 ml).

It is designed for operation in ...

[read more](#)