

## TR-5/30, rack

### DESCRIPTION

Test tube rack for  $\varnothing 30$ mm tubes, capacity 5 tubes. WxDxH: 155x90x112mm

### SPECIFICATIONS

Capacity	5
Tube size	$\varnothing 30$ mm



### CAT. NUMBER

BS-010406-KK	Test tube rack for $\varnothing 30$ mm tubes, capacity 5 tubes
--------------	--



**WB-4MS**  
With base BP-1 and lid  
Stirred water bath

Water bath-thermostat **WB-4MS** is designed for chemical, pharmaceutical, medical and biological laboratory research.

**WB-4MS** provides increased temperature stabilization (up to 0.1°C) due to built-in magnetic stirrer (regulated speed 250–1,000 rpm).

...

[read more](#)