

FB


Serana Secoll™ Lymphocyte Separation Solution

PRODUCT LAUNCH

Separation of Cells & Sub-Cellular Components

ANIMAL COMPONENT FREE

- Secoll™ Lymphocyte separation solution is a sterile, ready-to-use density gradient medium designed for the efficient isolation of viable lymphocytes from whole blood, bone marrow, or other biological samples.
- This solution exploits the principle of density gradient centrifugation to selectively separate lymphocytes from erythrocytes & granulocytes.



Catalog No.	Description	Vol	Price (ex GST)
RLL-001-100ML	Secoll™ Lymphocyte Separation Solution Density 1.077 g/ml, Sterile Filtered	100mL	\$40
RLL-001-500ML	Secoll™ Lymphocyte Separation Solution Density 1.077 g/ml, Sterile Filtered	500mL	\$155

Terms and Conditions

All prices are quoted in AUD and exclude GST.
Offers exclude servicing, spare parts
and freight.