

## P-5, shelf for pipettes

### DESCRIPTION

P-5, shelf for pipettes


### SPECIFICATIONS

Overall dimensions (W×D×H)

230 × 140 × 50 mm

Capacity

5



### CAT. NUMBER

BS-040104-CK

P-5

### ACCESSORIES



#### UVC/T-M-AR

With built in socket and inlet for power cords  
PCR cabinet

The Biosan **UVC/T-M-AR** is a benchtop DNA/RNA UV-cleaner workstation designed to create a clean, contamination-free environment for PCR setup and other nucleic acid-related workflows.

This model features a dual UV ...

[read more](#)



#### UVT-B-AR

With built in socket and inlet for power cords  
PCR cabinet

The Biosan **UVT-B-AR** is a benchtop DNA/RNA UV-cleaner workstation designed to create a clean, contamination-free environment for PCR setup and other nucleic acid-related workflows.

This model features a dual UV ...

[read more](#)



#### UVT-S-AR

BS-040107-A07  
PCR cabinet

The Biosan **UVT-S-AR** is a benchtop DNA/RNA UV-cleaner workstation designed to create a clean, contamination-free environment for PCR setup and other nucleic acid-related workflows.

This model features a dual UV ...

[read more](#)



#### UVC/T-M-AR

With built in socket and inlet for power cords  
DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box **UVC/T-M-AR** is designed for clean operations with DNA samples. UV-cleaner box provide protection against contamination.

Model is a bench-top type, made of metal framework, glass walls, working ...

[read more](#)



#### UVT-B-AR

With built in socket and inlet for power cords  
DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box **UVT-B-AR** is designed for clean operations with DNA samples. UV-Cabinet provide protection against contamination.

Model is a bench-top type, made of metal framework, working surface made of ...

[read more](#)



#### UVT-S-AR

BS-040107-AAA  
DNA/RNA UV-cleaner box

DNA/RNA UV-cleaner box **UVT-S-AR** is designed for clean operations with DNA samples. UV-cleaner box provide protection against contamination.

Model is a bench-top type, made of metal framework, glass walls, working ...

[read more](#)