

Remote grade 1 water dispenser

DESCRIPTION

Ultrapure water dispenser is designed to dispense ultrapure water that complies with ISO 3696 Grade I water requirements.

- The control of the dispenser is fulfilled using the control panel and a dispense button on the water dispenser.
- Dispenser is connected to water purification system with two 1/4" JG tubes.
- The dispenser is mounted on a stem which allows rotation around the axis by 360 degrees and the height adjustment.
- If desired, the dispenser can be removed from the stem for easier water dispense.



CAT. NUMBER

BS-070110-AK

Remote grade 1 water dispenser

ACCESSORIES



Labaqua Trace
BS-070105-A02
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

Any configuration of a Labaqua ultrapure system produces both ultrapure ...

[read more](#)



Labaqua HPLC
BS-070104-A02
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

Labaqua HPLC produces water with very low organic carbon (TOC) ...

[read more](#)



Labaqua Bio
30l storage tank
ultrapure water system

Labaqua ultrapure systems are multi-purpose water purification systems. The Labaqua systems produce ultrapure and pure water directly from tap water.

Labaqua Bio system produces water with very low organic and ...

[read more](#)