


P-EX, platform exchange set

DESCRIPTION

Platform Exchange Set allows for faster replacement of flat platforms. Thanks to the platform sliding system, it is now possible to install 4 × 2l flasks in ES-20/60, ES-20/80, ES-20/80C shakers-incubators. The platform exchange set can accommodate UP-168, P-6/1000, P-9/500, P-15/250, P-30/100, HSP-6/1000, HSP-9/500, HSP-15/250, HSP-30/100 platforms.



CAT. NUMBER

BS-010173-CK

P-EX (Slider + Tray)

ACCESSORIES



P-30/100
 BS-010135-BK
 platform

Platform with clamps for flasks.



P-16/250
 BS-010135-CK
 platform

Platform with clamps for flasks.



P-9/500
 BS-010135-AK
 platform

Platform with clamps for flasks.



P-6/1000
 BS-010135-DK
 platform

Platform with clamps for flasks.



HSP-6/1000
 BS-010167-LK
 platform

Platform with tight fit clamps for flasks.



HSP-9/500
BS-010167-NK
platform

Platform with tight fit clamps for flasks.



HSP-16/250
BS-010167-MK
platform

Platform with tight fit clamps for flasks.



HSP-30/100
BS-010167-KK
platform

Platform with tight fit clamps for flasks.



UP-168
BS-010135-JK
platform

Universal platform with clamps can accommodate flasks or bottles of different volume sizes. Clamps are not included with platform and need to be ordered separately.



ES-20/60
Without platform
Orbital Shaker-Incubator

Orbital Shaker-Incubator ES-20/60 for biotechnological and pharmaceutical laboratories is a professional category equipment designed for cultivation of microorganisms and eukaryotic cells including animal, plant and insect cells. It is also ...

[read more](#)



ES-20/80
Software included, without platform
Orbital Shaker-Incubator

ES-20/80 shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

[read more](#)



ES-20/80C
Software included, without platform
Orbital Shaker-Incubator with cooling

ES-20/80C shaker-incubator for biotechnological and pharmaceutical laboratories is a professional category equipment. The typical applications include - microbial and cell culture cultivation, protein expression, solubility studies, general mixing, as well ...

[read more](#)