



SHIMAS[®]

Glass Bead Sterilisers

STERIBEAD 250P

TISSUE CULTURE LAB EQUIPMENT

Glass bead sterilisers suitable for sterilising small instruments used in plating, like inoculating loops & spreaders

\$529 ex GST

- Also used in plant tissue culture for sterilising tweezers & needles
- Replaces Bunsen Burners & hot air sterilisers which can cause air flow displacement in BSCs & LAF cabinets
- Crucible Size: 50 mm Dia x 90mm Depth
- Sterilisation temperature is factory preset at 250-265oC
- Stainless Steel construction



This offer is valid whilst stocks last

Catalog No.	Description	Price (ex GST)
STERIBEAD250P	Glass Bead Steriliser, 250-265°C, 50x90mm Crucible	\$529

