



Sorenson Bioscience™

## Gel Loading Tips

**BULK BAGGED**

MiniFlex Tips feature an elongated ultra-thin microcapillary that is ideal for loading gels

**IDEAL FOR SDS-GELS**

- 0.5 µL - 200 µL volume range
- Graduated at 5 µL for quick volume check
- Fits standard universal pipettors
- Available sterile 83mm length for access to the bottom of narrow vessels and gel wells
- Round version is ideal for SDS-PAGE gels



Catalog No.	Description	Vol	Price (ex GST)
SOR-28480-1	Sorenson Clear Gel Loading Tip, Round, 0.57mm Outside Diameter	1000 Tips/Bag	\$53