

OneTouch Tips

A truly universal tip that fits virtually any pipettor – OneTouch tips provide an innovative, ergonomic and environmentally friendly liquid handling solution









Innovative – Pipetting like you’ve never felt it

A revolutionary tip that is manufactured from two resins, OneTouch tips provide a soft color-coded top collar and a rigid, ultra-clear polypropylene body.

Ergonomic – Soft tip requires less pressure to create an optimal seal

Tips reduce the necessary sealing and reciprocal ejection force felt with traditional single resin polypropylene tips. Tips are excellent for multichannel pipettors, and reduce pipettor wear.



	Cat. No.	Vol. Range	Type	Graduations	Packaging	Pack Type	Quantity
10 µl Barrier							
	10320	1-20 µl	Barrier	1 µl, 5 µl, 10 µl	Sterile, 96 Tips/Rack	 Cell Rack	10 Racks/Pack, 5 Packs/Case
200 µl Barrier							
	10350	1-330 µl	Barrier	10 µl, 50 µl, 100 µl	Sterile, 96 Tips/Rack	 Cell Rack	10 Racks/Pack, 5 Packs/Case
1000 µl Barrier							
	10380	50-1250 µl	Barrier	200 µl, 500 µl, 1000 µl	Sterile, 96 Tips/Rack	 Hinged Rack	8 Racks/Pack, 5 Packs/Case

FB
AUSTRALIA

Australian Distributors:

FB AUSTRALIA

Free Call: 1800 066 077

Email: info@fisherbiotec.com

Web: fisherbiotec.com.au