

# Axygen® MultiRack Tip System

## Filter & Reload Promo



The Axxygen MultiRack tip system has been engineered to fit most of the pipettors in the market.

It is an interchangeable, complete system offered in bulk, racked, filter, and reload packaging styles to meet the variety of needs in the lab setting.

The MultiRack tip system offers standard and extended-length tips allowing volumes from 0.1 µL to 1250 µL.

Axygen MultiRack tip system products are made in the USA and meet ISO standards. All tips are free of detectable RNase, DNase, DNA, and pyrogens.

### Features

- ▶ Reusable rack with reload systems
- ▶ Small footprint and stackable design for maximum bench space usage
- ▶ Extended tips offered in 10, 200, 300, and 1000 µL sizes
- ▶ Color-coded racks for easy volume identification

### Ordering Information

#### Axygen® MultiRack Maxymum Recovery® Pipet Tips

##### Racked, Filter, Sterile, Maxymum Recovery Surface

Cat. No.	Description	Quantity	Price (ex GST)
AX-MRF-10XT-L-R-S-1	10 µL MultiRack extended length filtered pipet tip, Maxymum Recovery surface, racked, sterile	960/pk	\$69
AX-MRF-20-L-R-S-1	20 µL MultiRack filtered pipet tip, Maxymum Recovery surface, racked, sterile	960/pk	\$69
AX-MRF-100-L-R-S-1	100 µL MultiRack filtered pipet tip, Maxymum Recovery surface, racked, sterile	960/pk	\$69
AX-MRF-200NX-L-R-S-1	200 µL NX MultiRack filtered pipet tip, Maxymum Recovery surface, racked sterile	960/pk	\$69
AX-MRF-200XT-L-R-S-1	200 µL MultiRack extended length filtered pipet tip, Maxymum Recovery surface, racked, sterile	768/pk	\$69
AX-MRF-300XT-L-R-S-1	300 µL MultiRack extended length filtered pipet tip, Maxymum Recovery surface, racked, sterile	768/pk	\$69
AX-MRF-1000XT-L-R-S-1	1000 µL MultiRack extended length, filtered pipet tip, Maxymum Recovery surface, racked, sterile	768/pk	\$69

#### Axygen MultiRack Reload System

Axygen MultiRack Reload System is designed to fit the Axxygen MultiRack System without tedious re-alignment. This design not only saves lab space but also reduces environmental impact. Sturdy racks resist accidental tipping, and the easy-close lid lock prevents tips from falling out of the rack while closed.

Cat. No.	Description	Quantity	Price (ex GST)
MR-10XT-RL-1	10 µL MultiRack extended length pipet tip, reload	960/pk	\$36
MR-200-RL-1	200 µL MultiRack pipet tip, reload	960/pk	\$36
MR-300-RL-1	300 µL MultiRack pipet tip, reload	960/pk	\$36
MR-1000XT-RL-1	1000 µL MultiRack extended length pipet tip, reload	768/pk	\$36

#### Accessories

Cat. No.	Description	Quantity	
MR-R-A	Empty rack for MultiRack 10 µL extended length, 200 µL, 200 µL NX, 300 µL	5/pk, 1 pk/cs	\$18
MR-R-B	Empty rack for MultiRack 200 µL extended length, 300 µL extended length, and 1000 µL extended length	5/pk, 1 pk/cs	\$20

Please note: Introductory promotion price valid until 31st Dec 2025



Australian Distributors:

FB AUSTRALIA

Free Call: 1800 066 077

Email: [info@fisherbiotec.com](mailto:info@fisherbiotec.com)

Web: [fisherbiotec.com.au](http://fisherbiotec.com.au)